



orthopedic **ESSENTIALS**  
product catalog

[www.deroyal.com](http://www.deroyal.com)

## TABLE OF CONTENTS:

Dear Valued Customer:

For over forty years, we have been privileged to stand quietly behind you in your emergency departments, orthopedic hospital floors, healthcare clinics, and physical therapy/rehab centers. We are proud to be one of the largest suppliers of orthopedic softgoods to hospitals and clinics in the nation. The diversity of our product line, the quality inherent in our manufacturing process and our committed staff enable us to offer our customers exceptional products at reasonable prices.

Our commitment continues - a commitment to providing patients with the best products available today. Using our own raw materials and many modern manufacturing techniques, we make comfortable, affordable products which improve your patient's quality of life.


We intend for this catalog to provide you with the information necessary to order and use our products. For more information, please contact your local DeRoyal Sales Representative or visit us online at [www.deroyal.com](http://www.deroyal.com). To place an order, fax us at 800.543.2182, e-mail [customerservice@deroyal.com](mailto:customerservice@deroyal.com) or call Customer Service directly at 800.251.9864.

Finally, as in any successful partnership, please let us know how we can best serve you. We strive for continuous improvement based on customer input, and we very much look forward to hearing from you.



Pete DeBusk

Founder and Chairman, DeRoyal



## SPINE

<b>Head/Neck</b>	<b>5</b>
CERVICAL COLLARS   Soft Foam   Rigid   Post-Op	
CERVICAL TRACTION	
<b>Back</b>	<b>13</b>
BACK BRACING	
LUMBAR SACRAL SUPPORTS	

## UPPER EXTREMITY

<b>Torso</b>	<b>19</b>
ABDOMINAL BINDERS	
RIB BELTS	
<b>Clavicle /Shoulder</b>	<b>23</b>
CLAVICLE SUPPORTS	
SHOULDER IMMOBILIZERS   Sling & Strap   Sling & Swathe   Velpeau   Shoulder Abductors	
<b>Arm/Elbow</b>	<b>31</b>
ARM SLINGS/ELEVATORS   Arm Slings   Elevators	
ELBOW BRACES & SUPPORTS   Post-Op Elbow   Supports & Sleeves	
<b>Wrist/Hand</b>	<b>41</b>
Wrist & Wrist/Forearm   Wrist & Thumb   Aluminum Splints	

## LOWER EXTREMITY


<b>Knee/Leg</b>	<b>57</b>
KNEE IMMOBILIZERS   Tri-Panel   Cutaway   Sized	
POST-OP KNEE	
KNEE SUPPORTS	
BUCKS TRACTION BOOTS	
<b>Ankle/Foot &amp; Heel</b>	<b>71</b>
ANKLE BRACES & SUPPORTS   Stirrups	
AMUBULATORY AIDS	

## PATIENT CARE

PRESSURE ULCER PREVENTION	<b>87</b>
FOAM POSITIONERS	
SPECIALTY NURSING	
<b>Temperature Therapy</b>	<b>101</b>
JETSTREAM & COOLJET	
WRAPS / BLANKETS	
CHILL-N-WRAP / HOT AND COLD PACK	



# SPINE | Head/Neck



## CONTOURED CERVICAL COLLAR soft foam



- Medium density foam construction
- Breathable cotton stockinette liner
- Additional cover included for better patient hygiene
- Deep contour holds the head in a neutral position

Size	Length	2½"	3"	4"
Universal	22"	1061-00	1057-00	1056-00
X-Small	13"	1061-01	1057-01	—
Small	15"	1061-02	1057-02	1056-02
Medium	17"	1061-03	1057-03	1056-03
Large	18"	1061-04	1057-04	1056-04
X-Large	20"	1061-05	1057-05	1056-05

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

## MEDIUM DENSITY CERVICAL COLLAR soft foam



- Medium density foam construction
- Breathable cotton stockinette liner
- Additional cover included for better patient hygiene

SIZE	LENGTH	2½"	3"	3½"	4"
Universal	22"	1000251	1000301	1000351	1000401
Small	16"	1010251	1010301	1010351	1010401
Medium	18"	1020251	1020301	1020351	1020401
Large	20"	1030251	1030301	1030351	1030401
X-Large	22"	1040251	1040301	1040351	1040401

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

## MEDIUM/FIRM DENSITY CERVICAL COLLAR soft foam



- Medium/firm density foam construction
- Breathable cotton stockinette liner
- Additional cover included for better patient hygiene
- Contour holds the head in a neutral position

SIZE	LENGTH	2½"	3"	3½"	4"
Universal	22"	1050-00	1050-10	1050-20	1050-30
Small	16"	1050-01	1050-11	1050-21	1050-31
Medium	18"	1050-02	1050-12	1050-22	1050-32
Large	20"	1050-03	1050-13	1050-23	1050-33
X-Large	22"	1050-04	1050-14	1050-24	1050-34

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with and without radiculopathy

## COMFO-EZE CERVICAL COLLAR soft foam



- Highest quality medium density foam construction
- Breathable synthetic stockinette liner
- Additional cover included for better patient hygiene
- Available in wide and narrow styles

SIZE	DIMENSIONS	WIDE
X-Small	3½" x 19½"	1001-01
Small	4" x 20½"	1001-02
Medium	4" x 22½"	1001-03
Large	4" x 23½"	1001-04
X-Large	4½" x 25½"	1001-05
Universal	4¾" x 22"	1010-01

SIZE	DIMENSIONS	REGULAR
X-Small	3" x 19½"	1002-01
Small	3½" x 20½"	1002-02
Medium	3½" x 22½"	1002-03
Large	3½" x 23½"	1002-04
X-Large	3½" x 25½"	1002-05
Universal	3½" x 22"	1010-02

SIZE	DIMENSIONS	NARROW
X-Small	2" x 19½"	1003-01
Small	2" x 20½"	1003-02
Medium	2" x 22½"	1003-03
Large	2" x 23½"	1003-04
X-Large	2" x 25½"	1003-05

SIZE	DIMENSIONS	PRODUCT #
Pediatric Small	2½" x 14"	1004-01
Pediatric Large	3½" x 17"	1004-02

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with or without radiculopathy

## SERPENTINE CERVICAL COLLAR soft foam



- Firm density foam construction
- Breathable synthetic stockinette liner
- Additional cover included for better patient hygiene
- Contoured serpentine design for patient comfort

SIZE	DIMENSIONS	PRODUCT #
Small	3½" x 17"	1006-02
Small, X-Wide	3¾" x 18"	1006-01
Small, X-Long	3¼" x 24½"	1006-05
Medium	4" x 20"	1006-03
Large	4" x 24"	1006-04

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with or without radiculopathy

## SOFT FOAM CERVICAL COLLAR soft foam



- Soft density foam construction
- Breathable cotton stockinette liner
- Packaged 6/Case

SIZE	DIMENSIONS	PRODUCT #
X-Small	3" x 19½"	1018-01
Small	3½" x 20½"	1018-02
Medium	3½" x 22½"	1018-03
Large	3½" x 23½"	1018-04
X-Large	3½" x 25½"	1018-05
Universal	3" x 22"	1018-10
Universal Narrow	3½" x 22"	1018-11
Universal Wide	4¼" x 22"	1018-12

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with or without radiculopathy

## ECONOMY FOAM COLLAR soft foam



- Soft density foam construction
- Breathable cotton stockinette liner

SIZE	DIMENSIONS	PRODUCT #
Universal	2" x 24"	A113000
Universal	3" x 24"	A115000
Universal	3.5" x 24"	A116000
Universal	4" x 24"	A117000
Universal	4.5" x 24"	A118000

SIZE	LENGTH	2"	3"	3.5"	4"
Pediatric	16"	A113002	—	—	—
Small	16"	—	A115005	A116005	A117005
Medium	18"	—	A115006	A116006	A117006
Large	20"	—	A115007	A116007	A117007
X-Large	22"	—	A115008	A116008	A117008

### Indications:

Cervical facet syndrome; Cervical sprains and strains; Discogenic pain; Degenerative disc disease with or without radiculopathy

## DEROYAL® RIGID COLLAR rigid



- 2-piece design restricts flexion, extension and rotation while supporting the occipital and mandibular areas
- Constructed of Plastazote® foam
- Large trachea opening allows quick access for carotid pulse monitoring, emergency airway procedures and visualization
- Posterior ventilation holes help reduce heat and moisture
- Not made with natural rubber latex

SIZE	CIRCUMFERENCE	1.25"	1.75"	2.25"	3.25"	4.25"	5.25"
Infant	6"-8"	A9917-01	—	—	—	—	—
Pediatric	8"-11"	—	A9917-00	—	—	—	—
Small	10"-13"	—	—	A9918-11	A9919-11	A9916-11	A9920-11
Medium	13"-16"	—	—	A9918-12	A9919-12	A9916-12	A9920-12
Large	16"-19"	—	—	A9918-13	A9919-13	A9916-13	A9920-13
X-Large	19" & up	—	—	A9918-14	A9919-14	A9916-14	A9920-14

### Indications:

Provides superior immobilization for both extrication and rehabilitation; Post-operative discectomy and cervical fusion



## EMT SELECT COLLAR extrication



- 1-piece flat design allows efficient storage where space is limited
- Universal/adjustable to four adult and three pediatric sizes
- Locking tabs, slide adjustment and reflective sizing indicators for maximum hold
- Closed cell foam liner to help prevent skin breakdown
- Large trachea opening allows easy access for carotid pulse monitoring, emergency airway procedures and visualization

SIZE	FRONT	CIRCUMFERENCE	QTY	PRODUCT #
Adult	3½"	11" - 23"	1 Each	1044-31
			20/Case	1044-30
			100/Case	1044-90
Pediatric	1½"	8" - 18"	1 Each	1044-41
			20/Case	1044-40

### Indications:

Provides superior immobilization for both extrication and rehabilitation;  
Cervical trauma

## ONE-PIECE EXTRICATION COLLAR extrication



- 1-piece flat design allows efficient storage where space is limited
- Four adult and two pediatric sizes ensure accurate fit
- Color-coded for quick size reference
- Large trachea opening allows easy access for carotid pulse monitoring, emergency airway procedures and visualization

SIZE	FRONT	CIRCUMFERENCE	PRODUCT #
Infant	1½"	8" - 18"	1044-80
Pediatric	1¾"	8" - 18"	1044-81
Low	1¾"	11" - 23"	1044-82
Short	2½"	11" - 23"	1044-83
Regular	3"	11" - 23"	1044-84
Tall	3½"	11" - 23"	1044-85

### Indications:

Provides superior immobilization for both extrication and rehabilitation;  
Cervical trauma

## XTW EXTENDED WEAR COLLAR post-op



- Designed to help restrict flexion, extension, lateral bending and axial rotation
- High density polyethylene frame for maximum rigidity and support
- Shoulder and mandible contours for increased stability
- Anterior and posterior flex tabs for individual comfort
- Open-cell liner helps prevent skin breakdown
- Five adult and five child sizes to ensure best patient fit
- Replacement pad set available upon request
- Not made with natural rubber latex

SIZE	COLLAR	COLLAR & PAD SET
Adult X-Short	11201000	11202000
Adult Short	11211000	11212000
Adult Regular	11221000	11222000
Adult Tall	11231000	11232000
Adult X-Tall	11241000	11242000
Infant Short	11251000	11252000
Infant Regular	11261000	11262000
Child Tot	11271000	11272000
Child Short	11281000	11282000
Child Regular	11291000	11292000

### Indications:

Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

## CAPITAL™ CERVICAL COLLAR post-op



- Designed to help restrict flexion, extension, lateral bending and axial rotation
- Reinforced sternal padding for optimum support and comfort
- Shoulder and mandible contours for increased stability
- Sheer finished open-cell liner is specifically designed to help prevent skin breakdown and provide sustained comfort during long-term use
- Replacement pad set available upon request
- Not made with natural rubber latex

SIZE	COLLAR	COLLAR & PAD SET	PAD SET
Adult Super Short	11421000	11422000	11423000
Adult Stout	11431000	11432000	11433000
Adult X-Small	11441000	11442000	11443000
Adult Short	11451000	11452000	11423000
Adult Regular	11461000	11462000	11423000
Adult Large	11471000	11472000	11423000
Infant 0-6 Months	11481000	11482000	11483000
Infant 6 Months – 2 Years	11481001	11482001	11483001
Child 2-6 Years	11481002	11482002	11483002
Child 6-12 Years	11481003	11482003	11483003

### Indications:

Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

## CAPITAL™ COLLAR ENHANCED post-op POWERED BY BOA® TECHNOLOGY



- Patented design featuring a micro-adjustable dial powered by Boa® technology that allows collar height modification in order to fit most patients
- Includes compliance ring to promote patient compliance
- Reinforced sternal padding for optimum support and comfort
- Designed to help restrict flexion, extension, lateral bending and axial rotation
- Shoulder, mandible and occipital contours for increased stability
- Sheer finished open-cell liner is specifically designed to help prevent skin breakdown and provide sustained comfort during long-term use
- Not made with natural rubber latex

SIZE	COLLAR	COLLAR & PAD SET	PAD SET ONLY
Universal	11311000	11312000	11313000

### Indications:

Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

## DELUXE HEAD HALTER cervical traction



- Constructed of perforated foam material with tricoot lining
- Designed to provide anterior tilt cervical traction
- Halter is easily adjustable with hook and loop closures
- Spreader bar not included

SIZE	PRODUCT #
Small	6101-01
Medium	6101-02
Large	6101-03

### Indications:

Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

## REPRO TYPE UNIVERSAL HEAD HALTER cervical traction



- Made of durable canvas material with tricoot lining
- Pulls from the occipital area to eliminate chin pressure
- Halter tilts head slightly during traction
- Chin straps have buckle closures for easy adjustment

SIZE	w/o SPREADER BAR
Universal	6102-00

### Indications:

Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

## DISCARD HEAD HALTER cervical traction



- Available in cotton or flannel material
- Product can be applied for hyperextension or flexion
- Spreader bar not included
- Packaged 12/Case

SIZE	COTTON	FLANNEL
Universal	6104-00	6105-00

### Indications:

Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

## STANDARD OVER-DOOR TRACTION UNIT cervical traction



- Unit comes with standard discard flannel head halter
- 12" spreader bar, traction cord and vinyl water bag included
- Fits easily on any door

SIZE	W/ 10 LB. BAG	W/ 20 LB. BAG
Universal	7102-00	7102-22

### Indications:

Degenerative disc disease; Cervical facet syndrome; Herniated nucleus pulposus with or without radiculopathy; Chronic cervical strains

## TRACTION CORD cervical traction



- Inner core constructed of polyester
- Polypropylene outer weave
- 600 pound tensile strength
- ¼" diameter

LENGTH	PRODUCT #
100'	7116-01

### Indications:

Traction

## WEIGHT BAG cervical traction




- Durable vinyl bag
- Graduated marks for weight adjustment

SIZE	PRODUCT #
20 lb.	7117-02

### Indications:

Traction

# SPINE | Back





## SOLACE® LOW PROFILE back bracing



- Soft, breathable material provides cool, comfortable fit increasing patient compliance
- Adjustable pulley style closure system is easy to secure and provides maximum compression
- Anterior and posterior plastic inserts provide control and support
- Low profile, lightweight design fits easily under patient's clothing
- Plastic D-rings for ease of application and removal

SIZE	WAIST		PRODUCT #
	CIRCUMFERENCE		
Small	26" - 34"		13400005
Medium	32" - 40"		13400006
Large	38" - 46"		13400007
X-Large	44" - 51"		13400008

### Indications:

Chronic low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic lumbar pain

## PROLIGN® LUMBAR ORTHOSIS back bracing



- Constructed of hook receivable material with soft breathable liner that wicks away moisture to help prevent skin breakdown
- Patented dual dynamic closure system provides instant full range adjustment and maximum compression
- Rigid anterior insert provides support to intra-abdominal and viscera region helping to unload the spine under gravitational load
- Available in 15° and 25° spot moldable posterior plastic inserts
- Plastic D-rings for ease of application and removal

SIZE	WAIST		
	CIRCUMFERENCE	15°	25°
X-Small	25" - 30"	13501004	13500004
Small	30" - 35"	13501005	13500005
Medium	35" - 40"	13501006	13500006
Large	40" - 45"	13501007	13500007
X-Large	45" - 50"	13501008	13500008
XX-Large	50" - 55"	13501009	13500009

### Indications:

Post-operative lumbar immobilization; Chronic low back pain; Lumbar muscle weakness; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; IDET procedure; Spinal blocks

## PROLIGN® EXT LUMBAR SACRAL ORTHOSIS back bracing



- Soft, breathable material to provide a cool, comfortable fit
- Dual, dynamic closure system provides infinite adjustment and maximum compression
- Rigid anterior insert provides support to the intra-abdominal and viscera region helping to unload the spine under gravitational load
- Posterior plastic insert (spot moldable) helps to maintain neutral sagittal alignment while maximizing support of anatomical structures
- Plastic D-rings for ease of application and removal
- Extended posterior panel for coverage from sacroccocygeal junction to the T9 vertebrae

SIZE	WAIST		SIZE	WAIST	
	CIRCUMFERENCE	25°		CIRCUMFERENCE	25°
X-Small	25" - 30"	155-01	X-Large	45" - 50"	155-05
Small	30" - 35"	155-02	XX-Large	50" - 55"	155-06
Medium	35" - 40"	155-03	XXX-Large	55" - 60"	155-07
Large	40" - 45"	155-04	XXXX-Large	60" - 65"	155-08

### Indications:

Post-operative lumbar immobilization; Chronic low back pain; Lumbar muscle weakness; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; IDET procedure; Spinal blocks

# ULTRALIGN®+ LUMBAR SACRAL ORTHOSIS back bracing



- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows infinite adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 0°, 15° or 25° lordosis
- Non-Tapered (male) and tapered (female) versions available for best patient fit
- Mid and low profile versions for short-waisted patients
- LSO meets requirements for most L1 to L5 pathologies

## Indications:

Post-operative lumbar fusion; Laminectomy or discectomy; Lumbar compression fractures; Degenerative disc disease; Osteoporosis; Spondylolisthesis; Spondylolysis; Spinal stenosis; Facet syndrome

## Sizing Guide

Sizing determined by three circumferential measurements:

1. 1" below tip of breast bone or xyphoid process
2. Around navel or waist
3. Around fullest area of trochanter or hips

NON-TAPERED	XYPHOID	WAIST	TROCHANTER
Small	31" - 35"	31" - 35"	31" - 35"
Medium	34" - 38"	32" - 38"	34" - 39"
Large	37" - 41"	35" - 41"	38" - 42"
X-Large	38" - 45"	38" - 45"	38" - 45"
XX-Large	45" - 53"	45" - 53"	45" - 53"

TAPERED	XYPHOID	WAIST	TROCHANTER
Small	25" - 31"	25" - 31"	27" - 33"
Medium	28" - 35"	29" - 36"	32" - 38"
Large	32" - 38"	32" - 38"	35" - 41"
X-Large	35" - 42"	36" - 42"	38" - 46"
XX-Large	40" - 46"	39" - 47"	43" - 49"

PRODUCT #	0°	15°	25°
<b>NON-TAPERED</b>			
Small	UPL1330-0	UPL1330-15	UPL1330-25
Medium	UPL1430-0	UPL1430-15	UPL1430-25
Large	UPL1530-0	UPL1530-15	UPL1530-25
X-Large	UPL1630-0	UPL1630-15	UPL1630-25
XX-Large	UPL1730-0	UPL1730-15	UPL1730-25

TAPERED	0°	15°	25°
Small	UPL2330-0	UPL2330-15	UPL2330-25
Medium	UPL2430-0	UPL2430-15	UPL2430-25
Large	UPL2530-0	UPL2530-15	UPL2530-25
X-Large	UPL2630-0	UPL2630-15	UPL2630-25
XX-Large	UPL2730-0	UPL2730-15	UPL2730-25

LOW PROFILE TAPERED	15°	25°
Small	UPLP2330-15	UPLP2330-25
Medium	UPLP2430-15	UPLP2430-25
Large	UPLP2530-15	UPLP2530-25
X-Large	UPLP2630-15	UPLP2630-25
XX-Large	UPLP2730-15	UPLP2730-25

MID PROFILE NON-TAPERED	15°
Small	UPMP1330-15
Medium	UPMP1430-15
Large	UPMP1530-15
X-Large	UPMP1630-15
XX-Large	UPMP1730-15

# ULTRALIGN®+ THORACIC LUMBAR SACRAL ORTHOSIS back bracing



**SPK**



**DLE**



Low Profile

- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows infinite adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 15° or 25° lordosis
- Non-tapered (male) and tapered (female) versions available for best patient fit
- DLE (Dorsal Lumbar Extension) or SPK (Sternal Pad Kit) available for maximum control
- Low profile version for short waisted patients
- TLSO (Thoracic Lumbar Sacral Orthosis) meets requirements for most T6 to L5 pathologies

## Indications:

Post-operative thoracic/lumbar fusion; Laminectomy or discectomy; Compression fractures; Degenerative disc disease; Osteoporosis; Single column instability immobilization; Facet syndrome

## Sizing Guide

Sizing determined by three circumferential measurements:

1. 1" below tip of breast bone or xiphoid process
2. Around navel or waist
3. Around fullest area of trochanter or hips

NON-TAPERED	XYPHOID	WAIST	TROCHANTER
Small	31" - 35"	31" - 35"	31" - 35"
Medium	34" - 38"	32" - 38"	34" - 39"
Large	37" - 41"	35" - 41"	38" - 42"
X-Large	38" - 45"	38" - 45"	38" - 45"
XX-Large	45" - 53"	45" - 53"	45" - 53"

TAPERED	XYPHOID	WAIST	TROCHANTER
Small	25" - 31"	25" - 31"	27" - 33"
Medium	28" - 35"	29" - 36"	32" - 38"
Large	32" - 38"	32" - 38"	35" - 41"
X-Large	35" - 42"	36" - 42"	38" - 46"
XX-Large	40" - 46"	39" - 47"	43" - 49"

## Product #

NON-TAPERED	15°	25°
Small	UPT1333-15	UPT1333-25
Medium	UPT1433-15	UPT1433-25
Large	UPT1533-15	UPT1533-25
X-Large	UPT1633-15	UPT1633-25
XX-Large	UPT1733-15	UPT1733-25

TAPERED	15°	25°
Small	UPT2333-15	UPT2333-25
Medium	UPT2433-15	UPT2433-25
Large	UPT2533-15	UPT2533-25
X-Large	UPT2633-15	UPT2633-25
XX-Large	UPT2733-15	UPT2733-25

LOW PROFILE TAPERED	15°	25°
Small	UPLT2333-15	UPLT2333-25
Medium	UPLT2433-15	UPLT2433-25
Large	UPLT2533-15	UPLT2533-25
X-Large	UPLT2633-15	UPLT2633-25
XX-Large	UPLT2733-15	UPLT2733-25

## ACCESSORIES

- Dorsal Lumbar Extension
- Sternal Pad Kit
- 3 1/2" Extension Panels
- 5 1/2" Extension Panels

## PRODUCT #

- UDLE
- USPK
- UEP2031
- UEP2051



## CRISS-CROSS LUMBO SACRAL SUPPORT lumbar supports



- Durable elastic construction with hook and loop closure
- Contoured waist for maximum comfort
- Two elastic side pulls for optimal compression
- Four adjustable posterior stays for stability and support
- 12" posterior heights
- Not made with natural rubber latex

SIZE	WAIST		SIZE	WAIST	
	CIRCUMFERENCE	PRODUCT #		CIRCUMFERENCE	PRODUCT #
Universal	26" - 46"	13850000	X-Large	40" - 44"	13850008
Small	24" - 28"	13850005	XX-Large	44" - 48"	13850009
Small/Medium	28" - 32"	13850056	XXX-Large	48" - 52"	13850010
Medium/Large	32" - 36"	13850067	XXXX-Large	52" - 56"	13850011
Large	36" - 40"	13850007			

### Indications:

Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic thoracolumbar pain

## ELASTIC LUMBO SACRAL SUPPORT lumbar supports



- Constructed of durable elastic with hook and loop closure
- Two elastic side pulls allow optimal compression
- Six flexible posterior stays for stability and support
- 10" posterior height
- Not made with natural rubber latex

SIZE	WAIST	
	CIRCUMFERENCE	PRODUCT #
Universal	26" - 46"	13810000

### Indications:

Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic thoracolumbar pain

## HEAT & MOLD LUMBO SACRAL SUPPORT lumbar supports



- Not made with natural rubber latex
- Posterior pocket accommodates moldable thermoplastic insert (sold separately)
- 8" posterior height

SIZE	WAIST	
	CIRCUMFERENCE	PRODUCT #
Small	24" - 32"	13800005
Medium	32" - 38"	13800006
Large	38" - 46"	13800007
X-Large	46" - 54"	13800008
Moldable Insert		13860900

### Indications:

Moderate to severe low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Lumbar instability

## ELASTIC/NEOPRENE BACK SUPPORT lumbar supports



- Constructed of beige or black elastic material
- Contoured shape for increased comfort in seated position
- Two dynamic side pulls for optimal compression
- Neoprene posterior pocket accommodates moldable thermoplastic insert (sold separately)
- 9" posterior height

SIZE	WAIST CIRCUMFERENCE		BEIGE	BLACK
Small	26" - 30"		NE7741-02	NE7741-12
Medium	30" - 34"		NE7741-03	NE7741-13
Large	34" - 38"		NE7741-04	NE7741-14
X-Large	38" - 42"		NE7741-05	NE7741-15
XX-Large	42" - 46"		NE7741-06	NE7741-16
XXX-Large	46" - 50"		NE7741-07	NE7741-17
Moldable Insert			13860900	

### Indications:

Mild to moderate low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Lumbar instability

## WORKERS BACK SUPPORT lumbar supports




- Constructed of high strength, breathable nylon material
- Five self-contouring stays
- Rubber lining on stay pockets help prevent migration
- 1½" suspension straps adjust easily and lay flat
- Hook and loop closure with vinyl covered side pulls for maximum support
- 9½" posterior height

SIZE	WAIST CIRCUMFERENCE	PRODUCT #
X-Small	28" - 32"	BS1000-10
Small	32" - 36"	BS1000-12
Medium	36" - 40"	BS1000-14
Large	40" - 44"	BS1000-16
X-Large	44" - 48"	BS1000-18
XX-Large	48" - 52"	BS1000-20

### Indications:

Mild low back pain; Lumbar muscle weakness; Lumbar sprains and strains

# UPPER EXTREMITY | Torso



## PREMIUM SIZED ABDOMINAL BINDER torso



- Constructed of knit elastic with soft lining for extended wear and comfort
- Three (9") or four (12") individual panel design for best patient fit
- Unique stitching prevents individual panels from rolling
- Hook and loop closure
- Not made with natural rubber latex

SIZE	MEASUREMENT*	9"	12"
Small/Medium	30" - 45"	13991000	13981000
Medium/Large	46" - 62"	13992000	13982000
Large	63" - 74"	13993000	13983000
X-Large	75" - 84"	13994000	13984000
XX-Large	85" - 94"	13995000	13985000

\*Measure circumference of male waist or female hips

### Indications:

Abdominal muscle strain; Soreness in abdominal region;  
Post-abdominal surgery

## PREMIUM UNIVERSAL ABDOMINAL BINDER torso



- Constructed of knit elastic
- Choose between three (9") or four (12") individual or solid (10") panel design
- Unique stitching design prevents individual panels from rolling
- Hook and loop closure
- Universal sizing reduces inventory needs
- Available in five panel (15") ilio bypass version
- Not made with natural rubber latex

Size	Measurement	9"	10"	12"	15"
Universal	27" - 48"	13970000	13950000	13960000	—
X-Long	up to 54"	—	—	—	13871000
XX-long	up to 70"	—	—	—	13872900

\*Measure circumference of male waist or female hips

### Indications:

Abdominal muscle strain; Soreness in abdominal region;  
Post-abdominal surgery

## SIZED ABDOMINAL BINDER torso



- Not made with natural rubber latex
- Three (9") or four (12") panel design
- Hook and loop closure

SIZE	MEASUREMENT*	9"	12"
Small/Medium	30" - 45"	13661056	13651056
Medium/Large	46" - 62"	13662067	13652067
Large	63" - 74"	13663008	13653008
X-Large	75" - 84"	13664009	13654009
XX-Large	85" - 94"	13665010	13655010

\*Measure circumference of male waist or female hips

### Indications:

Abdominal muscle strain; Soreness in abdominal region;  
Post-abdominal surgery

## UNIVERSAL ABDOMINAL BINDER torso



- Choose between three (9") or four (12") individual or solid (10") panel design
- Hook and loop closure
- Universal sizing reduces inventory needs
- Not made with natural rubber latex

SIZE	MEASUREMENT	9"	10"	12"
Universal	27" - 48"	13660000	13640000	13650000
X-Long	up to 54"	—	13600000	—
XX-Long	up to 70"	—	13590000	13630000

\*Measure circumference of male waist or female hips

### Indications:

Abdominal muscle strain; Soreness in abdominal region; Post-abdominal surgery

## ECONOMY ABDOMINAL BINDER torso



- Three (9"), four (12") or solid (10") panel design
- Hook and loop closure
- Not made with natural rubber latex

Size	Measurement*	9"	10"	12"
Universal	27"-48"	—	A133000	—
Small/Medium	30"-45"	A131056	—	A132056
Medium/Large	46"-62"	A131067	—	A132067

\*Measure circumference of male waist or female hips

### Indications:

Abdominal muscle strain; Soreness in abdominal region; Post-abdominal surgery

## PREMIUM 6" RIB BELT torso



- Constructed of knit elastic with soft inner liner
- Two panel design with unique stitching prevents material from rolling
- Hook and loop closure for ease of application
- Male and female (contoured) styles for best fit
- Available in universal and sized options
- Not made with natural rubber latex

Size	Dimensions	Male	Female
Universal	6" x 40"	13780300	13780400
XX-Large	6" x 50"	13790309	13790409

### Indications:

Fractured or bruised ribs; Intercostal strains; Costovertebral sprains

## UNIVERSAL 6" RIB BELT torso




- Constructed of knit elastic
- Solid panel design
- Hook and loop closure for ease of application
- Available in universal, gender specific styles for best fit
- Not made with natural rubber latex

SIZE	DIMENSIONS	QTY	MALE	FEMALE
Universal	6" x 40"	12/Case	13740300	13740400
Universal X-Long	6" x 43"	12/Case	13749300	13749400

### Indications:

Fractured or bruised ribs; Intercostal strains; Costovertebral sprains

# UPPER EXTREMITY | Clavicle /Shoulder



## PREMIUM CLAVICLE STRAP clavicle



- Synthetic stockinette over foam padding does not absorb water and can be worn while bathing
- Serpentine stitch design to prevent migration
- Available in buckle or hook and loop closure

SIZE	BUCKLE CLOSURE	HOOK & LOOP CLOSURE
Infant	3001-01	3002-01
Child	3001-02	3002-02
Small	3001-06	3002-06
Medium	3001-03	3002-03
Large	3001-04	3002-04
X-Large	3001-05	3002-05

### Indications:

Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations

## HEAVY-DUTY CLAVICLE STRAP clavicle



- Constructed of foam padding with stockinette covering
- Extra long straps for easy adjustment
- Available in buckle or hook and loop closure

SIZE	BUCKLE CLOSURE	HOOK & LOOP CLOSURE
Infant	3004-01	3004-11
Child	3004-02	3004-12
Small	3004-06	3004-16
Medium	3004-03	3004-13
Large	3004-04	3004-14
X-Large	3004-05	3004-15

### Indications:

Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations

## 3-WAY CLAVICLE STRAP clavicle



- Foam construction with cotton stockinette cover
- 3-way posterior vector helps prevent migration
- Available in buckle or hook and loop closure

SIZE	BUCKLE CLOSURE	HOOK & LOOP CLOSURE
Child	3015-01	3015-11
Small	3015-02	3015-12
Medium	3015-03	3015-13
Large	3015-04	3015-14
X-Large	3015-05	3015-15

### Indications:

Clavicle fractures; Grade II and III acromioclavicular separations; Osteolysis of the clavicle; AC degenerative joint pain; Mechanical cervical disease; Postural dysfunction; Sternoclavicular dislocations



## SHOULDER IMMOBILIZER W/ CANVAS SWATHE shoulder



- Sling with attached swathe constructed of navy canvas material
- Shoulder strap includes O-ring/hook and loop closure
- Shoulder pad for added comfort
- Thumb loop helps prevent wrist drop
- Universal right/left

SIZE	POUCH SIZE	PRODUCT #
X-Small	6" x 12½"	9001-01
Small	7" x 13½"	9001-02
Medium	8" x 16½"	9001-03
Large	9" x 19½"	9001-04
X-Large	9" x 22½"	9001-05

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## CANVAS SHOULDER IMMOBILIZER W/ WAIST STRAP shoulder



- Constructed of navy canvas sling with web shoulder and waist straps
- Shoulder strap includes O-ring/hook and loop closure
- Adjustable buckle on waist strap
- Thumb loop helps prevent wrist drop
- Universal right/left

SIZE	POUCH SIZE	PRODUCT #
X-Small	6" x 12½"	9002-01
Small	7" x 13½"	9002-02
Medium	8" x 16½"	9002-03
Large	9" x 19½"	9002-04
X-Large	9" x 22½"	9002-05

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## MUSLIN SHOULDER IMMOBILIZER W/ WAIST STRAP shoulder



- Constructed of lightweight muslin sling with web shoulder and waist straps
- Adjustable buckle closures
- Closed pouch design
- Thumb loop helps prevent wrist drop
- Universal right/left

SIZE	POUCH SIZE	PRODUCT #
X-Small	6" x 12½"	9003-01
Small	7" x 13½"	9003-02
Medium	8" x 16½"	9003-03
Large	9" x 19½"	9003-04
X-Large	9" x 22½"	9003-05

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## SHOULDER IMMOBILIZER W/ FOAM STRAPS shoulder



- Constructed of canvas or lightweight tietex material with foam shoulder and waist straps
- Adjustable straps feature O-ring/hook and loop closure
- Thumb loop helps prevent wrist drop
- Universal right/left

SIZE	POUCH SIZE	CANVAS	TIETEX
X-Small	6" x 12½"	9023-00	A111012
Small	7" x 13½"	9023-01	A111015
Medium	8" x 16½"	9023-02	A111016
Large	9" x 19½"	9023-03	A111017
X-Large	9" x 22½"	9023-04	A111018

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## UNIVERSAL FOAM SHOULDER IMMOBILIZER shoulder



- Wrist cuff, humeral cuff and shoulder strap are adjustable and can be placed anywhere on belt to achieve desired fit
- Provides shoulder and upper arm immobilization following dislocation or post-surgical application
- Universal right/left

SIZE	BELT SIZE	WRIST CUFF SIZE	PRODUCT #
Universal	6 ½" x 52"	3 ½" x 13"	A110001

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## UNIVERSAL SLING STYLE SHOULDER IMMOBILIZER shoulder



- Perforated foam construction ensures patient comfort
- Sling style wrist cuff for optimal comfort and support
- Wrist cuff, humeral cuff and padded shoulder strap are adjustable and can be placed anywhere on belt to achieve desired fit
- Provides shoulder and upper arm immobilization following dislocation or post-surgical application
- Universal right/left

SIZE	BELT SIZE	WRIST CUFF SIZE	PRODUCT #
Adult	5" x 51"	5½" x 15"	9012-00
Pediatric	4¾" x 34½"	2" x 9¾"	9011-00

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## DELUXE SLING & SWATHE shoulder



- Lightweight perforated foam material provides patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Open or closed elbow designs
- Swathe attaches at any point on sling
- Available in pediatric and extra long versions
- Universal right/left

SIZE	SWATHE SIZE	POUCH SIZE	PRODUCT #
Universal	5" x 54"	13" x 19"	9006-00
Universal w/ X-Long Swathe	5" x 58"	13" x 19"	9006-50
Universal Enclosed Elbow	5" x 54"	17" x 19"	9014-00
Universal Pediatric	4½" x 44"	9" x 16½"	9007-00
Swathe Only	5" x 54"	—	9006-10

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## FOAM SLING & SWATHE shoulder



- Constructed of perforated foam material
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Available in universal and sized options
- Universal right/left

SIZE	SWATHE SIZE	POUCH SIZE	PRODUCT #
Universal	5" x 54"	13" x 19"	TX4442-00
Small	5" x 46"	10" x 16½"	TX4442-02
Medium	5" x 48"	11" x 18½"	TX4442-03
Large	5" x 51"	12" x 19"	TX4442-04

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## UNIVERSAL FOAM SLING & SWATHE shoulder



- Constructed of perforated foam material for patient comfort
- Adjustable shoulder straps attach to hook receivable sling to ensure accurate fit
- Swathe attaches at any point on sling
- Extra-long swathe size available
- Universal right/left

SIZE	SWATHE SIZE	POUCH SIZE	PRODUCT #
Universal	5" x 54"	13" x 19"	4442-00
Universal w/ X-Long Swathe	5" x 58"	13" x 19"	4442-10

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## VELPEAU SHOULDER IMMOBILIZER shoulder



- Constructed of heavy duty elastic material
- Durable foam padded wrist and humeral cuffs
- Contoured design with hook and loop closure for patient comfort
- Universal right/left

SIZE	BELT SIZE	PRODUCT #
X-Small	6" x 29"	9008-00
Small	6" x 32"	9008-01
Medium	6" x 35"	9008-02
Large	6" x 38"	9008-03
X-Large	6" x 41"	9008-04
XX-Large	6" x 44"	9008-05
XXX-Large	6" x 47"	9008-06

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## ELASTIC SHOULDER IMMOBILIZER shoulder



- Constructed of lightweight elastic material
- Foam wrist and humeral cuffs
- Low-profile design is contoured for patient comfort
- Hook and loop closure
- Universal right/left

SIZE	BELT SIZE	PRODUCT #
X-Small	6" x 29"	1000159
Small	6" x 32"	1010159
Medium	6" x 35"	1020159
Large	6" x 38"	1030159
X-Large	6" x 41"	1040159
XX-Large	6" x 44"	1050159

### Indications:

S/P glenohumeral acromioplasty; Surgical repairs of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint

## ECONOMY ELASTIC SHOULDER IMMOBILIZER shoulder



- Constructed of lightweight elastic
- Soft foam wrist and humeral cuffs
- Low profile design is contoured for patient comfort
- Hook and loop closure
- Universal right/left

SIZE	BELT SIZE	PRODUCT #
X-Small	6" x 29"	A111004
Small	6" x 32"	A111005
Medium	6" x 35"	A111006
Large	6" x 38"	A111007
X-Large	6" x 41"	A111008
XX-Large	6" x 44"	A111009

### Indications:

S/P glenohumeral acromioplasty; Surgical repair of the rotator cuff; Anterior glenohumeral capsular reconstruction; Anterior dislocation of the G-H joint



## SHOULDER P.A.D. abductors



- Constructed of lightweight canvas material
- Foam block provides 35° of abduction
- Shoulder and waist straps can be trimmed to desired length
- Thumb loop helps prevent wrist drop
- Allows arm to be placed in an optimal neutral position for effective healing
- Exercise ball packaged with product
- Universal right/left

SIZE	POUCH SIZE	PRODUCT #
Small	5½" x 14"	11640005
Medium	6½" x 17"	11640006
Large	7" x 19"	11640007
X-Large	7½" x 20½"	11640008

### Indications:

Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

## SHOULDER ABDUCTION POSITIONER abductors



- Constructed of hook receivable foam material for infinite adjustability
- Foam block provides 55° of abduction
- Adjustable wrist and humeral cuffs for proper positioning
- Foam shoulder pad provides additional patient comfort
- Designed for abduction and immobilization of the shoulder
- Universal right/left

SIZE	PILLOW DIMENSIONS	PRODUCT #
Small	25¼" x 23¾"	11650005
Medium	29½" x 26¾"	11650006
Large	30½" x 36"	11650007

### Indications:

Acute shoulder instability; Surgical rotator cuff repair; Anterior glenohumeral capsular reconstruction; Acute rotator cuff tendinopathy; Dislocation/subluxation of G-H joint; Acromioclavicular joint dislocations; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

## WARRIOR™ WING abductors



- Universal Sling Design = one size fits most
- Aluminum waist frame is easily formed for accurate positioning & comfort
- Airplane hinge for optimal abduction and shoulder recovery
- Allows both abduction, rotation, and internal or external fixation
- Includes three triangle pillows for increased abduction and/or fixation (90°, 45°, 30°, 15°)
- Easy enough to be self-applied


SIZE	DESCRIPTION	PRODUCT #
Universal	Warrior™ Wing, Left	11631000
Universal	Warrior™ Wing, Right	11632000

### Indications:

For non-surgical or post-op shoulder conditions that require diagnosis-specific arm positioning; Rotator cuff repairs; Bankart lesions; SLAP lesions; Glenohumeral dislocations / subluxation; Posterior / Anterior capsule repairs; Shoulder instabilities; Muscle and tendon repair; Joint reconstruction



# UPPER EXTREMITY | Arm/Elbow



## SPECIALTY ARM SLING arm sling/elevators



- Durable navy canvas construction with deep pocket design
- Infant and child sizes feature cartoon print
- Choice of foam or vinyl shoulder pad for patient comfort
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with O-rings

SIZE	POUCH SIZE	STRAP LENGTH	FOAM PAD	VINYL PAD
Infant	5" x 9"	1" x 21"	8004-01	8004-11
Child	6" x 11"	1" x 27"	8004-02	8004-12
Small	7" x 13"	2" x 34"	8004-03	8004-13
Medium	8" x 17"	2" x 37"	8004-04	8004-14
Large	9" x 20"	2" x 37"	8004-05	8004-15
X-Large	9" x 23"	2" x 42"	8004-06	8004-16

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## ENVELOPE ARM SLING arm sling/elevators



- Constructed of lightweight navy canvas material
- Deep pocket pouch design
- Foam shoulder pad for added comfort
- Thumb loop helps prevent wrist drop
- Buckle closure

SIZE	POUCH SIZE	STRAP LENGTH	PRODUCT #
Child	6" x 11"	1" x 27"	8001-01
Small	7" x 13"	1" x 34"	8001-02
Medium	8" x 17"	1" x 37"	8001-03
Large	9" x 20"	1" x 42"	8001-04
X-Large	9" x 23"	1" x 45"	8001-05

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## LIGHT BLUE ARM SLING arm sling/elevators



- Constructed of light blue cotton material with deep pocket design
- Foam shoulder pad provides additional comfort
- Thumb loop helps prevent wrist drop
- Hook and loop closure

SIZE	POUCH SIZE	STRAP LENGTH	PRODUCT #
Infant	5" x 9"	1" x 21"	10101
Child	6" x 11"	1" x 27"	10102
Small	7" x 13"	2" x 34"	10103
Medium	8" x 17"	2" x 37"	10104
Large	9" x 20"	2" x 37"	10105
X-Large	9" x 23"	2" x 42"	10106

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations



## NARROW POUCH ARM SLING arm sling/elevators



- Lightweight cotton material with narrow, closed pouch design
- Buckle closure
- Available in each or 12/Case quantity
- Available with or without shoulder pad

SIZE	POUCH SIZE	STRAP LENGTH	EACH	CASE	W/OUT PAD
Child	5" x 10"	1" x 28"	—	8017-00	A119901
Small	5" x 13"	1" x 34"	8017-51	8017-01	A119902
Medium	6" x 15"	1" x 38"	8017-52	8017-02	A119903
Large	6" x 18"	1" x 42"	8017-53	8017-03	A119904
X-Large	7" x 20"	1" x 46"	8017-54	8017-04	A119905

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## DEEP POCKET ARM SLING arm sling/elevators



- Deep pocket, open pouch design constructed of breathable tietex material
- Foam shoulder pad provides additional comfort
- Thumb loop helps prevent wrist drop
- Hook and loop closure with O-rings and slide buckle adjustment
- Available in each or 100/Case quantity

SIZE	POUCH SIZE	STRAP LENGTH	EACH	CASE
Infant	5" x 9"	1" x 21"	TX9902-01	—
Child	6" x 11"	1" x 27"	TX9902-02	TX9902-02/100
Small	7" x 13"	2" x 37"	TX9902-03	TX9902-03/100
Medium	8" x 17"	2" x 37"	TX9902-04	TX9902-04/100
Large	9" x 20"	2" x 37"	TX9902-05	TX9902-05/100
X-Large	9" x 22"	2" x 42"	TX9902-06	TX9902-06/100

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## ARM SLING W/ FOAM STRAP arm sling/elevators



- Constructed of canvas or lightweight tietex material
- Deep pocket, closed pouch design
- Foam shoulder strap provides maximum comfort
- Hook closure fastens anywhere on foam strap for superior adjustment

SIZE	POUCH SIZE	STRAP LENGTH	CANVAS	TIETEX
Small	7" x 13"	2" x 37"	8023-01	TX8023-21
Medium	8" x 17"	2" x 37"	8023-02	TX8023-22
Large	9" x 20"	2" x 37"	8023-03	TX8023-23
X-Large	9" x 22"	2" x 40"	8023-04	TX8023-24
XX-Large	10" x 23"	2" x 40"	8023-05	—

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## UNIVERSAL ARM SLING arm sling/elevators



- Features easy fold-back universal design
- Foam shoulder pad provides additional comfort
- Hook and loop closure with O-rings

SIZE	POUCH SIZE	STRAP LENGTH	CANVAS	COTTON	TIETEX
Adult	8" x 19"	2" x 37"	8020-00	A112050	TX8020-00
Child	5" x 9"	1" x 31"	8020-01	—	—

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## UNIVERSAL ARM SLING W/ FOAM STRAP arm sling/elevators



- Foam shoulder strap provides maximum comfort
- Hook and loop closure with O-rings

SIZE	PRODUCT #
Universal	A112000

### Indications:

Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

## PERFORATED ARM ELEVATOR/SLING arm sling/elevators



- Constructed of hook receivable perforated foam
- Converts easily to sling with shoulder pad included
- Can be trimmed for best patient fit
- Hook and loop closures for easy application

SIZE	POUCH SIZE	PRODUCT #
Universal	20" x 20½"	8788-50

### Indications:

Post upper extremity surgery; Upper extremity edema

## CONVOLUTED FOAM ARM ELEVATOR/SLING arm sling/elevators



- Constructed of convoluted foam inner with hook receivable outer lining
- Converts easily to sling with shoulder pad included
- Slide buckle for adjustment of shoulder strap
- Hook and loop closures for easy application

SIZE	POUCH SIZE	PRODUCT #
Universal	10" x 16"	8788-03

### Indications:

Post upper extremity surgery; Upper extremity edema

## ARM ELEVATOR/SLING W/ ICE POCKETS arm sling/elevators



- Heavy duty canvas construction
- Includes two ice bags for cold therapy
- Converts easily to sling

SIZE	POUCH SIZE	PRODUCT #
Small	10" x 14"	8011-01
Medium	10" x 16"	8011-02
Large	12" x 17"	8011-03

### Indications:

Post upper extremity surgery; Upper extremity edema

## WARRIOR™ POST-OP ELBOW elbow bracing



- Universal sizing (left & right specific)
- Single hinge lateral application for easy donning and doffing
- Hinge is easily adjusted in 10-degree increments (Flexion: 0° to 120°, Extension: 0° to 90°)
- Length can be adjusted from 16" to 24"
- Shoulder strap included for extra support and immobilization

SIZE	DESCRIPTION	PRODUCT #
Universal	Warrior™ Elbow, Left	EB6000-01
Universal	Warrior™ Elbow, Right	EB6000-02

### Indications:

Post-operative rehabilitation

## SLIMLINE ELBOW BRACE elbow bracing



- Hinge features range of motion settings in 10-degree increments
- Contoured softgoods at biceps and forearm are trimmable to provide precise fit
- Universal size fits most patients

SIZE	PRODUCT #
Universal	EB2040-00

### Indications:

Post-operative rehabilitation

## WRIST & ELBOW SPLINT elbow supports



- Constructed of durable leatherette material
- Provides immobilization of elbow at 90°
- Aluminum frame provides protection to the elbow

SIZE	BICEPS	FOREARM	PRODUCT #
	CIRCUMFERENCE	CIRCUMFERENCE	
Small	11½"	9¾"	6001-01
Medium	12½"	10"	6001-02
Large	13½"	11"	6001-03

### Indications:

Olecranon fractures; Radial head fractures; Hairline fracture of the radius, humerus and ulna; Ligamentous reconstruction procedures of the elbow

## PREMIUM TENNIS ELBOW SUPPORT elbow supports



- Made of lightweight, ventilated material
- Support is low-profile increasing patient compliance
- Inner pad can be strategically placed for maximum pain relief
- D-ring closure for easy application
- Universal size fits most patients

SIZE	MEASUREMENT	PRODUCT #
Universal	2" x 16"	12870000

### Indications:

Medial and lateral epicondylitis

## PNEUGEL® TENNIS ELBOW SUPPORT elbow supports



- Constructed of soft foam material for patient comfort
- Concentration of pressure allows compression to be applied where needed while not constricting uninvolved muscles or blood flow
- Air/gel removable bladder can be placed in freezer to provide cold therapy
- Inflation bulb included to adjust bladder
- Universal size fits most patients

SIZE	MEASUREMENT	PRODUCT #
Universal	2" x 20"	12840000

### Indications:

Medial and lateral epicondylitis

## UNIVERSAL TENNIS ELBOW SUPPORT elbow supports



- Perforated foam material is breathable for patient comfort
- Counterforce strap with loop lock contact closure for adjustable compression
- Universal size fits most patients

SIZE	MEASUREMENT	PRODUCT #
Universal	3" x 13"	6005-00

### Indications:

Medial and lateral epicondylitis

## SIZED TENNIS ELBOW SUPPORT elbow supports



- Constructed of surgical webbing with foam lining
- Counterforce strap with loop lock contact closure for adjustable compression
- Sized to ensure a precise fit

SIZE	MEASUREMENT	PRODUCT #
Small	3" x 9"	6004-01
Medium	3" x 10"	6004-02
Large	3" x 11"	6004-03
X-Large	3" x 12"	6004-04

### Indications:

Medial and lateral epicondylitis

## UNIVERSAL NEOPRENE TENNIS ELBOW SUPPORT elbow supports



- 1/8" neoprene construction provides compression, warmth and support
- Adjustable felt buttress can be placed just below the epicondyle for pain relief
- Loop lock contact closure for adjustable compression and application
- Universal size fits most patients

SIZE	MEASUREMENT	PRODUCT #
Universal	3" x 16"	NE7730-70

### Indications:

Medial and lateral epicondylitis

## NEOPRENE TENNIS ELBOW SUPPORT W/ PAD elbow supports



- 1/8" neoprene construction provides compression, warmth and support
- Incorporates 1/4" pad to alleviate stress over the lateral or medial epicondyle
- Elastic counterforce strap and loop lock contact closure allow for maximum compression
- Sized to ensure a precise fit

SIZE	FOREARM CIRCUMFERENCE	PRODUCT #
Small	9" - 10"	NE7731-72
Medium	10" - 11"	NE7731-73
Large	11" - 12"	NE7731-74
X-Large	12" - 14"	NE7731-75

### Indications:

Medial and lateral epicondylitis



## NEOPRENE ELBOW SLEEVE elbow sleeves




- 1/8" neoprene construction provides compression, warmth and support
- Anatomically contours to the elbow
- Available with or without elbow strap

SIZE	FOREARM		w/o STRAP	w/ STRAP
	CIRCUMFERENCE			
Small	9" - 10"		NE7727-72	NE7728-72
Medium	10" - 11"		NE7727-73	NE7728-73
Large	11" - 12"		NE7727-74	NE7728-74
X-Large	12" - 14"		NE7727-75	NE7728-75
XX-Large	14" - 16"		NE7727-76	NE7728-76

### Indications:

Tendonitis; Osteoarthritis; Arthritis; Olecranon bursitis; (Medial and lateral epicondylitis w/ strap)

## ELASTIC ELBOW SUPPORT elbow sleeves



- Pull on cotton elastic design
- Provides snug fit to achieve maximum comfort and compression

SIZE	FOREARM		PRODUCT #
	CIRCUMFERENCE		
Small	7" - 10"		6003-01
Medium	10" - 13"		6003-02
Large	13" - 15"		6003-03
X-Large	15" - 17"		6003-04


### Indications:

Minor elbow sprains and strains; Edema; Osteoarthritis; Olecranon bursitis





# UPPER EXTREMITY | Wrist/Hand



## SUEDE LEATHERETTE WRIST & WRIST/FOREARM SPLINT wrist/hand



- Suede leatherette material helps reduce migration
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

SIZE	LENGTH	WRIST		LENGTH	WRIST/FOREARM	
		RIGHT	LEFT		RIGHT	LEFT
Pediatric	—	—	—	6"	12991202	12991102
Child	—	—	—	9"	12992202	12992102
Small	8"	12982205	12982105	11"	12993205	12993105
Medium	8"	12982206	12982106	11"	12993206	12993106
Large	8"	12982207	12982107	11"	12993207	12993107
X-Large	8"	12982208	12982108	11"	12993208	12993108

### Indications:

Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

## CANVAS WRIST & WRIST/FOREARM SPLINT wrist/hand



- Constructed of sturdy blue canvas material with foam tricot lining
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

SIZE	LENGTH	WRIST		LENGTH	WRIST/FOREARM	
		RIGHT	LEFT		RIGHT	LEFT
Pediatric	4"	1010115	1060115	6"	1010025	1060025
Small	8"	1020115	1070115	11"	1020025	1070025
Medium	8"	1030115	1080115	11"	1030025	1080025
Large	8"	1040115	1090115	11"	1040025	1090025
X-Large	8"	1050115	1100115	11"	1050025	1100025

### Indications:

Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

## TIETEX WRIST & WRIST/FOREARM SPLINT wrist/hand



- Tietex material construction is lightweight and comfortable for patient
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past the distal palmar crease allowing full finger range of motion

SIZE	LENGTH	WRIST		LENGTH	WRIST/FOREARM	
		RIGHT	LEFT		RIGHT	LEFT
X-Small	4"	A123204	A123104	6"	A124204	A124104
Small	8"	A123205	A123105	11"	A124205	A124105
Medium	8"	A123206	A123106	11"	A124206	A124106
Large	8"	A123207	A123107	11"	A124207	A124107
X-Large	8"	A123208	A123108	11"	A124208	A124108

### Indications:

Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

## UNIVERSAL FOAM WRIST & WRIST/FOREARM wrist/hand



- Perforated foam material allows skin to breathe
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Moveable dorsal stay makes product universal in size
- D-ring closures for easy application and removal
- Available without binding around edges for trimming

SIZE	LENGTH	w/ BINDING RIGHT	w/ BINDING LEFT
Universal	8"	5066-82	5066-83
Universal	10"	5066-80	5066-81
SIZE	LENGTH	w/o BINDING RIGHT	w/o BINDING LEFT
Universal	8½"	5005-03	5005-04
Universal	11½"	5005-01	5005-02

### Indications:

Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

## UNIVERSAL BLACK FOAM WRIST & WRIST/FOREARM SPLINT



- Comfortable black foam material with tricot lining
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Moveable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

SIZE	LENGTH	RIGHT	LEFT
Universal	8"	BF5066-82	BF5066-83
Universal	10"	BF5066-80	BF5066-81

### Indications:

Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

## UNIVERSAL TWO-PANEL WRIST & WRIST/FOREARM SPLINT



- Blue canvas construction with foam tricot lining
- Malleable palmar and dorsal stays limit dorsiflexion motion
- Loop lock closure at wrist provides optimal immobilization and compression
- Universal right/left and sized to fit majority of patients

SIZE	LENGTH	PRODUCT #
Universal	8"	1000275
Universal	11"	1000285

### Indications:

Sprains or strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

## UNIVERSAL CUTAWAY WRIST & WRIST/FOREARM SPLINT wrist/hand



- Perforated foam construction allows skin to breath
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Loop lock closure at wrist provides compression and support
- Foam can be trimmed for desired fit

SIZE	LENGTH	RIGHT	LEFT
Universal	8"	5040-03	5040-04
Universal	11"	5040-01	5040-02

### Indications:

Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

## BLACK WRIST SPLINT wrist/hand



- Constructed of comfortable black foam material
- Malleable palmar stay provides anatomical fit and rigid support
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	8" RIGHT	8" LEFT	10" RIGHT	10" LEFT
Pediatric	5073-11	5073-12	-	-
X-Small	5073-01	5073-06	5074-01	5074-06
Small	5073-02	5073-07	5074-02	5074-07
Medium	5073-03	5073-08	5074-03	5074-08
Large	5073-04	5073-09	5074-04	5074-09
X-Large	5073-05	5073-10	5074-05	5074-10

### Indications:

Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## PREMIUM WRIST SPLINT wrist/hand



- Breathable material with four way stretch provides patient comfort
- Removable/malleable palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	LENGTH	RIGHT	LEFT
X-Small	8"	351XSR	351XSL
Small	8"	351SR	351SL
Medium	8"	351MR	351ML
Large	8"	351LR	351LL
X-Large	8"	351XLR	351XLL
XX-Large	8"	351XXLR	351XXLL

### Indications:

Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## LACE UP SUEDE LEATHERETTE WRIST SPLINT<sup>wrist/hand</sup>



- Soft perforated suede leatherette material with nylon tricot lining
- Removable/malleable palmar stay
- Single pull lace closure for compression and easy adjustment
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	LENGTH	RIGHT	LEFT
X-Small	8"	5023-01	5023-06
Small	8"	5023-02	5023-07
Medium	8"	5023-03	5023-08
Large	8"	5023-04	5023-09
X-Large	8"	5023-05	5023-10

### Indications:

Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## PREMIUM UNIVERSAL WRIST SPLINT<sup>wrist/hand</sup>



- Constructed of comfortable foam material that helps keep patient's skin cool and dry
- Dynamic closure allows easy adjustment for an intimate fit
- Removable palmar and dorsal stays make product universal and allows for proper support
- Universal right/left fits most patients

SIZE	LENGTH	PRODUCT #
Universal	8"	5068-00

### Indications:

Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## UNIVERSAL WRIST SPLINT<sup>wrist/hand</sup>



- Constructed of comfortable foam material
- Malleable palmar and dorsal stays provides anatomical fit and rigid support
- Movable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

SIZE	LENGTH	RIGHT	LEFT
Universal	7"	5067-01	5067-02

### Indications:

Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## ELASTIC WRIST SPLINT wrist/hand



- Cotton elastic construction provides compression
- Removable/malleable palmar stay
- Hook and loop closure for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	6" RIGHT	6" LEFT	8" RIGHT	8" LEFT
Pediatric	5015-01	5015-06	—	—
Small	A121205	A121105	A122205	A122105
Medium	A121206	A121106	A122206	A122106
Large	A121207	A121107	A122207	A122107
X-Large	A121208	A121108	A122208	A122108
XX-Large	—	—	5016-15	5016-20

### Indications:

Mild sprains and strains of the wrist; Tendonitis, Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## UNIVERSAL AMBIDEXTROUS WRISTLET wrist/hand



- 1/8" neoprene material provides compression, warmth and support
- Wraparound design for universal sizing
- Hook and loop closure for easy application
- Universal size fits most patients

SIZE	PRODUCT #
Universal	NE7740-70

### Indications:

Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

## ELASTIC WRIST SUPPORT wrist/hand



- Cotton elastic material with thumb opening for precise fit
- Provides snug fit for maximum comfort and compression
- Universal right/left

Size	Product #
Small	5019-01
Medium	5019-02
Large	5019-03
X-Large	5019-04

### Indications:

Wrist sprains and strains

## PREMIUM WRIST & THUMB SPLINT wrist/thumb



- Breathable material with four way stretch provides patient comfort
- Removable/adjustable thumb and palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	LENGTH	RIGHT	LEFT
X-Small	8"	350XSR	350XSL
Small	8"	350SR	350SL
Medium	8"	350MR	350ML
Large	8"	350LR	350LL
X-Large	8"	350XLR	350XLL
XX-Large	8"	350XXLR	350XXLL

### Indications:

De Quervain's Syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

## BLACK FOAM WRIST & THUMB wrist/thumb



- Constructed of comfortable foam material
- Removable/adjustable thumb and palmar stay
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

SIZE	LENGTH	RIGHT	LEFT
X-Small	8"	A125204	A125104
Small	8"	A125205	A125105
Medium	8"	A125206	A125106
Large	8"	A125207	A125107
X-Large	8"	A125208	A125108
XX-Large	8"	A125209	A125109

### Indications:

De Quervain's Syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

## PREMIUM UNIVERSAL WRIST & THUMB SPLINT wrist/thumb



- Constructed of comfortable foam material that helps keep patient's skin cool and dry
- Dynamic closure allows easy adjustment for an intimate fit
- Removable palmar and dorsal stays make product universal and allows for proper support
- Adjustable thumb stay ensures desired thumb positioning
- Universal right/left and size fits most patients

Size	Length	Product #
Universal	8"	5069-00

### Indications:

De Quervain's Syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries



## UNIVERSAL AMBIDEXTROUS WRISTLET W/ ABDUCTED THUMB wrist/thumb



- 1/8" neoprene construction provides warmth, compression and support
- Wraparound design for universal sizing
- Designed to hold thumb in abducted position
- Hook and loop closure for easy application

SIZE	PRODUCT #
Universal	NE7738-70

### Indications:

Arthritis; Thumb sprains; Tendonitis

## UNIVERSAL LEATHERETTE THUMB SPLINT wrist/thumb



- Constructed of durable perforated leatherette material
- Can be used alone or in combination with wrist and forearm splint for thumb support
- Universal size fits most patients

SIZE	PRODUCT #
Universal	5026-00

### Indications:

Thumb intrinsic muscle injuries; Post-thumb fracture;  
Collateral ligament injuries of the thumb

## WRIST/THUMB & THUMB SPICA wrist/hand



- Lightweight splint constructed of perforated foam
- Thumb spica supports thumb while allowing IP joint to remain free
- Wrist & thumb spica supports the wrist and IP joint of the thumb
- Soft inner lining wicks away perspiration for comfort and hygiene
- Stretch loop straps designed for a precise anatomical fit allowing hand and finger function
- May be trimmed for a custom fit
- Thumb version available in white or black foam

### THUMB - WHITE FOAM

SIZE	MEASUREMENT*	RIGHT	LEFT
Small/Medium	2½" - 3½"	340S/MR	340S/ML
Medium/Large	3½" - 4¼"	340M/LR	340M/LL

### THUMB - BLACK FOAM

SIZE	MEASUREMENT*	RIGHT	LEFT
Small/Medium	2½" - 3½"	343S/MR	343S/ML
Medium/Large	3½" - 4¼"	343M/LR	343M/LL

### WRIST/THUMB - WHITE FOAM

SIZE	MEASUREMENT*	RIGHT	LEFT
Small/Medium	2½" - 3½"	345S/MR	345S/ML
Medium/Large	3½" - 4¼"	345M/LR	345M/LL

\*Measure width of hand at MP joints

### Indications:

De Quervain's Syndrome; CMC arthritis; Post-thumb fractures; Tendon injury/repair



## BASEBALL FINGER SPLINTS aluminum splints



- Lightweight aluminum can be easily adjusted for proper fit
- Foam padding for patient comfort
- Easily moldable to maintain desired position
- Packaged 12/Case

SIZE	LENGTH	PRODUCT #
Small	4 ¼"	12101
Medium	4 ¾"	12102
Large	5 ¾"	12103

### Indications:

Fractures of the distal phalanx

## TOAD FINGER SPLINTS aluminum splints



- Lightweight aluminum can be easily adjusted for proper fit
- Foam padding for patient comfort
- Malleable arms fold to hold finger in desired alignment
- Packaged 12/Case

SIZE	LENGTH	PRODUCT #
Small	2 ¼"	11101
Medium	2 ¾"	11102
Large	2 7/8"	11103

### Indications:

Distal and proximal interphalangeal joint positioning

## FOLD OVER FINGER SPLINTS aluminum splints




- Lightweight aluminum can be adjusted for proper fit
- Foam padding for patient comfort
- Packaged 12/Case

SIZE	LENGTH	PRODUCT #
Small	1"	12201
Medium	2"	12202
Large	3"	12203
X-Large	3 ½"	9111-07

### Indications:

Protects injured nails and fingertips; Baseball or transverse fractures of the finger

## FOUR PRONG FINGER SPLINTS aluminum splints



- Lightweight aluminum can be adjusted for proper fit
- Foam padding for patient comfort
- Packaged 12/Case

SIZE	LENGTH	PRODUCT #
Small	2"	12301
Medium	3"	12302
Large	3 1/2"	12303

### Indications:

Provides protection and immobilization to injured fingers

## GUTTER SPLINTS aluminum splints



- Lightweight aluminum can be easily adjusted for proper fit
- Foam padding for patient comfort
- Features slight angle at tip to address a variety of phalangeal injuries
- Packaged 12/Case

MEASUREMENT	PRODUCT #
1 1/2" x 7/8"	11201
3" x 7/8"	11202
4" x 7/8"	11205
5 1/2" x 7/8"	11203
7" x 7/8"	11204

### Indications:

Phalangeal and metacarpal fractures; Finger dislocations

## PADDED FINGER SPLINT ASSORTMENT KIT aluminum splints



- Includes an assortment of 57 splints packaged in convenient storage container
- Constructed of lightweight aluminum with foam padding

DESCRIPTION	QUANTITY	PRODUCT #
Complete Kit	1 Kit	9100-00
Kit Includes:		
2 3/4" Toad Splints	12 Each	
2 7/8" Toad Splints	12 Each	
1 1/2" Gutter Splints	6 Each	
3" Gutter Splints	6 Each	
2" Fold Over Finger Splints	6 Each	
3" Fold Over Finger Splints	6 Each	
3 1/2" Fold Over Finger Splints	9 Each	

(All sizes sold individually as 12/Case)

## FOAM PADDED FINGER STRIPS aluminum splints



- Constructed of aluminum with foam padding
- Easy to form and can be cut to desired length with scissors
- Packaged 12/Case

MEASUREMENT	PRODUCT #
½" x 9"	11001
½" x 18"	11002
¾" x 18"	11003
1" x 18"	11004
1 ½" x 18"	9115-05
2" x 18"	9115-06

### Indications:

Phalangeal fractures

## FOAM PADDED SPLINTING MATERIAL aluminum splints



- Constructed of aluminum with foam padding
- Can be trimmed with scissors to achieve desired length for comfortable immobilization
- Packaged 3/Case

MEASUREMENT	PRODUCT #
2" x 16"	3302D39
3" x 16"	3302D59
4" x 16"	3302D79

### Indications:

Fractures and dislocations

## FINGER SPLINTS aluminum splints



- Lightweight aluminum foam construction can be easily adjusted for proper fit
- Universal right/left
- Packaged 6/Case

SIZE	PRODUCT #
Small	9125-01
Medium	9125-02
Large	9125-03
X-Large	9125-04

### Indications:

Phalangeal fractures

## FINGER DRESSING/FINGER CAGE aluminum splints



- 25 yard roll of tubular bandage
- Perfect for all types of dressings

SIZE	WIDTH	PRODUCT #
Small Fingers & Toes	5/8"	FD1
Fingers & Toes, Infants Hands & Feet	1"	FD2
Small Adult Hands, Arms & Feet, Infant & Child Arms & Feet	1 1/2"	FD3

- Cage constructed of durable aluminum
- Assists in application of tubular finger dressing
- Reusable

SIZE	DIAMETER	Length	Product #
Small	1"	6"	9150-01
Large	1 3/4"	6"	9150-02

## PLASTALUME SPLINTS aluminum splints



- Includes an assortment of 54 splints packaged in convenient storage container
- Constructed of blue plastic attached to white foam padding

DESCRIPTION	QTY	PRODUCT #
Complete Kit	1 Kit	9120-00
Kit Includes:		
1 1/4" Length	6 Each	9119-01
2 1/4" Length	6 Each	9119-02
3 1/4" Length	6 Each	9119-03
4 1/4" Length	6 Each	9119-04
3 1/4" Length, circular base	6 Each	9119-05
4 1/4" Length, circular base	6 Each	9119-06
5 1/4" Length, circular base	6 Each	9119-07
6 1/4" Length, circular base	6 Each	9119-08
7 1/4" Length, circular base	6 Each	9119-09

(All sizes sold individually as 12/Case)

## STAX FINGER SPLINT KIT aluminum splints



- Includes an assortment of 30 splints packaged in convenient storage container
- Protects distal phalanx and extends the DIP joint
- Open tip permits examination of fingertip and nailbed

DESCRIPTION	QTY	PRODUCT #
Complete Kit	1 Kit	9122-00
Kit Includes:		
Size 1	5 Each	9121-01
Size 2	5 Each	9121-02
Size 3	5 Each	9121-03
Size 4	5 Each	9121-04
Size 5	5 Each	9121-05
Size 5 1/2	2 Each	9121-07
Size 6	2 Each	9121-06
Size 7	1 Each	9121-08

(All sizes sold individually)

### Indications:

Fingertip injuries; Nailbed injuries; Mallet finger

## WRIST & FOREARM SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Foam padding for patient comfort

SIZE	PRODUCT #	
	RIGHT	LEFT
Child	11301	11304
Youth	11302	11305
Adult	11303	11306

### Indications:

Post-operative rehabilitation; Severe sprains and strains of the wrist; Wrist and distal forearm immobilization

## COLLES' SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design for air circulation
- Foam padding for patient comfort

SIZE	w/ FOAM	
	RIGHT	LEFT
Infant	9105-07	9105-08
Child	11701	11704
Youth	11702	11705
Adult	11703	11706

### Indications:

Colles' injury; Wrist immobilization

## PLASTIC COLLES' SPLINT aluminum splints



- Pre-molded plastic design for an anatomical fit and proper support
- Perforated design for air circulation
- Foam padding for patient comfort

SIZE	RIGHT	LEFT
Child	9105-24	9105-20
Youth	9105-25	9105-21
Adult	9105-26	9105-22

### Indications:

Colles' injury; Wrist immobilization

## SHORT METACARPAL SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Contoured design to fit in palm of hand
- Perforated design for air circulation
- Foam padding for patient comfort

SIZE	RIGHT	LEFT
Child	9136-01	9136-02
Youth	9134-01	9134-02
Adult	9135-01	9135-02

### Indications:

Phalangeal and metacarpal fractures

## LONG METACARPAL SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Perforated design promotes air circulation
- Foam padding for patient comfort

SIZE	RIGHT	LEFT
Child	9148-01	9148-04
Youth	9148-02	9148-05
Adult	9148-03	9148-06

### Indications:

Post-operative rehabilitation; Severe sprains and strains of the wrist; Wrist and distal forearm immobilization

## MASON ALLEN HAND SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Available with or without foam padding
- Universal right/left

SIZE	w/ FOAM	w/o FOAM
Small	9133-01	9133-04
Universal	9133-03	9133-06

### Indications:

Immobilization of entire hand and wrist



## ELBOW SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Aluminum bar can be bent to adjust angle
- Perforated design for air circulation
- Foam padding for patient comfort
- Universal right/left

SIZE	PRODUCT #
Child	9118-01
Youth	9118-02
Adult	9118-03

### Indications:

Elbow immobilization

## POSTERIOR LEG SPLINT aluminum splints



- Pre-molded aluminum design for an anatomical fit and proper support
- Open front for easy application
- Short or full leg length options
- Foam padding for patient comfort

SHORT LEG	LENGTH	PRODUCT #
Child	12"	9127-01
Youth	14"	9127-02
Adult	16"	9127-03


FULL LEG	LENGTH	PRODUCT #
Child	19 ½"	9128-01
Youth	25"	9128-02
Adult	31"	9128-03

### Indications:

**Short:** Transverse fractures of the tibia/fibula; Ankle fractures; Lower leg immobilization; **Full:** Shaft fractures of the femur; Knee immobilization



# LOWER EXTREMITY | Knee/Leg



## DELUXE TRI-PANEL KNEE IMMOBILIZER knee immobilizers



- Constructed of foam body with durable blue canvas panels
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application
- Also available in black foam construction

LENGTH	CANVAS	FOAM
12"	1120917	BF1120917
14"	1140917	—
16"	1160917	BF1160917
19"	1190917	BF1190917
22"	1220917	BF1220917
24"	1240917	BF1240917

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## SUPERLITE TRI-PANEL KNEE IMMOBILIZER knee immobilizers



- Constructed of foam body with white perforated canvas panels
- T-bars stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application

LENGTH	PRODUCT #
12"	7091-50
16"	7091-56
19"	7091-51
22"	7091-55
24"	7091-52
26"	7091-57
28"	7091-58

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## FOAM TRI-PANEL KNEE IMMOBILIZER knee immobilizers



- Constructed of blue perforated foam for patient comfort
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Trimmable foam allows for a more customized fit

LENGTH	PRODUCT #
12"	7091-03
16"	7091-04
19"	7091-05
22"	7091-06
24"	7091-07
26"	7091-08

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## TIETEX TRI-PANEL KNEE IMMOBILIZER knee immobilizers



- Constructed of black foam body with durable tietex panels
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral panels allow for universal sizing
- Patella strap assists with application

Length	Product #
12"	A141000
16"	A142000
19"	A143000
22"	A144000
24"	A145000
26"	A146000

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## KNEE IMMOBILIZER W/ ELASTIC STRAPS knee immobilizers



- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Adjustable medial and lateral stays allow for universal sizing
- Elastic straps provide compression helping to reduce edema

LENGTH	BLUE	BLACK
12"	—	BF7072-01
16"	7072-02	BF7072-02
20"	7072-03	BF7072-03
22"	7072-04	BF7072-04
24"	7072-05	BF7072-05
26"	—	BF7072-06

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## KNEE IMMOBILIZER W/ ELASTIC STRAPS knee immobilizers



- Perforated foam construction for patient comfort and breathability
- Available with or without t-bar stays
- Adjustable medial and lateral stays allow for universal sizing
- Elastic straps provide compression helping to reduce edema

LENGTH	W/O T-BARS	W/ T-BARS	40" CIRCUMFERENCE
13"	7071-01	7071-11	—
16"	7071-02	7071-12	—
18"	—	—	7071-55
19"	7071-03	7071-13	—
23"	7071-05	7071-15	7071-57
28"	7071-07	7071-17	7071-59

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## CUTAWAY KNEE IMMOBILIZER W/ PATELLA STRAP knee immobilizers



- Perforated foam construction for patient comfort and breathability
- Available with or without t-bar stays
- Contoured posterior stays for correct anatomical fit
- Adjustable medial and lateral stays allow for universal sizing
- Elastic patella strap helps reduce knee edema
- Also available in black foam construction (w/o t-bars)

LENGTH	W/ T-BARS	W/O T-BARS
13"	—	9307-01
14"	7066-14	—
16"	7066-16	9307-07
18"	7066-18	9307-02
20"	7066-20	9307-03
22"	7066-22	9307-06
24"	7066-24	9307-04
26"	7066-26	—
28"	7066-28	9307-05

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## CUTAWAY KNEE IMMOBILIZER knee immobilizers



- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Available with straight or contoured posterior stays
- Adjustable medial and lateral stays allow for universal sizing
- Quick fit tabs assist with application

LENGTH	STRAIGHT	CONTOURED
12"	7041-01	7041-11
16"	7041-02	7041-22
20"	7041-03	7041-33
24"	7041-04	7041-44
28"	7041-05	7041-55

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## CUTAWAY KNEE IMMOBILIZER knee immobilizers



- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Available with straight or contoured posterior stays
- Adjustable medial and lateral stays allow for universal sizing
- Loop lock straps affix posterior to encircle entire leg
- Also available in black foam construction

LENGTH	T-BARS	CONTOURED	STRAIGHT	BLACK FOAM
12"	—	—	—	BF4433-12
16"	—	4433-10	—	BF4433-16
19"	4433-11	4433-01	4433-02	BF4433-19
22"	4433-22	4433-05	—	BF4433-22
24"	4433-33	4433-03	4433-04	BF4433-24
26"	—	—	—	BF4433-26

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures



## BLACK CUTAWAY KNEE IMMOBILIZER knee immobilizers



- Constructed of lightweight black foam
- T-bar stays provide optimal strength
- Adjustable medial and lateral stays allow for universal sizing

LENGTH	PRODUCT #
12"	A140001
16"	A147000
19"	A148000
22"	A149000
24"	A140000

\*Universal size fits thigh circumference 12" - 24"

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## SIZED CANVAS KNEE IMMOBILIZER knee immobilizers



- Constructed of durable blue canvas
- Available with or without t-bar stays
- Contoured posterior stays for correct anatomical fit
- Alternating web straps with loop lock closure

W/O T-BAR	MEASUREMENT*	12"	16"	19"	22"	24"
Small	14" - 16"	1010127	1010167	1010197	1010227	1010247
Medium	16" - 18"	1020127	1020167	1020197	1020227	1020247
Large	18" - 20"	1030127	1030167	1030197	1030227	1030247
X-Large	20" - 22"	1040127	1040167	1040197	1040227	1040247
XX-Large	22" - 24"	1050127	1050167	1050197	1050227	1050247

\*Measure thigh circumference 6" above mid-patella

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## SIZED FOAM KNEE IMMOBILIZER knee immobilizers



- Perforated foam construction for patient comfort and breathability
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide desired compression

SIZE	MEASUREMENT*	16"	20"	22"	24"
Small	14" - 16"	4461-02	4462-02	4464-02	4463-02
Medium	16" - 18"	4461-03	4462-03	4464-03	4463-03
Large	18" - 20"	4461-04	4462-04	4464-04	4463-04
X-Large	20" - 22"	4461-05	4462-05	4464-05	4463-05
XX-Large	22" - 24"	—	4462-06	—	4463-06

\*Measure thigh circumference 6" above mid-patella

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## SIZED BLACK FOAM KNEE IMMOBILIZER knee immobilizers



- Constructed of black foam with tricot liner
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide desired compression

SIZE	MEASUREMENT*	10"	12"	16"	20"
Small	14" - 16"	—	BF7006-01	BF7007-01	BF7008-01
Medium	16" - 18"	—	BF7006-02	BF7007-02	BF7008-02
Large	18" - 20"	—	BF7006-03	BF7007-03	BF7008-03
X-Large	20" - 22"	—	BF7006-04	BF7007-04	BF7008-04
XX-Large	22" - 24"	—	BF7006-05	BF7007-05	BF7008-05
XXX-Large	24" - 26"	—	BF7006-06	BF7007-06	BF7008-06
Pediatric Small	8" - 10"	14970005			
Pediatric Medium	10" - 12"	14970006			
Pediatric Large	12" - 14"	14970007			

\*Measure thigh circumference 6" above mid-patella

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## SIZED TIETEX KNEE IMMOBILIZER knee immobilizers



- Constructed of durable tietex foam
- T-bar stays provide optimal strength
- Contoured posterior stays for correct anatomical fit
- Alternating straps with loop lock closure provide desired compression

SIZE	MEASUREMENT*	16"	20"	22"	24"
Small	14" - 16"	A146005	A142005	A143005	A144005
Medium	16" - 18"	A146006	A142006	A143006	A144006
Large	18" - 20"	A146007	A142007	A143007	A144007
X-Large	20" - 22"	A146008	A142008	A143008	A144008
XX-Large	22" - 24"	A146009	A142009	A143009	A144009

\*Measure thigh circumference 6" above mid-patella

### Indications:

Mild to severe knee injuries; Post-operative immobilization; Grade II and III collateral ligament sprain; Patellar dislocations; Osteochondral fractures

## WARRIOR® RECOVERY KNEE BRACE post-op



- Easy-to-use, lever-action, micro-adjustable telescoping struts for length adjustments from 18" to 29"
- Four convenient preset drop lock settings, 0°, 15°, 30°, 45°
- Enhanced range-of-motion (ROM) hinge interface for easy setting and patient compliance
- ROM is set in 10° increments
- Soft, circumferential and trimmable straps plus removable Y-tab feature optimizes individual patient fit and comfort
- Inner pad panel adjustment (2mm hex wrench, not included) to avoid incision sites and provide a more customized fit
- Customizable anti-migration strap pads for superior fixture and reduced slippage
- Integrated "bend zone" on lower hinge arms allows for customization and accommodation to varying patient anatomy
- Laminate reinforced padding to support prolonged use and maintain pad integrity
- Extended padding at malleoli to prevent irritation and improve comfort during use

SIZE	DESCRIPTION	ADJUSTMENT LENGTH	PRODUCT #
Universal	Warrior® Recovery Post-Op Knee Brace	18" - 29"	KB9000-01

### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management.

## TRANSITION KNEE BRACE post-op



- Drop-lock design offers easy access for pin adjustment
- Range of motion can be set in 10° increments
- Superior and inferior score zones for length adjustment
- Dial can be locked to encourage patient compliance
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

DESCRIPTION	LENGTH ADJUSTMENT	PRODUCT #
Full Foam	23" - 17"	KB5000-00
Cool Foam	23" - 17"	KB5000-01
Full Foam, X-Long	26" - 20"	KB5000-02

### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

## SLIMLINE II KNEE BRACE post-op



- Snap-lock hinge cover design with simple pin adjustment
- Range of motion can be set at 5°, 15°, 30°, 60° or 90° of flexion and 0°, 15°, 30°, 60° or 90° of extension
- Superior and inferior score zones for length adjustment
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

DESCRIPTION	LENGTH ADJUSTMENT	PRODUCT #
Full Foam	23" - 17"	KB6000-00
Cool Foam	23" - 17"	KB6000-01

### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

## SLIMLINE KNEE BRACE post-op



- Posterior offset hinge for correct anatomical fit
- Hinge can be easily adjusted in 10° degree increments of flexion and extension without tools
- Paddles allow for easy length adjustment
- Unique calf strap helps prevent migration
- Dial can be locked to encourage patient compliance

DESCRIPTION	LENGTH ADJUSTMENT	PRODUCT #
Telescopic Paddles	27" - 22"	KB2040-00
Fixed Paddles	22"	KB2040-10
No Paddles	21"	KB2040-40

### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

## KEY KNEE BRACE post-op



- Posterior offset hinge for best anatomical fit
- Hinge can be adjusted in 10 degree increments of flexion and extension
- Length is adjustable with easy-to-use key
- Unique calf strap helps prevent migration

DESCRIPTION	LENGTH ADJUSTMENT	PRODUCT #
Telescopic Paddles	27" - 22"	KB3001-00

### Indications:

Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

## OPEN PATELLA KNEE SUPPORT knee supports



- 1/8" neoprene construction provides warmth, compression and support
- Available with oval pad for additional patella protection
- Open patella for patient comfort

SIZE	MEASUREMENT*	W/O PAD	W/ PAD
X-Small	13½" - 15½"	NE7702-71	NE7702-01
Small	15½" - 18"	NE7702-72	NE7702-02
Medium	18" - 20½"	NE7702-73	NE7702-03
Large	20½" - 23"	NE7702-74	NE7702-04
X-Large	23" - 25½"	NE7702-75	NE7702-05
XX-Large	25½" - 28"	NE7702-76	NE7702-06
XXX-Large	28" - 30½"	NE7702-77	NE7702-07
XXXX-Large	30½" - 33"	—	NE7702-08

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

## UNIVERSAL KNEE SUPPORT knee supports



- 1/8" neoprene construction provides warmth, compression and support
- Anterior wraparound strapping above and below knee helps prevent migration
- Universal adult and pediatric sizes available

SIZE	MEASUREMENT*	PRODUCT #
Universal	up to 25"	NE7710-70
Pediatric	up to 12"	NE7710-71

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Quadriceps tendonitis; Plica syndrome

## FAN KNEE STABILIZER knee supports



- 1/8" neoprene construction provides warmth, compression and support
- Features four criss-cross elastic straps for additional patella support
- Strapping above and below the knee helps prevent migration
- Medial/lateral stays provide stability
- Enclosed tubular buttress and adjustable felt buttress included

SIZE	MEASUREMENT*	PRODUCT #
Small	15½" - 18"	NE7707-72
Medium	18" - 20½"	NE7707-73
Large	20½" - 23"	NE7707-74
X-Large	23" - 25½"	NE7707-75
XX-Large	25½" - 28"	NE7707-76
XXX-Large	28" - 30½"	NE7707-77

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Patellar dislocation; Quadriceps tendonitis; Plica syndrome

## KNEE SUPPORT W/TRIMMABLE BUTTRESS knee supports



- 1/8" neoprene construction provides warmth, compression and support
- Anterior strapping above and below knee helps prevent migration while medial/lateral stays provide additional stability
- Enclosed tubular buttress can be trimmed for desired patella position

SIZE	MEASUREMENT*	PRODUCT #
X-Small	13½" - 15½"	NE7717-71
Small	15½" - 18"	NE7717-72
Medium	18" - 20½"	NE7717-73
Large	20½" - 23"	NE7717-74
X-Large	23" - 25½"	NE7717-75
XX-Large	25½" - 28"	NE7717-76
XXX-Large	28" - 30½"	NE7717-77

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Quadriceps tendonitis; Plica syndrome

## DELUXE KNEE SUPPORT W/TRIMMABLE BUTTRESS knee supports



- 1/8" neoprene construction provides warmth, compression and support
- Anterior strapping above and below knee helps prevent migration while medial/lateral stays provide additional stability
- Enclosed tubular buttress can be trimmed for desired patella position

SIZE	MEASUREMENT*	PRODUCT #
X-Small	13½" - 15½"	NE7718-71
Small	15½" - 18"	NE7718-72
Medium	18" - 20½"	NE7718-73
Large	20½" - 23"	NE7718-74
X-Large	23" - 25½"	NE7718-75
XX-Large	25½" - 28"	NE7718-76
XXX-Large	28" - 30½"	NE7718-77

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Patellar dislocation; Quadriceps tendonitis; Plica syndrome

## ELASTIC KNEE SUPPORT knee supports



- Pull up cotton elastic design
- Provides snug fit for maximum compression
- Available in open and closed patella styles

SIZE	MEASUREMENT*	CLOSED PATELLA	OPEN PATELLA
Small	12" - 15"	7042-01	7043-01
Medium	15" - 18"	7042-02	7043-02
Large	18" - 21"	7042-03	7043-03
X-Large	21" - 24"	7042-04	7043-04
XX-Large	24" - 27"	7042-05	7043-05

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

## VISCO ELASTIC KNEE SUPPORT knee supports



- Visco elastic construction with four way stretch for maximum compression
- Silicone buttress for patella stabilization
- Spiral stays provide moderate medial/lateral stability
- Closed patella design

SIZE	MEASUREMENT*	PRODUCT #
X-Small	13 ½" – 15 ½"	14801004
Small	15 ½" – 18"	14801005
Medium	18" – 20 ½"	14801006
Large	20 ½" – 23"	14801007
X-Large	23" – 25 ½"	14801008
XX-Large	25 ½" – 28"	14801009
XXX-Large	28" – 30 ½"	14801010

\*Measure circumference 6" above mid-patella

### Indications:

Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

## ARTHROSCOPIC FOAM KNEE WRAP knee supports



- Close fitting foam wrap provides compression to help control edema
- Optional pressure-sensitive adhesive foam pads can be placed in desired position after lateral patellar retinacular release or as adjunct for patellar stabilization
- Hook tab closures attach anywhere on wrap


SIZE	PRODUCT #
Universal	FF4061-00

### Indications:

Post-arthroscopic or arthroplasty knee surgery; Moderate to severe joint effusion



## TOTAL UNIVERSAL BUCKS TRACTION BOOT bucks traction




- Constructed of convoluted foam and laminated hook receivable material
- Medial and lateral stays provide stability and are removable for ease of application
- Full-foot design helps prevent foot drop
- Trimmable foam body and loop lock closure for more customized fit
- Spreader bar attached

SIZE	PRODUCT #
Universal	8140-00

### Indications:

Hip fractures; Hip dislocations; Femur fractures

## FULL FOOT BUCKS TRACTION BOOT bucks traction



- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- Full foot design helps prevent foot drop
- 1" ventilation holes increase airflow
- Trimmable foam and loop lock closures
- Spreader bar attached

SIZE	PRODUCT #
Universal	8130-00
Pediatric	8106-01
Small	8106-03
Medium	8106-04
Large	8106-05
X-Large	8106-06

### Indications:

Hip fractures; Hip dislocations; Femur fractures

## UNIVERSAL BUCKS TRACTION BOOT bucks traction



- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable foam body for more customized fit
- Universal size fits most patients
- Spreader bar attached

SIZE	PRODUCT #
Universal	8121-00

### Indications:

Hip fractures; Hip dislocations; Femur fractures

## SIZED BUCKS TRACTION BOOT bucks traction




- Non-convoluted foam construction helps reduce pressure points
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable foam and loop lock closures
- Spreader bar attached

SIZE	PRODUCT #
Universal	8120-00
Pediatric	8105-01
X-Small	8105-02
Small	8105-03
Medium	8105-04
Large	8105-05
X-Large	8105-06

### Indications:

Hip fractures; Hip dislocations; Femur fractures

## BUCKS TRACTION BOOT W/TRIANGLE BUCKLE bucks traction




- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- 1" ventilation holes increase airflow
- Trimmable foam and hook closures
- Traction cord attaches to triangle buckle

SIZE	PRODUCT #
Pediatric	19401
X-Small	19402
Small	19403
Medium	19404
Large	19405
X-Large	19406

### Indications:

Hip fractures; Hip dislocations; Femur fractures

## ECONOMY BUCKS TRACTION BOOT bucks traction



- Convoluted foam construction helps prevent migration
- Medial and lateral stays provide additional stability
- Trimmable foam and loop lock closures
- Spreader bar attached


SIZE	PRODUCT #
Universal	A150000
Small	A155005
Medium	A155006
Large	A155007
X-Large	A155008

### Indications:

Hip fractures; Hip dislocations; Femur fractures



# LOWER EXTREMITY | Ankle/Foot & Heel



## CONFOR® ANKLE STIRRUP ankle



- Confor® open cell polyurethane foam provides a superior fit
- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Elasticized stretch loop straps provide compression to help reduce edema
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left
- Available in each or 6/Case quantity

SIZE	HEIGHT	EACH	CASE
Standard	10"	AB2340-00	AB2340-06
Small	8½"	AB2341-00	AB2341-06

### Indications:

Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

## AIR ANKLE STIRRUP ankle



- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air bladder allows for simple adjustment of air with inflation bulb
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left
- Kit features: Standard air ankle stirrup, ice pack, elastic bandage and exercise band

SIZE	HEIGHT	PRODUCT #
Standard	10"	AB2360-00
Small	8½"	AB2361-00
Pediatric	6"	AB2362-00
Kit		SC1000-00

### Indications:

Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

## AIR/GEL ANKLE STIRRUP ankle



- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air/gel bladder can be placed in freezer for cold therapy
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits easily in shoe
- Universal right/left

SIZE	HEIGHT	PRODUCT #
Standard	10"	A152000
Small	8.5"	A152005

### Indications:

Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

## DEROYAL® FUNCTIONAL ANKLE BRACE ankle



- Anatomically designed uprights for a secure fit
- Hinged stirrup maintains support throughout range of motion
- Padded foam inlay for high energy absorption during activity
- Polyurethane foam padding for patient comfort and shock absorption
- Stretch loop straps for maximum support and edema control
- Universal right/left

SIZE	MALE	FEMALE	PRODUCT #
	SHOE SIZE	SHOE SIZE	
X-Small	—	5 - 7	AB2000-10
Small	5 - 8	7½ - 10	AB2000-12
Medium	8½ - 12	10½ - 13	AB2000-14
Large	12½ - 16	13½ +	AB2000-16

### Indications:

Acute and chronic ankle sprains; Chronic instability; Prophylactic talocrural support

## DEROYAL® SPORTS ORTHOSIS ankle



- Constructed of durable nylon for long lasting wear
- Figure 8 strapping mimics athletic taping
- Lace-up design for a secure fit
- Elastic cuff closure for added support and compression
- Low Profile design allows use in athletic shoe
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
X-Small	10"-11"	AB2900-10
Small	11"-12"	AB2900-12
Medium	12"-13"	AB2900-14
Large	13"-14"	AB2900-16
X-Large	14"-15"	AB2900-18
XX-Large	15"-16"	AB2900-20

\*Measure circumference from the base of heel around the ankle

### Indications:

Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

## DEROYAL® SPORTS ANKLE BRACE 2 ankle



- Constructed of durable nylon for long lasting wear
- Universal stabilization strap forces ankle into inversion or eversion
- Cuff closure for added support and compression
- Lace-up design for secure fit and primary stabilization
- Low profile design allows use in athletic shoe
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
X-Small	10" - 11"	15600004
Small	11" - 12"	15600005
Medium	12" - 13"	15600006
Large	13" - 14"	15600007
X-Large	14" - 15"	15600008
XX-Large	15" - 16"	15600009

\*Measure circumference from the base of heel around the ankle

### Indications:

Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

## CANVAS LACE-UP ANKLE BRACE ankle



- Constructed of durable canvas with flannel lining
- Contoured design for best anatomical fit
- Removable medial and lateral metal stays provide stability
- Lace-up design for a secure fit
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
Small	7" - 8"	8761-01
Medium	8" - 10"	8761-02
Large	10" - 12"	8761-03
X-Large	12" - 14"	8761-04

\*Measure circumference from the base of heel around the ankle

### Indications:

Acute and chronic ankle sprains; Ankle rehabilitation

## FIGURE 8 WRAP ANKLE SUPPORT ankle



- Constructed of a durable, breathable mesh sock material
- Elastic strap provides compression and edema control
- Figure 8 wrap with hook and loop closure provides additional support
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
Small	7" - 8"	4016-01
Medium	8" - 10"	4016-02
Large	10" - 12"	4016-03
X-Large	12" - 14"	4016-04

\*Measure circumference from the base of heel around the ankle

### Indications:

Mild ankle sprains and strains; Ankle edema; Peroneal and anterior tibialis tendonitis; Osteoarthritis of the talocrural and subtalar joints

## ELASTIC ANKLE SLEEVE ankle



- Pull on cotton elastic design
- Provides snug fit for maximum comfort and compression
- Low profile design fits easily in shoe
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
Small	7" - 8"	4004-01
Medium	8" - 10"	4004-02
Large	10" - 12"	4004-03
X-Large	12" - 14"	4004-04

\*Measure circumference from the base of heel around the ankle

### Indications:

Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints



## CLOSED HEEL ELASTIC ANKLE SLEEVE ankle



- Closed heel, elastic construction
- Provides snug fit for maximum comfort and compression
- Pull-on design for ease of application
- Low profile design fits easily in shoe
- Universal right/left

SIZE	ANKLE CIRCUMFERENCE*	PRODUCT #
Small	7" - 8"	4005-01
Medium	8" - 10"	4005-02
Large	10" - 12"	4005-03
X-Large	12" - 14"	4005-04

\*Measure circumference from the base of heel around the ankle

### Indications:

Mild ankle sprains and strains; Ankle edema; Osteoarthritis of the talocrural and subtalar joints

## DELUXE ANKLE FOOT ORTHOSIS ankle



- Closed heel support provides additional security
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits easily in shoe
- Available in right/left design for best fit

SIZE	MALE	FEMALE	RIGHT	LEFT
	SHOE SIZE	SHOE SIZE		
Small	6 - 7½	7½ - 9	AF1001-01	AF1001-04
Medium	8 - 9 ½	9½ - 11	AF1001-02	AF1001-05
Large	10 - 11½	11½ - 13	AF1001-03	AF1001-06
X-Large	12 +	13½ +	AF1001-07	AF1001-08

### Indications:

Patients suffering from neurological or orthopedic conditions;  
Lower leg instabilities

## ANKLE FOOT ORTHOSIS ankle



- Open back and heel design
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits easily in shoe
- Universal in size

DESCRIPTION	MALE	FEMALE
Right	AF1000-01	AF1000-03
Left	AF1000-02	AF1000-04

### Indications:

Patients suffering from neurological or orthopedic conditions;  
Lower leg instabilities

## FULL SHELL WALKER 2 ankle



- Built in medial and lateral air pumps allow easy operation and instant air bladder control
- Air bladders accommodate swelling and bulky dressings
- Wide foot bed accommodates varying patient foot anatomies
- Textured tread, non-skid sole provides maximum shock absorption and grip
- Rocker bottom mimics natural gait
- Unique adjustable toe guard for increased protection
- Additional pads included for extra comfort where needed

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD	SHORT
X-Small	Up to 4	Up to 5	15500041	15500046
Small	4 – 7	5 – 8	15500042	15500047
Medium	7 – 10	8 – 11	15500043	15500048
Large	10 – 13	11 – 15	15500044	15500049
X-Large	13 +	15 +	15500045	15500050

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

## INLINE PNEUMATIC WALKER 2 ankle



- High quality rigid outer skeleton for superior protection
- Dense foam liner to provide comfort and support
- Shock absorbing embedded sole tread improves patient comfort and helps prevent slips
- Rocker bottom mimics natural gait
- Plastic uprights help reduce ankle rotation and range of motion
- Large air bladder creates a larger surface for contact pressure to help control swelling and ensure a snug, comfortable fit

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD	SHORT
X-Small	Up to 4 ½	Up to 5 ½	15500031	15500036
Small	4 ½ - 7 ½	5 ½ - 8	15500032	15500037
Medium	7 ½ - 10 ½	8 - 11 ½	15500033	15500038
Large	10 ½ - 12 ½	11 ½ - 13 ½	15500034	15500039
X-Large	12 ½ +	13 ½ +	15500035	15500040

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

## PACESETTER™ AIR WALKER ankle



- Inflatable pump design for compression and stabilization
- Unique non-welded air bladder reduces chance of over-inflation to the point of bursting or tearing
- Wide foot plate for patient comfort and stability
- Rocker sole helps promote natural gait
- Heel strap helps prevent foot migration
- Offered in both standard and low top versions

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD OPEN HEEL	LOW TOP OPEN HEEL
Youth	-	-	15780004	-
Small	3 - 6½	4 - 7½	15780005	15800005
Medium	7 - 8½	8 - 10	15780006	15800006
Large	9 - 10½	10½ - 12½	15780067	15800067
X-Large	11 - 16	13 +	15780008	15800008

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles Tendonitis

## PACESETTER™ II WALKER ankle



- Rocker sole helps promote rolling gait pattern
- Double laminated foam wrap is breathable to prevent skin breakdown
- Wide foot plate for patient comfort and stability
- Available in open or closed heel design
- Offered in both standard and low top versions

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD OPEN HEEL	STANDARD CLOSED HEEL	LOW TOP OPEN HEEL
Small	3 - 6½	4 - 7½	CSG53510-A	CSG53550-A	CSG53520-A
Medium	7 - 8½	8 - 10	CSG53510-B	CSG53550-B	CSG53520-B
Large	9 - 10½	10½ - 12½	CSG53510-D	CSG53550-D	CSG53520-D
X-Large	11 - 16	13 +	CSG53510-C	CSG53550-C	CSG53520-C

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles Tendonitis

## INLINE WALKER ankle



- Heavy duty flared uprights are color-coded for easy size identification
- Low profile rocker bottom promotes natural gait
- Wide foot plate for patient comfort and stability
- Treaded sole for patient safety
- Available in open or closed heel design
- Offered in both standard and low top versions

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD OPEN HEEL	STANDARD CLOSED HEEL	LOW TOP OPEN HEEL	LOW TOP CLOSED HEEL
Small	3½ - 7	5 - 8½	15520005	15530005	15540005	15550005
Medium	7½ - 8½	9 - 10	15520006	15530006	15540006	15550006
Large	9 - 10½	10½ - 12	15520067	15530067	15540067	15550067
X-Large	11 +	12½ +	15520008	15530008	15540008	15550008
Pediatric	8 - 12 (Child)	—	15520002	15530002	15540002	15550002

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles Tendonitis

## TRACKER EX WALKER ankle



- Lightweight design for patient comfort
- Low profile rocker bottom sole helps promote natural gait
- Wide foot plate for patient comfort and stability
- Closed heel design helps prevent heel migration
- Extra long wrap provides full coverage
- Available in standard and low top versions

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD	LOW TOP
X-Small	up to 5½	up to 4½	15590002	15580002
Small	6 - 8½	5 - 8½	15590005	15580005
Medium	9 - 10½	9 - 11½	15590006	15580006
Large	11 - 13	12 - 13½	15590067	15580067
X-Large	13½ +	14 +	15590008	15580008

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles Tendonitis

## ONE-PIECE SAMSON WALKER ankle



- Low profile anatomical design promotes natural gait pattern
- Two layer memory foam inlay conforms to foot
- Adjustable straps for even compression
- Wrap material allows skin to breathe and helps to prevent skin breakdown
- Available in standard and low top designs

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD	LOW TOP
X-Small	up to 5	up to 6	SL8001-01	—
Small	5 - 7	6 - 8	SL8001-03	SL8001-33
Medium	7½ - 10	8½ - 11	SL8001-05	SL8001-35
Large	10½ - 12	11½ - 13	SL8001-07	SL8001-37
X-Large	12½ +	13½ +	SL8001-09	—

### Indications:

Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles Tendonitis

## ORIGINAL DARBY SHOE foot & heel



- Semi-rigid metatarsal shank provides controlled forefoot motion
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

SIZE	MALE SHOE SIZE	PRODUCT #	FEMALE SHOE SIZE	PRODUCT #
Small	6 - 8	DS1014-10	4 - 6	DS1013-10
Medium	8½ - 10	DS1014-12	6½ - 8	DS1013-12
Large	10½ - 12	DS1014-14	8½ - 10	DS1013-14
X-Large	12½ - 14	DS1014-16	—	—
Pediatric	12 - 1 (Child)	DS1012-00	—	—

### Indications:

Post-operative foot protection and support; Minor foot fractures

## AEGIS POST-OP SHOE foot & heel



- Semi-rigid rubber sole for support
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

SIZE	MALE SHOE SIZE	PRODUCT #	FEMALE SHOE SIZE	PRODUCT #
Small	6 - 8	2044-01	4 - 6	2044-07
Medium	8½ - 10	2044-02	6½ - 8	2044-08
Large	10½ - 12	2044-03	8½ - 10	2044-09
X-Large	12½ - 14	2044-04	—	—

### Indications:

Post-operative foot protection and support; Minor foot fractures

## BLACK POST-OP SHOE foot & heel



- Semi-rigid rubber sole for support
- Black nylon upper for comfort
- Loop-Lock closure for easy application
- Universal right/left

SIZE	MALE SHOE SIZE	PRODUCT #	FEMALE SHOE SIZE	PRODUCT #
Small	6 - 8	A152305	4 - 6	A152405
Medium	8½ - 10	A152306	6½ - 8	A152406
Large	10½ - 12	A152307	8½ - 10	A152407
X-Large	12½ - 14	A152308	—	—

### Indications:

Post-operative foot protection and support; Minor foot fractures

## CANVAS POST-OP SHOE foot & heel



- Rigid rubber sole provides controlled forefoot motion
- Canvas upper for durability
- Hook and loop closure
- Universal right/left

SIZE	MALE SHOE SIZE	PRODUCT #	FEMALE SHOE SIZE	PRODUCT #
Small	6 - 8	2041-01	4 - 6	2041-07
Medium	8½ - 10	2041-02	6½ - 8	2041-08
Large	10½ - 12	2041-03	8½ - 10	2041-09
Pediatric	6 - 8 (Child)	2041-00	—	

### Indications:

Post-operative foot protection and support; Minor foot fractures

## DELUXE WOODEN POST-OP SHOE foot & heel



- Rigid wood sole with rubber bottom for support
- Perforated suede upper for comfort and durability
- Available in hook and loop or lace up closure design
- Universal right/left

SIZE	MALE SHOE SIZE	HOOK/LOOP CLOSURE	LACE UP CLOSURE
Child	2 - 4½	2007-01	2006-01
Small	5 - 6½	2007-02	2006-02
Medium	7 - 8½	2007-03	2006-03
Large	9 - 11	2007-04	2006-04
X-Large	11½ - 13	2007-05	2006-05
SIZE	FEMALE SHOE SIZE	HOOK/LOOP CLOSURE	LACE UP CLOSURE
Child	2 - 3½	2007-06	2006-06
Small	4 - 6	2007-07	2006-07
Medium	6½ - 8	2007-08	2006-08
Large	8½ - 10	2007-09	2006-09
X-Large	10½ - 12	2007-10	2006-10

### Indications:

Post-operative foot protection and support; Minor foot fractures

## CANVAS WOODEN POST-OP SHOE foot & heel



- Rigid wood sole with rubber bottom for support
- Blue canvas upper
- Available in hook and loop or lace up closure design
- Universal right/left

SIZE	MALE SHOE SIZE	HOOK/LOOP CLOSURE	LACE UP CLOSURE
Small	5 - 6½	2020-11	2020-01
Medium	7 - 8½	2020-12	2020-02
Large	9 - 11	2020-13	2020-03
X-Large	11½ - 13	2020-14	2020-04
SIZE	FEMALE SHOE SIZE	HOOK/LOOP CLOSURE	LACE UP CLOSURE
Small	4 - 6	2021-11	2021-01
Medium	6½ - 8	2021-12	2021-02
Large	8½ - 10	2021-13	2021-03
X-Large	10½ - 12	2021-14	2021-04

### Indications:

Post-operative foot protection and support; Minor foot fractures

## DARCO® MEDSURG™ SHOE foot & heel



- MetaShank™ protection provides rigid control under metatarsal heads
- Strapless forefoot closure expands to accommodate bulky dressings
- Square toe design serves as a bumper for additional protection
- Rocker sole reduces plantar pressure on forefoot and heel
- Ankle strap helps eliminate heel migration
- Universal right/left

SIZE	MALE SHOE SIZE	PRODUCT #	FEMALE SHOE SIZE	PRODUCT #
Small	6 - 8	DS1016-10	4 - 6	DS1015-10
Medium	8½ - 10	DS1016-12	6½ - 8	DS1015-12
Large	10½ - 12	DS1016-14	8½ - 10	DS1015-14
X-Large	12½ - 14	DS1016-16	—	—

### Indications:

Post-operative foot protection and support; Minor foot fractures

## DARCO® ORTHOWEDGE™ SHOE foot & heel



- Helps to reduce weight bearing pressure from the forefoot with 15° dorsi-angle
- Removable 10 mm insole allows customization for off-loading specific areas
- Square toe design serves as a bumper for additional protection
- Hook and loop closure for easy application
- Universal right/left

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	PRODUCT #
X-Small	—	4 - 7	HS1005-08
Small	6 - 8	7½ - 10	HS1005-10
Medium	8½ - 10	10½ - 13	HS1005-12
Large	10½ - 12	13½ +	HS1005-14
X-Large	12½ - 14	—	HS1005-16

### Indications:

Post-operative foot protection and support; Minor foot fractures



## COMFO-EZE ROCKER CAST SHOE foot & heel



- Constructed of durable navy canvas material
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

SIZE	LENGTH	WIDTH	PRODUCT #
Pediatric	5 ½"	2 ½"	2004-00
X-Small	7"	3 ½"	2004-01
Small	8 ½"	4 ½"	2004-02
Medium	10"	4 ¾"	2004-03
Large	11 ½"	5 ¼"	2004-04
X-Large	13"	5 ½"	2004-05

### Indications:

Lower extremity cast protection

## CANVAS CAST SHOE foot & heel



- Constructed of navy canvas material
- Padded canvas footbed for maximum comfort
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

SIZE	LENGTH	WIDTH	PRODUCT #
Child	5 ½"	2 ½"	2035-01
Small	8 ½"	4 ½"	2035-02
Medium	10"	4 ¾"	2035-03
Large	11 ½"	5 ¼"	2035-04
X-Large	13"	5 ½"	2035-05

### Indications:

Lower extremity cast protection

## BLACK CANVAS CAST SHOE foot & heel



- Constructed of durable black canvas material
- Non-slip rubber sole for patient safety
- Rocker bottom promotes rolling gait pattern
- Heel strap helps prevent foot migration
- Hook and loop closure for easy application

SIZE	LENGTH	WIDTH	PRODUCT #
X-Small	7"	3 ½"	2033-01
Small	8 ½"	4 ½"	2033-02
Medium	11 ½"	5 ¼"	2033-03
Large	13"	5 ½"	2033-04

### Indications:

Lower extremity cast protection

## VINYL CAST BOOT foot & heel



- Constructed of durable vinyl material
- Rocker bottom promotes rolling gait pattern
- Closed heel design to prevent foot migration
- Hook and loop closure for easy application

SIZE	LENGTH	WIDTH	PRODUCT #
X-Small	7"	3 ½"	2027-01
Small	8 ½"	4 ½"	2027-02
Medium	10"	4 ¾"	2027-03
Large	11 ½"	5 ¼"	2027-04
X-Large	13"	5 ½"	2027-05

### Indications:

Lower extremity cast protection

## DARCO® SLIMLINE CAST BOOT foot & heel



- Square toe design serves as a bumper for protection and provides room for larger casts
- Stretchable upper material helps boot conform to cast allowing precise fit
- Rocker sole promotes rolling gait pattern
- Padded insole to cushion patient's foot
- Strapless forefoot closure adjusts easily to cast size

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	PRODUCT #
X-Small	—	3 - 5	2039-01
Small	—	5½ - 8	2039-02
Medium	6 - 8	8½ - 10	2039-03
Large	8½ - 11	10½ - 13	2039-04
X-Large	11½ - 14	13+	2039-05
SIZE	CHILD SHOE SIZE		PRODUCT #
Child/Small	8 - 10½		2039-98
Child/Medium	11 - 1½		2039-99
Child/Large	2 - 4		2039-00

### Indications:

Lower extremity cast protection

## DEROYAL® NIGHT SPLINT foot & heel



- Dorsiflexion design applies stretch for effective pain relief
- Lightweight and low profile design for patient comfort
- Adjustable straps to vary amount of stretch
- Removable foam liner for laundering
- Optional foam wedges provided for additional stretch
- Kit features: night splint, ice pack, exercise band, heel cups, foot control strap and self-management guide

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	SPLINT	KIT
Small	4 - 6	5 - 7	NS2000-10	SC2000-10
Medium	6½ - 9	7½ - 10	NS2000-12	SC2000-12
Large	9½ +	10½ +	NS2000-14	SC2000-14

### Indications:

Plantar fasciitis; Achilles tendonitis; Tendonosis

## DORSAL NIGHT SPLINT foot & heel



- Thermoplastic design holds foot in 90° neutral position
- Neoprene straps with hook and loop closure
- Plastic free bottom and non-skid transfer surface
- Low profile design allows for a cooler, more comfortable fit

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	PRODUCT #
Small/Medium	4 - 9	5 - 10	NS2000-20
Large/X-Large	9½ +	10½ +	NS2000-21

### Indications:

Plantar fasciitis; Achilles tendonitis; Tendonosis

## SILICONE HEEL CUPS foot & heel



- Constructed of silicone for long-lasting support and comfort
- Soft blue spot helps relieve heel strike fatigue and discomfort through shock absorption
- Packaged 1/Pair

SIZE	PRODUCT #
Small	HC2000-10
Medium	HC2000-12
Large	HC2000-14

### Indications:

Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

## GREY SPOT HEEL CUPS foot & heel



- Constructed of viscoelastic polymers to dissipate heel strike energy
- “U-shaped” design helps divert weight bearing forces away from plantar fascia
- Packaged 1/Pair

---

SIZE	PRODUCT #
Regular	HC1000-01
Large	HC1000-02

### Indications:

Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

**CRUTCHES** ambulatory aids



- Light-weight and sturdy aluminum construction
- Rubber tips for maximum surface traction
- Padded underarm rests and hand grips
- Push-button height adjustment

SIZE	HEIGHT	QTY	PRODUCT #
Short	4'8" – 5'2"	10 Pr/Cs	CR4000-LS
Regular	5'2" – 5'10"	10 Pr/Cs	CR2000-LS
Tall	5'10" – 6'6"	10 Pr/Cs	CR3000-LS

**FOLDING WALKERS** ambulatory aids



- Light-weight and sturdy aluminum construction
- Rubber tips for maximum surface grip
- Padded hand grips
- Push-button height adjustment on each leg
- Push-button folding for easy storing/carrying
- Optional 3" wheels

SIZE	QTY	PRODUCT #
Standard	2/Cs	WLK02
Standard w/wheels	2/Cs	WLK02-W


**CANES** ambulatory aids



- Light-weight and sturdy aluminum construction
- Rubber tip for maximum surface traction
- Padded hand grip
- Push-button height adjustment

SIZE	QTY	PRODUCT #
Standard	20/Cs	CA298

**PATIENT CARE** | Pressure Ulcer Prevention | Foam Positioners |  
Specialty Nursing



## ANKLE CONTRACTURE BOOT pressure ulcer prevention



Fleece



- Static ankle foot orthosis used for positioning and pressure reduction
- Structured to completely suspend the heel
- Available in fleece or foam material with or without boot sole
- Softgoods completely cover splint leaving no exposed metal against skin
- "No strap" closures make it easy for patients to apply and remove the boot
- Anti-rotation bar, padded toe piece and mesh laundry bag included
- Accessories include hip abduction bar, boot sole and extra softgoods
- Universal right/left

DESCRIPTION	FOAM	FLEECE	ACCESSORIES	PRODUCT #
Ankle Contracture Boot	4302	4307	Extra Liner - Foam	4300K
Ankle Contracture Boot w/ Boot Sole	4301	4306	Extra Liner - Fleece	4305K
Ankle Contracture Boot w/ Boot Sole & Extra Liner	4300	4305	Hip Abduction Bar	4300/05AB
			Boot Sole	4300/05WK

### Indications:

Hip rotation; Foot drop; Not recommended for severe plantar flexion contractures, open wounds, severe spasticity or thrombosis

### Sizing Guide

Add letter to end of product # to indicate size

Measure calf circumference 8" up from bottom of heel

- A\* = 8" - 12"
- B = 10" - 14"
- C = 14" - 16"
- D = 16" - 21"
- E = 16" - 21"

Measure foot circumference around the widest part of foot

- A\* = 6" - 8"
- B = 7" - 9"
- C = 9" - 11"
- D = 11" - 13"
- E = 13" +

\*Size A available in 4301 style only

## PRUVENTOR™ HEEL OFF-LOADING DEVICE pressure ulcer prevention



- Soft, fluid and stain resistant material can be cleaned if soiled
- Pressure-absorbing filling creates fluid-like interior and absorbs pressure away from bony prominences
- Soft, elastic straps with hook tabs for easy and secure placement (Long boot comes with additional strap for added security around lower leg)
- Anti-rotation wedge easily attaches to either side of the boot
- Compatible with sequential and intermittent compression devices (SCD/ICD) with easy exits for SCD/ICD tubing (Long boot has additional exits on either side of the boot)

DESCRIPTION	LENGTH	QTY	PRODUCT #
Universal Short	9"	1 Each	M70-100S
Universal Long	13"	1 Each	M70-100L
PRUventor XL	13"	1 Each	M70-100B

### Indications:

Heel decubitus ulcer prevention; Foot drop; Foot rotation



## HEEL PILLOW pressure ulcer prevention



- Soft, cotton padding provides maximum heel cushioning to help protect against skin breakdown
- Cotton strap(s) with hook and loop closure for easy application and adjustment
- Available in two-strap, criss-cross version for more security (M3042)

SIZE	QTY	PRODUCT #
Universal	1 Pair	M3040
Universal	1 Pair	M3042

### Indications:

Heel decubitus ulcer prevention; Skin shearing and/or skin friction

## DEROYAL® HEEL SUSPENSION BOOT pressure ulcer prevention



- Smooth, breathable foam construction
- Sloped wedge is precisely placed inside boot to completely suspend the heel
- Ventilation holes help promote air circulation and dissipate heat
- Soft straps with D-ring/hook closures
- Posterior panel helps reduce friction against bed surface
- Extra foam block with adhesive can be attached to side of boot to help prevent foot rotation

DESCRIPTION	QTY	PRODUCT #
Smooth w/Attached Wedge	6/Case	M60-300

### Indications:

Heel decubitus ulcer prevention; Foot drop; Foot rotation

## BUNNY BOOT pressure ulcer prevention




- Resilient, breathable polyfiber liner allows air circulation for increased comfort
- Adjustable foam wedge at base of frame
- Rigid plastic frame holds foot at a 90° angle and can act as a foot rotator with internal or external rotation of 10° or 20° arcs
- Soft straps with hook and loop/D-ring closures

SIZE	QTY	PRODUCT #
Universal	1 Each	M0501D01
Universal	24/Case	M0501D05

### Indications:

Toe and heel pressure reduction; Foot drop; Static stretching of foot

## FOAM HEEL & ANKLE PROTECTORS pressure ulcer prevention



- Breathable, convoluted foam construction for patient comfort
- Available in high or low-density foam
- Foam strap with hook and loop closure for easy application and adjustment

SIZE	DESCRIPTION	QTY	PRODUCT #
Universal	High-density	1 Pair	M190FF
Universal	High-density	24 Pair/Case	M190FF-24
Universal	Low-density	20 Pair/Case	M190FFE

### Indications:

Decubitus ulcer prevention; Skin shearing and/or skin friction

## HEEL & ELBOW PROTECTORS - FULL FOOT DESIGN pressure ulcer prevention



- Deluxe or standard, resilient polyfiber construction
- Ventilation holes help promote air circulation (heel only)
- Foam strap with hook and loop closure for easy application and adjustment

DELUXE	SIZE	QTY	PRODUCT #
Heel	Universal	1 Pair	M3014D
Heel	Universal	6 Pair/Case	M3014D-6
Elbow	Universal	1 Pair	M3006

STANDARD	SIZE	QTY	PRODUCT #
Heel	Universal	1 Pair	M3015

### Indications:

Decubitus ulcer prevention; Skin shearing and/or skin friction

## GELBO® HEEL & ELBOW PROTECTORS pressure ulcer prevention



- Nylon/spandex sleeve conforms around the heel or elbow to hold product in place
- Lined with natural mineral oil-based gel formula for shock absorption and cushioning
- Available in standard and large sizes

SIZE	QTY	PRODUCT #
Standard (>110 lbs)	1 Pair	M3022
Large (<110 lbs)	1 Pair	M3022-L

### Indications:

Decubitus ulcer prevention; Skin shearing and/or skin friction

## PADDED HEEL & ELBOW PROTECTORS pressure ulcer prevention



- Knitted sleeve conforms around heel or elbow to help hold product in place
- Breathable foam or resilient polyfiber pad sewn into cushioned sleeve for maximum pressure protection with no irritating seam


---

SIZE	QTY	FOAM PAD
Universal	12 Pair/Case	M3000U
Universal	36 Pair/Case	M3001U
Bariatric	12 Pair/Case	M3000UB

### Indications:

Decubitus ulcer prevention; Skin shearing and/or skin friction

## CONCAVE HIP ABDUCTION PILLOW - HIGH-DENSITY foam positioners



- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Long foam straps with hook closure for easy application and adjustment
- Available in compressed packaging for efficient storage

SIZE	DIMENSIONS				QTY	COMPRESSED
	A	B	C	D		
Infant	4½"	5"	2"	2"	10/Case	M60-076-I
X-Small	9"	10"	3"	4"	5/Case	M60-076-XS
Small	14"	18"	7"	5"	4/Case	M60-076-S
Medium	16½"	21"	7½"	6"	4/Case	M60-076-M
Large	18"	24"	9"	6"	4/Case	M60-076-L


  

SIZE	DIMENSIONS				QTY	NON-COMPRESSED
	A	B	C	D		
Small	14"	18"	7"	5"	2/Case	M60-078-S
Medium	16½"	21"	7½"	6"	2/Case	M60-078-M
Large	18"	24"	9"	6"	2/Case	M60-078-L

### Indications:

Post surgical hip abduction

## CONCAVE HIP ABDUCTION PILLOW - WIDE foam positioners




- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Wider angle for increased abduction
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

SIZE	DIMENSIONS				QTY	PRODUCT #
	A	B	C	D		
Small	16"	16"	10"	7"	5/Case	M60-024-S
Medium	18"	20"	10"	9"	5/Case	M60-024-M
Large	20"	24"	10"	9"	5/Case	M60-024-L

### Indications:

Post surgical hip abduction

## CONCAVE HIP ABDUCTION PILLOW - WIDE W/ LINER foam positioners




- Breathable, high-density foam construction
- Concave sides hold legs in secure position
- Wider angle for increased abduction
- Polyfiber liner provides additional patient comfort
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

SIZE	DIMENSIONS				QTY	PRODUCT #
	A	B	C	D		
Small	16"	16"	10"	7"	5/Case	M60-025-S
Medium	18"	20"	10"	9"	5/Case	M60-025-M
Large	20"	24"	10"	9"	5/Case	M60-025-L

### Indications:

Post surgical hip abduction

## FLAT HIP ABDUCTION PILLOW - NARROW foam positioners




- Breathable, high-density foam construction
- Long foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

SIZE	DIMENSIONS				QTY	PRODUCT #
	A	B	C	D		
Small	10"	16"	6"	7"	5/Case	M60-031-S
Medium	12"	21"	6"	7"	5/Case	M60-031-M
Large	14"	23"	6"	7"	5/Case	M60-031-L

### Indications:

Post surgical hip abduction

## HIP ABDUCTION PILLOW - ADJUSTABLE foam positioners




- Breathable, high-density foam construction
- 4" tear-off blocks on both ends for length adjustment
- Long, perforated foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

SIZE	DIMENSIONS				QTY	PRODUCT #
	A	B	C	D		
Universal	10"	18" - 27"	6"	5"	4/Case	M60-224

### Indications:

Post surgical hip abduction

## ARM CRADLE foam positioners




- Breathable, high-density foam construction
- Stabilizes arm in upright position
- Ventilation holes promote air circulation, dissipate heat and allow use of ice packs to reduce edema
- Compressed packaging for efficient storage

SIZE	DIMENSIONS			QTY	PRODUCT #
	A	B	C		
Universal	12½"	9"	15½"	5/Case	M60-034

### Indications:

Arm elevation; Post-operative stabilization

## FOOT CRADLE foam positioners




- Breathable, high-density foam construction
- Stabilizes foot in upright position
- Ventilation holes promote air circulation, dissipate heat and allow use of ice packs to reduce edema
- Compressed packaging for efficient storage

SIZE	DIMENSIONS		QTY	PRODUCT #
	A	B		
Universal	14"	9"	10/Case	M60-045

### Indications:

Foot drop support; Post-operative stabilization

## FOOT DROP SUPPORT foam positioners




- Breathable, high-density foam construction
- Adjustable anti-rotation strap holds foot in neutral position and keeps heel suspended
- Hook receivable foam strap
- Compressed packaging for efficient storage

SIZE	DIMENSIONS		QTY	PRODUCT #
	A	B		
Universal	8"	10"	5/Case	M60-018

### Indications:

Decubitus ulcer prevention; Foot drop support

## LEG POSITIONER foam positioners




- Breathable, high-density standard or flat foam construction
- Available with or without 45° angle
- Stabilizes knee and lower leg while promoting better circulation
- Foam strap with hook closure for easy application and adjustment (concave only)
- Compressed packaging for efficient storage

SIZE	DIMENSIONS			QTY	PRODUCT #
	A	B	C		
Standard	11"	29"	11"	2/Case	M60-015A
Standard w/ 45° angle	11"	28"	10"	1 Each	M60-015B
Flat w/ 45° angle	8"	31"	15"	2/Case	M60-015CD

### Indications:

Pressure relief of heel, ankle, mid-foot and mid-calf; Decubitus ulcer prevention; Static back pain relief

## LEG POSITIONER W/ FOOT DROP STRAP foam positioners




- Breathable, high-density foam construction
- Hook receivable foot drop strap holds foot in neutral position
- Stabilizes knee and lower leg while promoting better circulation
- Foam straps with hook closure for easy application and adjustment
- Compressed packaging for efficient storage

SIZE	DIMENSIONS			QTY	PRODUCT #
	A	B	C		
Universal	15"	32"	11"	4/Case	M60-015

### Indications:

Pressure relief of heel, ankle, mid-foot and mid-calf; Decubitus ulcer prevention; Static back pain relief

## BODY ALIGNER foam positioners



- Breathable, high-density or low density (M60-032E) foam construction
- Helps maintain proper positioning while reducing pressure from sacrum
- Available in compressed packaging for efficient storage or non-compressed option
- Available in covered version (M60-032C, M60-032D) with non-skid bottom that can be cleaned and re-used

SIZE	A	B	C	QTY	COMPRESSED	COMPRESSED
Universal	24"	11"	8"	5/Case	M60-032	—
Universal	24"	11"	8"	10/Case	M60-032B	—
Universal	26"	27"	12"	2/Case	M60-033	—
Universal	21"	10½"	8"	6/Case	—	M60-032E
Universal	24"	11"	8"	2/Case	—	M60-032C
Universal	24"	11"	6"	4/Case	—	M60-032D

### Indications:

Decubitus ulcer prevention; Positioning



**SURGICAL BRA** specialty nursing



- Constructed of Lycra® material
- Front hook and loop closure for easy access to dressings and application
- Wide adjustable shoulder straps for maximum patient comfort

SIZE	MEASUREMENT	PRODUCT #
Small	32" - 35"	M5001-S
Medium	36" - 39"	M5001-M
Large	40" - 43"	M5001-L
X-Large	44" - 47"	M5001-XL
XX-Large	48" - 51"	M5001-XXL
XXX-Large	52" - 55"	M5001-XXXL

**Indications:**

Post operative compression; Dressing holder

**TRACH TUBE HOLDER** specialty nursing




- Soft, foam holder secures most tracheostomy tubes
- Available with wide or narrow fasteners for secure fit (adult only)
- Adjustable to accommodate different neck circumferences
- 12/Case

SIZE	FASTENER	NECK CIRCUMFERENCE	PRODUCT #
Adult	Wide	up to 23"	M1150
Adult	Narrow	up to 20"	M1151
Pediatric	Narrow	up to 10"	M1152

**Indications:**

Patients with tracheostomy

**CATHETER STRAP** specialty nursing



- Soft, elastic construction for patient comfort
- Excellent substitute to messy, skin-irritating tape
- Hook and loop fasteners for easy application and placement
- Available in standard or extra long length
- Ideal for any area where foley catheters are applied

SIZE	DIMENSIONS	QTY	ELASTIC
Universal	22"L x 2"W	36/Case	M1133
Universal, X-Long	38"L x 2"W	24/Case	M1133-38

**Indications:**

Securing foley catheters; Catheter positioning

## ELBOW IMMOBILIZER specialty nursing




- Soft, perforated foam construction for patient comfort and breathability
- Metal T-bar stays for optimal immobilization and strength
- Prevents flexion of arm to prevent pulling NG tubes, tracheostomy tubes or I.V. lines
- Elastic straps with hook closures for fast and accurate placement (Adult only)

SIZE	DIMENSIONS	QTY	PRODUCT #
Adult	14"L x 17"W	1 Pair	M2035
Infant	11"L x 5"W	12 Eaches	M7035-XS
Child	12½"L x 7½"W	1 Pair	M7035-S
Child	12½"L x 7½"W	12 Eaches	M7035-SB
Child	12"L x 7½"W	12 Eaches	M7035-M
Child	17"L x 7½"W	1 Pair	M7035-L

### Indications:

Patients at risk of unplanned extubations; Patients who disrupt medical treatment

## CARDIAC CUSHION specialty nursing



- Heart diagram printed on pillow with permanent marker included for bedside education purposes
- Cotton cover with polyester fill provides superior cushioning
- Packaged 10/Case

SIZE	PRODUCT #
Adult	M13-100
Pediatric	M13-500

\*Contact your local DeRoyal sales representative for custom logo information.

### Indications:

Post-cardiac surgery; Support during expectoration; Sternal support to aid in healing; Helps reduce pain when coughing and deep breathing; Marketing tool\*

## DISPOSABLE TELEMETRY POUCH specialty nursing



- Durable, non-woven material
- Neck and torso tie straps for easy application and adjustment
- Open pouch design
- 50/Case

DIMENSIONS	PRODUCT #
7½"L x 6½"W	M4209D06

### Indications:

Patients with telemetry units

## POLYESTER TELEMETRY POUCH specialty nursing



- Durable polyester pouch with foam neck strap for patient comfort
- Hook receivable material for easy application and adjustability
- 50/Case

<b>DIMENSIONS</b>	<b>PRODUCT #</b>
6½" L x 5" W	M13-055

### Indications:

Patients with telemetry units

## BED RAIL PADS specialty nursing



- Foam padded, vinyl covered pads help prevent patient injury
- Vinyl covering wipes clean to help prevent cross contamination issues
- Available in 2-piece (side rails only) or 4-piece set (side rails and head/foot boards)
- Hook and loop closures (M3000D08) or elastic straps (M3002D09) for side rails and snap buttons for head/foot boards for easy application and removal

DESCRIPTION	SIDE RAIL DIMENSIONS	HEAD/FOOT BOARD DIMENSIONS	PRODUCT #
2 Piece Set	69" L x 12" H	—	M3000D08
2 Piece Set	32" L x 12" H	—	M3002D09
4 Piece Set	69" L x 12" H	32" L x 12" H	M3002D08

### Indications:

Bed side rail injury and bed exit

## DISPOSABLE FOAM BED RAIL PADS specialty nursing



- Foam construction provides protective padding to help prevent patient injury
- Hook and loop closures for easy application and removal

DESCRIPTION	SIDE RAIL DIMENSIONS	HEAD/FOOT BOARD DIMENSIONS	PRODUCT #
2 Piece Set	60" L x 18" H	—	M10-080

### Indications:

Bed side rail injury and bed exit

## CERVICAL TENSION PILLOW specialty nursing



- Constructed of soft foam with cotton cover
- Designed to help relieve pressure and tension on cervical spine
- Anatomically shaped to fit back of neck and head area

---

SIZE	PRODUCT #
Universal	9222-00

### Indications:

Cervical spine positioning and support; Neck tension relief

## CERV-O-REST PILLOW specialty nursing



- Constructed of soft foam with cotton cover
- Compresses to fill in space behind neck when in a supine position

---

SIZE	PRODUCT #
Universal	9221-00

### Indications:

Cervical spine positioning and support

## ORTHOPEDIC FELT specialty nursing




- Ideal for extra padding

---

MEASUREMENT	PRODUCT #
36" x 21" x ¼"	8829-01

## INVALID RING specialty nursing



- Soft foam construction provides pressure-free comfort while sitting
- Comes with cloth cover

SIZE	DIMENSIONS			PRODUCT #
	A	B	C	
Universal	14½"	12"	2½"	M10-070- C

### Indications:

Post-natal discomfort; Post-operative pain; Hemorrhoids; Bruised coccyx

## ELASTIC BANDAGES specialty nursing



- Durable elastic construction
- Available in standard length (5 yds) and double length (10 yds)
- Packaged individually
- Non-sterile
- Available in hook and loop or clip closure
- Premium option features hook closure on both ends of roll

SIZE	QUANTITY	PREMIUM		ECONO	
		CLIP CLOSURE	HOOK & LOOP Closure	CLIP CLOSURE	HOOK & LOOP Closure
2" x 5 yds	10 rolls/Case	EB2	EBV2	DEB2	DEBV2
3" x 5 yds	10 rolls/Case	EB3	EBV3	DEB3	DEBV3
4" x 5 yds	10 rolls/Case	EB4	EBV4	DEB4	DEBV4
6" x 5 yds	10 rolls/Case	EB6	EBV6	DEB6	DEBV6
4" x 10 yds	6 rolls/Case	HEB4	—	—	—
6" x 10 yds	6 rolls/Case	HEB6	—	—	—

## JETSTREAM® temperature therapy



- Patented Automatic Temperature Control (ATC™) measures and regulates temperature to provide therapeutic hot or cold therapy without manual adjustment
- Intuitive LED light lets user know when the unit is in the therapeutic zone and advises when they should add ice or hot water
- Portable system is suitable for any medical setting or treatment at home
- Lightweight, low profile design and Quick Start Guide are patient friendly
- Unit provides hot and cold therapy

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

### TEMPERATURE THERAPY UNIT T700

COMBO UNITS (WITH BLANKET)	NON-STERILE	STERILE
Hip/Lumbar	T725NS	T725ST
Knee/Shoulder	T735NS	T735ST
Knee/Shoulder (w/attached straps)	T736NS	T736ST
Shoulder/Knee (w/attached straps)	T737NS	T737ST

## COOLJET™ temperature therapy



- Patented ATC™ measures and regulates temperature to provide therapeutic cold therapy without manual adjustment
- Intuitive LED light lets user know when in therapeutic zone and when to add ice
- Portable system suitable for any medical setting or treatment at home
- Therapy Blankets optimize temperature distribution, insulation and comfort
- Blanket Wraps increase insulation and ensure precise placement of blankets
- Intuitive screw top opening

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

### TEMPERATURE THERAPY UNIT T800

COMBO UNITS (WITH BLANKET)	NON-STERILE	STERILE
Hip/Lumbar	T825NS	T825ST
Knee/Shoulder	T835NS	T835ST
Knee/Shoulder (w/attached straps)	T836NS	T836ST
Shoulder/Knee (w/attached straps)	T837NS	T837ST



## HOT/COLD BLANKETS temperature therapy



- Therapy blankets optimize temperature distribution, insulation and comfort
- Available in non-sterile or sterile versions
- Universal right/left

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma



DESCRIPTION	DIMENSIONS	QTY	NON-STERILE	STERILE
Arm/Leg	10" x 5"	1 Each 10/Case	T5015NS T5015NS-10	T5015ST T5015ST-10
Hip/Lumbar	10 ½" x 11 ½"	1 Each 10/Case	T5025NS T5025NS-10	T5025ST T5025ST-10
Knee/Shoulder	10 ½" x 11"	1 Each 10/Case	T5035NS T5035NS-10	T5035ST T5035ST-10
Knee/Shoulder – Long	12 ½" x 11"	1 Each 10/Case	T5035NS-L T5035NS-L-10	T5035ST-L T5035ST-L-10
Knee/Shoulder (w/attached straps)	10 ½" x 11"	1 Each 10/Case	T5036NS T5036NS-10	T5036ST T5036ST-10
Knee/Shoulder – Long (w/attached straps)	12 ½" x 11"	1 Each 10/Case	T5036NS-L T5036NS-L-10	T5036ST-L T5036ST-L-10
Shoulder/Knee (w/attached straps)	10 ½" x 11"	1 Each 10/Case	T5037NS T5037NS-10	T5037ST T5037ST-10
Shoulder/Knee – Long (w/attached straps)	12 ½" x 11"	1 Each 10/Case	T5037NS-L T5037NS-L-10	T5037ST-L T5037ST-L-10
Hand/Wrist	3 ½" x 5"	1 Each 10/Case	T5055NS T5055NS-10	T5055ST T5055ST-10
Cervical	10" x 21"	1 Each 10/Case	T5075NS T5075NS-10	T5075ST T5075ST-10
Foot/Ankle	10" x 11 ½"	1 Each 10/Case	T5085NS T5085NS-10	T5085ST T5085ST-10

## HOT/COLD THERAPY WRAPS temperature therapy



- Wrap keeps blanket well insulated and provides a simple way to ensure precise blanket placement

DESCRIPTION	DIMENSIONS	QTY	PRODUCT #
Chest/Lumbar	12" x 58"	1 Each	T6100
Knee	15" x 27"	1 Each	T6110
Shoulder	15" x 29"	1 Each	T6120
Hips	17" x 58"	1 Each	T6160



## HOT/COLD FOAM THERAPY WRAPS temperature therapy




- Laminated foam construction keeps hot/cold pack insulated
- Anatomically designed to hold hot/cold pack in place
- Foam is compressive and can be tightened to apply even compression and help reduce swelling
- Universal right/left

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

DESCRIPTION	PACKS	PRODUCT #
Shoulder	2	9381-00
	4	9381-10
Knee	2	9383-00
	4	9383-10
Knee w/Straps	4	9387-00
Back	1	9384-00
Ankle	1	9385-00
Wrist	1	9386-00

## CHILL-N-WRAP ICE BAG temperature therapy



- Constructed of soft outer cover with moisture absorbing middle layer
- Leak-proof plastic bladder prevents patient exposure
- Plastic closure clip is completely enclosed
- Comfortable opening holds shape while ice is being added
- Packaged 10/Case

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

SIZE	POUCH DIMENSIONS	STRAP LENGTH	PRODUCT #
Large	6" x 16"	16 1/2"	8726-50

## HOT/COLD PACK temperature therapy



- Reusable for hot or cold therapy
- Packaged each or 20/case

**Indications:** Pain and soft tissue swelling as a result of surgery or trauma

MEASUREMENT	EACH	CASE
4 1/2" x 10"	9380-00	9380-20

## PEDIATRIC INDEX

### Head/Neck

Product Name	Page#
Comfo-EZE Cervical Collar	7
Economy Foam Collar	8
DeRoyal Rigid Collar	8
EMT Select Extrication Collar	9
One-Piece Extrication Collar	9
XTW Extended Wear Collar	9
Capital™ Collar	10

### Clavicle/Shoulder

Product Name	Page#
Premium Clavicle Strap	24
Heavy-Duty Clavicle Strap	24
3-Way Clavicle Strap	24
Universal Sling Style Shoulder Immobilizer	26
Deluxe Sling & Swathe	27

### Arm/Elbow

Product Name	Page#
Specialty Arm Sling	32
Envelope Arm Sling	32
Light Blue Arm Sling	32
Narrow Pouch Arm Sling	33
Deep Pocket Arm Sling	33
Universal Arm Sling	34

### Wrist/Hand

Product Name	Page#
Suede Leatherette Wrist & Wrist/Forearm Splint	42
Canvas Wrist & Wrist/Forearm Splint	42
Black Wrist Splint	44
Elastic Wrist Splint	46
Finger Dressing/Finger Cage	52
Wrist & Forearm Splint	53
Colles' Splint	53
Plastic Colles' Splint	53
Short Metacarpal Splint	54
Long Metacarpal Splint	54
Elbow Splint	55
Posterior Leg Splint	55

### Knee/Leg

Product Name	Page#
Sized Black Foam Knee Immobilizer	62
Universal Knee Support	65
Full Foot Bucks Traction Boot	68
Sized Bucks Traction Boot	69
Bucks Traction Boot w/Triangle Buckle	69

### Ankle/Foot & Heel

Product Name	Page#
Air Ankle Stirrup	72
Pacesetter™ Air Walker	76
Inline Walker	77
Original Darby Shoe	79
Canvas Post-Op Shoe	80
Deluxe Wooden Post-Op Shoe	80
Comfo-EZE Rocker Cast Shoe	82
Canvas Cast Shoe	82
Darco® SlimLine Cast Boot	83

### Patient Care

Product Name	Page#
Concave Hip Abduction Pillow	92
Trach Tube Holder	96
Elbow Immobilizer	97
Cardiac Cushion	97

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
1000159	28	1020025	42	1050-20	6	11204	50
1000251	6	1020115	42	1050-21	6	11205	50
1000275	43	1020127	61	1050-22	6	1120917	58
1000285	43	1020159	28	1050-23	6	11211000	9
1000301	6	1020167	61	1050-24	6	11212000	9
1000351	6	1020197	61	1050-30	6	11221000	9
1000401	6	1020227	61	1050-31	6	11222000	9
1001-01	7	1020247	61	1050-32	6	11231000	9
1001-02	7	1020251	6	1050-33	6	11232000	9
1001-03	7	1020301	6	1050-34	6	11241000	9
1001-04	7	1020351	6	1050025	42	11242000	9
1001-05	7	1020401	6	1050115	42	11251000	9
1002-01	7	1030025	42	1050127	61	11252000	9
1002-02	7	1030115	42	1050159	28	11261000	9
1002-03	7	1030127	61	1050167	61	11262000	9
1002-04	7	1030159	28	1050197	61	11271000	9
1002-05	7	1030167	61	1050227	61	11272000	9
1003-01	7	1030197	61	1050247	61	11281000	9
1003-02	7	1030227	61	1056-00	6	11282000	9
1003-03	7	1030247	61	1056-02	6	11291000	9
1003-04	7	1030251	6	1056-03	6	11292000	9
1003-05	7	1030301	6	1056-04	6	11301	53
1004-01	7	1030351	6	1056-05	6	11302	53
1004-02	7	1030401	6	1057-00	6	11303	53
1006-01	7	1040025	42	1057-01	6	11304	53
1006-02	7	1040115	42	1057-02	6	11305	53
1006-03	7	1040127	61	1057-03	6	11306	53
1006-04	7	1040159	28	1057-04	6	11311000	10
1006-05	7	1040167	61	1057-05	6	11312000	10
1010-01	7	1040197	61	1060025	42	11313000	10
1010-02	7	1040227	61	1060115	42	1140917	58
1010025	42	1040247	61	1061-00	6	11421000	10
10101	32	1040251	6	1061-01	6	11422000	10
1010115	42	1040301	6	1061-02	6	11423000	10
1010127	61	1040351	6	1061-03	6	11431000	10
1010159	28	1040401	6	1061-04	6	11432000	10
1010167	61	1044-30	9	1061-05	6	11433000	10
1010197	61	1044-31	9	1070025	42	11441000	10
10102	32	1044-40	9	1070115	42	11442000	10
1010227	61	1044-41	9	1080025	42	11443000	10
1010247	61	1044-80	9	1080115	42	11451000	10
1010251	6	1044-81	9	1090025	42	11452000	10
10103	32	1044-82	9	1090115	42	11461000	10
1010301	6	1044-83	9	1100025	42	11462000	10
1010351	6	1044-84	9	11001	51	11471000	10
10104	32	1044-85	9	1100115	42	11472000	10
1010401	6	1044-90	9	11002	51	11481000	10
10105	32	1050-00	6	11003	51	11481001	10
10106	32	1050-01	6	11004	51	11481002	10
1018-01	8	1050-02	6	11101	49	11481003	10
1018-02	8	1050-03	6	11102	49	11482000	10
1018-03	8	1050-04	6	11103	49	11482001	10
1018-04	8	1050-10	6	11201	50	11482002	10
1018-05	8	1050-11	6	11201000	9	11482003	10
1018-10	8	1050-12	6	11202	50	11483000	10
1018-11	8	1050-13	6	11202000	9	11483001	10
1018-12	8	1050-14	6	11203	50	11483002	10

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
11483003	10	13500006	14	13983000	20	15540005	77
1160917	58	13500007	14	13984000	20	15540006	77
11631000	29	13500008	14	13985000	20	15540008	77
11632000	29	13500009	14	13991000	20	15540067	77
11640005	29	13501004	14	13992000	20	15550002	77
11640006	29	13501005	14	13993000	20	15550005	77
11640007	29	13501006	14	13994000	20	15550006	77
11640008	29	13501007	14	13995000	20	15550008	77
11650005	29	13501008	14	14801004	67	15550067	77
11650006	29	13501009	14	14801005	67	15580002	77
11650007	29	13590000	21	14801006	67	15580005	77
11701	53	13600000	21	14801007	67	15580006	77
11702	53	13630000	21	14801008	67	15580008	77
11703	53	13640000	21	14801009	67	15580067	77
11704	53	13650000	21	14801010	67	15590002	77
11705	53	13651056	20	14970005	62	15590005	77
11706	53	13652067	20	14970006	62	15590006	77
1190917	58	13653008	20	14970007	62	15590008	77
12101	49	13654009	20	155-01	14	15590067	77
12102	49	13655010	20	155-02	14	15600004	73
12103	49	13660000	21	155-03	14	15600005	73
12201	49	13661056	20	155-04	14	15600006	73
12202	49	13662067	20	155-05	14	15600007	73
12203	49	13663008	20	155-06	14	15600008	73
1220917	58	13664009	20	155-07	14	15600009	73
12301	50	13665010	20	155-08	14	15780004	76
12302	50	13740300	22	15500031	76	15780005	76
12303	50	13740400	22	15500032	76	15780006	76
1240917	58	13749300	22	15500033	76	15780008	76
12840000	37	13749400	22	15500034	76	15780067	76
12870000	37	13780300	21	15500035	76	15800005	76
12982105	42	13780400	21	15500036	76	15800006	76
12982106	42	13790309	21	15500037	76	15800008	76
12982107	42	13790409	21	15500038	76	15800067	76
12982108	42	13800005	17	15500039	76	19401	69
12982205	42	13800006	17	15500040	76	19402	69
12982206	42	13800007	17	15500041	76	19403	69
12982207	42	13800008	17	15500042	76	19404	69
12982208	42	13810000	17	15500043	76	19405	69
12991102	42	13850000	17	15500044	76	19406	69
12991202	42	13850005	17	15500045	76	2004-00	82
12992102	42	13850007	17	15500046	76	2004-01	82
12992202	42	13850008	17	15500047	76	2004-02	82
12993105	42	13850009	17	15500048	76	2004-03	82
12993106	42	13850010	17	15500049	76	2004-04	82
12993107	42	13850011	17	15500050	76	2004-05	82
12993108	42	13850056	17	15520002	77	2006-01	80
12993205	42	13850067	17	15520005	77	2006-02	80
12993206	42	13860900	17	15520006	77	2006-03	80
12993207	42	13860900	18	15520008	77	2006-04	80
12993208	42	13871000	20	15520067	77	2006-05	80
13400005	14	13872900	20	15530002	77	2006-06	80
13400006	14	13950000	20	15530005	77	2006-07	80
13400007	14	13960000	20	15530006	77	2006-08	80
13400008	14	13970000	20	15530008	77	2006-09	80
13500004	14	13981000	20	15530067	77	2006-10	80
13500005	14	13982000	20	15540002	77	2007-01	80

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
2007-02	80	2044-04	79	350SL	47	4462-02	61
2007-03	80	2044-07	79	350SR	47	4462-03	61
2007-04	80	2044-08	79	350SXL	47	4462-04	61
2007-05	80	2044-09	79	350XLL	47	4462-05	61
2007-06	80	3001-01	24	350XLR	47	4462-06	61
2007-07	80	3001-02	24	350XSR	47	4463-02	61
2007-08	80	3001-03	24	350XLL	47	4463-03	61
2007-09	80	3001-04	24	350XXLR	47	4463-04	61
2007-10	80	3001-05	24	351LL	44	4463-05	61
2020-01	80	3001-06	24	351LR	44	4463-06	61
2020-02	80	3002-01	24	351ML	44	4464-02	61
2020-03	80	3002-02	24	351MR	44	4464-03	61
2020-04	80	3002-03	24	351SL	44	4464-04	61
2020-11	80	3002-04	24	351SR	44	4464-05	61
2020-12	80	3002-05	24	351XLL	44	5005-01	43
2020-13	80	3002-06	24	351XLR	44	5005-02	43
2020-14	80	3004-01	24	351XSL	44	5005-03	43
2021-01	80	3004-02	24	351XSR	44	5005-04	43
2021-02	80	3004-03	24	351XXLL	44	5015-01	46
2021-03	80	3004-04	24	351XXLR	44	5015-06	46
2021-04	80	3004-05	24	4004-01	74	5016-15	46
2021-11	80	3004-06	24	4004-02	74	5016-20	46
2021-12	80	3004-11	24	4004-03	74	5019-01	46
2021-13	80	3004-12	24	4004-04	74	5019-02	46
2021-14	80	3004-13	24	4005-01	75	5019-03	46
2027-01	83	3004-14	24	4005-02	75	5019-04	46
2027-02	83	3004-15	24	4005-03	75	5023-01	45
2027-03	83	3004-16	24	4005-04	75	5023-02	45
2027-04	83	3015-01	24	4016-01	74	5023-03	45
2027-05	83	3015-02	24	4016-02	74	5023-04	45
2033-01	82	3015-03	24	4016-03	74	5023-05	45
2033-02	82	3015-04	24	4016-04	74	5023-06	45
2033-03	82	3015-05	24	4300	88	5023-07	45
2033-04	82	3015-11	24	4300/05AB	88	5023-08	45
2035-01	82	3015-12	24	4300/05WK	88	5023-09	45
2035-02	82	3015-13	24	4300K	88	5023-10	45
2035-03	82	3015-14	24	4301	88	5026-00	48
2035-04	82	3015-15	24	4302	88	5040-01	44
2035-05	82	3302D39	51	4305	88	5040-02	44
2039-00	83	3302D59	51	4305K	88	5040-03	44
2039-01	83	3302D79	51	4306	88	5040-04	44
2039-02	83	340M/LL	48	4307	88	5066-80	43
2039-03	83	340M/LR	48	4433-01	60	5066-81	43
2039-04	83	340S/ML	48	4433-02	60	5066-82	43
2039-05	83	340S/MR	48	4433-03	60	5066-83	43
2039-98	83	343M/LL	48	4433-04	60	5067-01	45
2039-99	83	343M/LR	48	4433-05	60	5067-02	45
2041-00	80	343S/ML	48	4433-10	60	5068-00	45
2041-01	80	343S/MR	48	4433-11	60	5069-00	47
2041-02	80	345M/LL	48	4433-22	60	5073-01	44
2041-03	80	345M/LR	48	4433-33	60	5073-02	44
2041-07	80	345S/ML	48	4442-00	27	5073-03	44
2041-08	80	345S/MR	48	4442-10	27	5073-04	44
2041-09	80	350LL	47	4461-02	61	5073-05	44
2044-01	79	350LR	47	4461-03	61	5073-06	44
2044-02	79	350ML	47	4461-04	61	5073-07	44
2044-03	79	350MR	47	4461-05	61	5073-08	44

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
5073-09	44	7066-24	60	8017-00	33	9007-00	27
5073-10	44	7066-26	60	8017-01	33	9008-00	28
5073-11	44	7066-28	60	8017-02	33	9008-01	28
5073-12	44	7071-01	59	8017-03	33	9008-02	28
5074-01	44	7071-02	59	8017-04	33	9008-03	28
5074-02	44	7071-03	59	8017-51	33	9008-04	28
5074-03	44	7071-05	59	8017-52	33	9008-05	28
5074-04	44	7071-07	59	8017-53	33	9008-06	28
5074-05	44	7071-11	59	8017-54	33	9011-00	26
5074-06	44	7071-12	59	8020-00	34	9012-00	26
5074-07	44	7071-13	59	8020-01	34	9014-00	27
5074-08	44	7071-15	59	8023-01	33	9023-00	26
5074-09	44	7071-17	59	8023-02	33	9023-01	26
5074-10	44	7071-55	59	8023-03	33	9023-02	26
6001-01	36	7071-57	59	8023-04	33	9023-03	26
6001-02	36	7071-59	59	8023-05	33	9023-04	26
6001-03	36	7072-02	59	8105-01	69	9100-00	50
6003-01	39	7072-03	59	8105-02	69	9105-07	53
6003-02	39	7072-04	59	8105-03	69	9105-08	53
6003-03	39	7072-05	59	8105-04	69	9105-20	53
6003-04	39	7091-03	58	8105-05	69	9105-21	53
6004-01	38	7091-04	58	8105-06	69	9105-22	53
6004-02	38	7091-05	58	8106-01	68	9105-24	53
6004-03	38	7091-06	58	8106-03	68	9105-25	53
6004-04	38	7091-07	58	8106-04	68	9105-26	53
6005-00	37	7091-08	58	8106-05	68	9111-07	49
6101-01	11	7091-50	58	8106-06	68	9115-05	51
6101-02	11	7091-51	58	8120-00	69	9115-06	51
6101-03	11	7091-52	58	8121-00	68	9118-01	55
6102-00	11	7091-55	58	8130-00	68	9118-02	55
6104-00	11	7091-56	58	8140-00	68	9118-03	55
6105-00	11	7091-57	58	8726-50	103	9119-01	52
7041-01	60	7091-58	58	8761-01	74	9119-02	52
7041-02	60	7102-00	12	8761-02	74	9119-03	52
7041-03	60	7102-22	12	8761-03	74	9119-04	52
7041-04	60	7116-01	12	8761-04	74	9119-05	52
7041-05	60	7117-02	12	8788-03	35	9119-06	52
7041-11	60	8001-01	32	8788-50	34	9119-07	52
7041-22	60	8001-02	32	8829-01	99	9119-08	52
7041-33	60	8001-03	32	9001-01	25	9119-09	52
7041-44	60	8001-04	32	9001-02	25	9120-00	52
7041-55	60	8001-05	32	9001-03	25	9121-01	52
7042-01	66	8004-01	32	9001-04	25	9121-02	52
7042-02	66	8004-02	32	9001-05	25	9121-03	52
7042-03	66	8004-03	32	9002-01	25	9121-04	52
7042-04	66	8004-04	32	9002-02	25	9121-05	52
7042-05	66	8004-05	32	9002-03	25	9121-06	52
7043-01	66	8004-06	32	9002-04	25	9121-07	52
7043-02	66	8004-11	32	9002-05	25	9121-08	52
7043-03	66	8004-12	32	9003-01	25	9122-00	52
7043-04	66	8004-13	32	9003-02	25	9125-01	51
7043-05	66	8004-14	32	9003-03	25	9125-02	51
7066-14	60	8004-15	32	9003-04	25	9125-03	51
7066-16	60	8004-16	32	9003-05	25	9125-04	51
7066-18	60	8011-01	35	9006-00	27	9127-01	55
7066-20	60	8011-02	35	9006-10	27	9127-02	55
7066-22	60	8011-03	35	9006-50	27	9127-03	55

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
9128-01	55	A115006	8	A125106	47	A155006	69
9128-02	55	A115007	8	A125107	47	A155007	69
9128-03	55	A115008	8	A125108	47	A155008	69
9133-01	54	A116000	8	A125109	47	A9916-11	8
9133-03	54	A116005	8	A125204	47	A9916-12	8
9133-04	54	A116006	8	A125205	47	A9916-13	8
9133-06	54	A116007	8	A125206	47	A9916-14	8
9134-01	54	A116008	8	A125207	47	A9917-00	8
9134-02	54	A117000	8	A125208	47	A9917-01	8
9135-01	54	A117005	8	A125209	47	A9918-11	8
9135-02	54	A117006	8	A131056	21	A9918-12	8
9136-01	54	A117007	8	A131067	21	A9918-13	8
9136-02	54	A117008	8	A132056	21	A9918-14	8
9148-01	54	A118000	8	A132067	21	A9919-11	8
9148-02	54	A119901	33	A133000	21	A9919-12	8
9148-03	54	A119902	33	A140000	61	A9919-13	8
9148-04	54	A119903	33	A140001	61	A9919-14	8
9148-05	54	A119904	33	A141000	59	A9920-11	8
9148-06	54	A119905	33	A142000	59	A9920-12	8
9150-01	52	A121105	46	A142005	62	A9920-13	8
9150-02	52	A121106	46	A142006	62	A9920-14	8
9221-00	99	A121107	46	A142007	62	AB2000-10	73
9222-00	98	A121108	46	A142008	62	AB2000-12	73
9307-01	60	A121205	46	A142009	62	AB2000-14	73
9307-02	60	A121206	46	A143000	59	AB2000-16	73
9307-03	60	A121207	46	A143005	62	AB2340-00	72
9307-04	60	A121208	46	A143006	62	AB2340-06	72
9307-05	60	A122105	46	A143007	62	AB2341-00	72
9307-06	60	A122106	46	A143008	62	AB2341-06	72
9307-07	60	A122107	46	A143009	62	AB2360-00	72
9380-20	103	A122108	46	A144000	59	AB2361-00	72
9381-00	103	A122205	46	A144005	62	AB2362-00	72
9381-10	103	A122206	46	A144006	62	AB2900-10	73
9383-00	103	A122207	46	A144007	62	AB2900-12	73
9383-10	103	A122208	46	A144008	62	AB2900-14	73
9384-00	103	A123104	42	A144009	62	AB2900-16	73
9385-00	103	A123105	42	A145000	59	AB2900-18	73
9386-00	103	A123106	42	A146000	59	AB2900-20	73
9387-00	103	A123107	42	A146005	62	AF1000-01	75
A110001	26	A123108	42	A146006	62	AF1000-02	75
A111004	28	A123204	42	A146007	62	AF1000-03	75
A111005	28	A123205	42	A146008	62	AF1000-04	75
A111006	28	A123206	42	A146009	62	AF1001-01	75
A111007	28	A123207	42	A147000	61	AF1001-02	75
A111008	28	A123208	42	A148000	61	AF1001-03	75
A111009	28	A124104	42	A149000	61	AF1001-04	75
A111012	26	A124105	42	A150000	69	AF1001-05	75
A111015	26	A124106	42	A152000	72	AF1001-06	75
A111016	26	A124107	42	A152005	72	AF1001-07	75
A111017	26	A124108	42	A152305	79	AF1001-08	75
A111018	26	A124204	42	A152306	79	BF1120917	58
A112000	34	A124205	42	A152307	79	BF1160917	58
A112050	34	A124206	42	A152308	79	BF1190917	58
A113000	8	A124207	42	A152405	79	BF1220917	58
A113002	8	A124208	42	A152406	79	BF1240917	58
A115000	8	A125104	47	A152407	79	BF4433-12	60
A115005	8	A125105	47	A155005	69	BF4433-16	60



## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
BF4433-19	60	DEB6	100	M0501D01	89	M60-034	93
BF4433-22	60	DEBV2	100	M0501D05	89	M60-045	94
BF4433-24	60	DEBV3	100	M10-070-C	99	M60-076-I	92
BF4433-26	60	DEBV4	100	M10-080	98	M60-076-L	92
BF5066-80	43	DEBV6	100	M1133	96	M60-076-M	92
BF5066-81	43	DS1012-00	79	M1133-38	96	M60-076-S	92
BF5066-82	43	DS1013-10	79	M1150	96	M60-076-XS	92
BF5066-83	43	DS1013-12	79	M1151	96	M60-078-L	92
BF7006-01	62	DS1013-14	79	M1152	96	M60-078-M	92
BF7006-02	62	DS1014-10	79	M13-055	98	M60-078-S	92
BF7006-03	62	DS1014-12	79	M13-100	97	M60-224	93
BF7006-04	62	DS1014-14	79	M13-500	97	M60-300	89
BF7006-05	62	DS1014-16	79	M190FF	90	M70-100B	88
BF7006-06	62	DS1015-10	81	M190FF-24	90	M70-100L	88
BF7007-01	62	DS1015-12	81	M190FFE	90	M70-100S	88
BF7007-02	62	DS1015-14	81	M2035	97	M7035-L	97
BF7007-03	62	DS1016-10	81	M3000D08	98	M7035-M	97
BF7007-04	62	DS1016-12	81	M3000U	91	M7035-S	97
BF7007-05	62	DS1016-14	81	M3000UB	91	M7035-SB	97
BF7007-06	62	DS1016-16	81	M3001U	91	M7035-XS	97
BF7008-01	62	EB2	100	M3002D08	98	NE7702-01	65
BF7008-02	62	EB2040-00	36	M3002D09	98	NE7702-02	65
BF7008-03	62	EB3	100	M3006	90	NE7702-03	65
BF7008-04	62	EB4	100	M3014D	90	NE7702-04	65
BF7008-05	62	EB6	100	M3014D-6	90	NE7702-05	65
BF7008-06	62	EB6000-01	36	M3015	90	NE7702-06	65
BF7072-01	59	EB6000-02	36	M3022	90	NE7702-07	65
BF7072-02	59	EBV2	100	M3022-L	90	NE7702-08	65
BF7072-03	59	EBV3	100	M3040	89	NE7702-71	65
BF7072-04	59	EBV4	100	M3042	89	NE7702-72	65
BF7072-05	59	EBV6	100	M4209D06	97	NE7702-73	65
BF7072-06	59	FD1	52	M5001-L	96	NE7702-74	65
BS1000-10	18	FD2	52	M5001-M	96	NE7702-75	65
BS1000-12	18	FD3	52	M5001-S	96	NE7702-76	65
BS1000-14	18	FF4061-00	67	M5001-XL	96	NE7702-77	65
BS1000-16	18	HC1000-01	85	M5001-XXL	96	NE7707-72	65
BS1000-18	18	HC1000-02	85	M5001-XXXL	96	NE7707-73	65
BS1000-20	18	HC2000-10	84	M60-015	95	NE7707-74	65
CA298	86	HC2000-12	84	M60-015A	94	NE7707-75	65
CR2000-LS	86	HC2000-14	84	M60-015B	94	NE7707-76	65
CR3000-LS	86	HEB4	100	M60-015CD	94	NE7707-77	65
CR4000-LS	86	HEB6	100	M60-018	94	NE7710-70	65
CSG53510-A	77	HS1005-08	81	M60-024-L	92	NE7710-71	65
CSG53510-B	77	HS1005-10	81	M60-024-M	92	NE7717-71	66
CSG53510-C	77	HS1005-12	81	M60-024-S	92	NE7717-72	66
CSG53510-D	77	HS1005-14	81	M60-025-L	92	NE7717-73	66
CSG53520-A	77	HS1005-16	81	M60-025-M	92	NE7717-74	66
CSG53520-B	77	KB2040-00	64	M60-025-S	92	NE7717-75	66
CSG53520-C	77	KB2040-10	64	M60-031-L	93	NE7717-76	66
CSG53520-D	77	KB2040-40	64	M60-031-M	93	NE7717-77	66
CSG53550-A	77	KB3001-00	64	M60-031-S	93	NE7718-71	66
CSG53550-B	77	KB5000-00	63	M60-032	95	NE7718-72	66
CSG53550-C	77	KB5000-01	63	M60-032B	95	NE7718-73	66
CSG53550-D	77	KB5000-02	63	M60-032C	95	NE7718-74	66
DEB2	100	KB6000-00	64	M60-032D	95	NE7718-75	66
DEB3	100	KB6000-01	64	M60-032E	95	NE7718-76	66
DEB4	100	KB9000-01	63	M60-033	95	NE7718-77	66

## INDEX BY PRODUCT#

Product #	Page#	Product #	Page#	Product #	Page#	Product #	Page#
NE7727-72	39	T5035NS-L-10	102	TX4442-03	27	UPLP2630-15	15
NE7727-73	39	T5035ST	102	TX4442-04	27	UPLP2630-25	15
NE7727-74	39	T5035ST-10	102	TX8020-00	34	UPLP2730-15	15
NE7727-75	39	T5035ST-L	102	TX8023-21	33	UPLP2730-25	15
NE7727-76	39	T5035ST-L-10	102	TX8023-22	33	UPLT2333-15	16
NE7728-72	39	T5036NS	102	TX8023-23	33	UPLT2333-25	16
NE7728-73	39	T5036NS-10	102	TX8023-24	33	UPLT2433-15	16
NE7728-74	39	T5036NS-L	102	TX9902-01	33	UPLT2433-25	16
NE7728-75	39	T5036NS-L-10	102	TX9902-02	33	UPLT2533-15	16
NE7728-76	39	T5036ST	102	TX9902-02/100	33	UPLT2533-25	16
NE7730-70	38	T5036ST-10	102	TX9902-03	33	UPLT2633-15	16
NE7731-72	38	T5036ST-L	102	TX9902-03/100	33	UPLT2633-25	16
NE7731-73	38	T5036ST-L-10	102	TX9902-04	33	UPLT2733-15	16
NE7731-74	38	T5037NS	102	TX9902-04/100	33	UPLT2733-25	16
NE7731-75	38	T5037NS-10	102	TX9902-05	33	UPMP1330-15	15
NE7738-70	48	T5037NS-L	102	TX9902-05/100	33	UPMP1430-15	15
NE7740-70	46	T5037NS-L-10	102	TX9902-06	33	UPMP1530-15	15
NE7741-02	18	T5037ST	102	TX9902-06/100	33	UPMP1630-15	15
NE7741-03	18	T5037ST-10	102	UDLE	16	UPMP1730-15	15
NE7741-04	18	T5037ST-L	102	UEP2031	16	UPT1333-15	16
NE7741-05	18	T5037ST-L-10	102	UEP2051	16	UPT1333-25	16
NE7741-06	18	T5055NS	102	UPL1330-0	15	UPT1433-15	16
NE7741-07	18	T5055NS-10	102	UPL1330-15	15	UPT1433-25	16
NE7741-12	18	T5055ST	102	UPL1330-25	15	UPT1533-15	16
NE7741-13	18	T5055ST-10	102	UPL1430-0	15	UPT1533-25	16
NE7741-14	18	T5075NS	102	UPL1430-15	15	UPT1633-15	16
NE7741-15	18	T5075NS-10	102	UPL1430-25	15	UPT1633-25	16
NE7741-16	18	T5075ST	102	UPL1530-0	15	UPT1733-15	16
NE7741-17	18	T5075ST-10	102	UPL1530-15	15	UPT1733-25	16
NS2000-10	84	T5085NS	102	UPL1530-25	15	UPT2333-15	16
NS2000-12	84	T5085NS-10	102	UPL1630-0	15	UPT2333-25	16
NS2000-14	84	T5085ST	102	UPL1630-15	15	UPT2433-15	16
NS2000-20	84	T5085ST-10	102	UPL1630-25	15	UPT2433-25	16
NS2000-21	84	T6100	102	UPL1730-0	15	UPT2533-15	16
SC1000-00	72	T6110	102	UPL1730-15	15	UPT2533-25	16
SC2000-10	84	T6120	102	UPL1730-25	15	UPT2633-15	16
SC2000-12	84	T6160	102	UPL2330-0	15	UPT2633-25	16
SC2000-14	84	T700	101	UPL2330-15	15	UPT2733-15	16
SL8001-01	78	T725NS	101	UPL2330-25	15	UPT2733-25	16
SL8001-03	78	T725ST	101	UPL2430-0	15	USPK	16
SL8001-05	78	T735NS	101	UPL2430-15	15	WLK02	86
SL8001-07	78	T735ST	101	UPL2430-25	15	WLK02-W	86
SL8001-09	78	T736NS	101	UPL2530-0	15		
SL8001-33	78	T736ST	101	UPL2530-15	15		
SL8001-35	78	T737NS	101	UPL2530-25	15		
SL8001-37	78	T737ST	101	UPL2630-0	15		
T5015NS	102	T800	101	UPL2630-15	15		
T5015NS-10	102	T825NS	101	UPL2630-25	15		
T5015ST	102	T825ST	101	UPL2730-0	15		
T5015ST-10	102	T835NS	101	UPL2730-15	15		
T5025NS	102	T835ST	101	UPL2730-25	15		
T5025NS-10	102	T836NS	101	UPLP2330-15	15		
T5025ST	102	T836ST	101	UPLP2330-25	15		
T5025ST-10	102	T837NS	101	UPLP2430-15	15		
T5035NS	102	T837ST	101	UPLP2430-25	15		
T5035NS-10	102	TX4442-00	27	UPLP2530-15	15		
T5035NS-L	102	TX4442-02	27	UPLP2530-25	15		

**Trademarks**

ATC™ is a trademark of DeRoyal Industries, Inc.  
Boa® is a registered trademark of Boa Technology, Inc.  
Capital™ is a trademark of DeRoyal Industries, Inc.  
Confor® is a registered trademark of Cabot Safety Intermediate LLC  
CoolJet™ is a trademark of DeRoyal Industries, Inc.  
Darco® is a registered trademark of Darco International, Inc.  
DeRoyal® is a registered trademark of DeRoyal Industries, Inc.  
Gelbo® is a registered trademark of Knit-Rite, Inc.  
JetStream® is a registered trademark of DeRoyal Industries, Inc.  
Hill-Rom® is a registered trademark of Hill-Rom Services, Inc.  
Lycra® is a registered trademark of Invista North America S.A.R.L.  
MedSurg™ is a trademark of Darco International, Inc.  
Metashank™ is a trademark of Darco International, Inc.  
Pacesetter™ is a trademark of DeRoyal Industries, Inc.  
Plastazote® is a registered trademark of Zotefoams PLC  
PneuGel® is a registered trademark of Remington Products Company  
PROlign® is a registered trademark of DeRoyal Industries, Inc.  
PRUventor™ is a trademark of DeRoyal Industries, Inc.  
Orthowedge™ is a trademark of Darco International, Inc.  
Solace® is a registered trademark of DeRoyal Industries, Inc.  
ULTRAlign® is a registered trademark of DeRoyal Industries, Inc.  
Warrior® is a registered trademark of DeRoyal Industries, Inc.  
Warrior™ is a trademark of DeRoyal Industries, Inc.





Improving Care. Improving Business.®

DeRoyal Industries, Inc.  
200 DeBusk Lane, Powell, TN 37849  
1-800-DeRoyal (inside U.S.)  
865-938-7828 (outside U.S.)  
[www.deroyal.com](http://www.deroyal.com)  
Reprint# 0-2098 | rev. 01/16

