

Signature
Product Line



 DeRoyal®

Improving Care. Improving Business.®



ORTHOPEDICS
AND SPORTS
MEDICINE



SIGNATURE
PRODUCT LINE



DeRoyal®

Improving Care. Improving Business.®



DeRoyal Industries, Inc. | 200 DeBusk Lane, Powell, TN 37849 | 800.DEROYAL or 865.938.7828 | www.deroyal.com

REV 1/16 REPRINT # 0-1678

■ Head/Neck/Spine	
Collars	4
Universal	6
SO	8
LSO	9
TLSO/OSO	11
■ Shoulder/Arm/Elbow/Wrist	
Abduction/Immobilizers	14
Arm Slings	16
Elbow Brace	18
Sleeves	20
Tennis Supports	21
Wrist	22
Wrist/Wrist/Forearm	24
Wrist/Thumb	25
■ Knee/Leg	
OA Bracing	28
ACL Bracing	29
Post-Op	30
Hinged Bracing	32
Patellafemoral	35
Supports/ Sleeves	37
■ Ankle/Foot	
Ankle Bracing	44
Ankle Stirrups	49
AFO	51
Walkers	52
■ Temperature Therapy	
Temperature Therapy Units	58
Therapy Blankets/Wraps	60
Thermal Supports	62
Accessories	63
■ Ambulatory Aids	
Fracture Bracing	66
Post-Op/Cast Shoes	67
Miscellaneous Ambulatory Aids	70





Collars

Universal

SO

LSO

TLSO/OSO

Head/Neck/Spine



HCPCS Code: L0174

Capital™ Collar Enhanced powered by Boa® Technology

- Patented design featuring a micro-adjustable dial powered by Boa® Technology that allows collar height modification in order to fit most patients
- Includes compliance ring to ensure patient compliance
- Reinforced sternal padding for optimum support and comfort
- Designed to help restrict flexion, extension, lateral bending and axial rotation
- Shoulder, mandible and occipital contours for increased stability
- Sheer finished open-cell liner is specifically designed to help prevent skin breakdown and provide sustained comfort during long-term use
- Not made with natural rubber latex

Indications: Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

Size	Collar	Collar & Pad Set	Pad Set Only
Universal	11311000	11312000	11313000

*Capital™ is a trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.*



HCPCS Code: L0174

Capital™ Cervical Collar

- Designed to help restrict flexion, extension, lateral bending and axial rotation
- Reinforced sternal padding for optimum support and comfort
- Shoulder and mandible contours for increased stability
- Sheer finished open-cell liner is specifically designed to help prevent skin breakdown and provide sustained comfort during long-term use
- Replacement pad set available upon request
- Not made with natural rubber latex

Indications: Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

Size	Collar	Collar & Pad Set	Pad Set Only
Adult Super Short	11421000	11422000	11423000
Adult Stout	11431000	11432000	11433000
Adult X-Small	11441000	11442000	11443000
Adult Short	11451000	11452000	11423000
Adult Regular	11461000	11462000	11423000
Adult Large	11471000	11472000	11423000
Infant 0 - 6 Months	11481000	11482000	11483000
Infant 6 Months - 2 Years	11481001	11482001	11483001
Child 2 - 6 Years	11481002	11482002	11483002
Child 6 - 12 Years	11481003	11482003	11483003



HCPCS Code: L0172


XTW Extended Wear Collar

- Designed to help restrict flexion, extension, lateral bending and axial rotation
- High density polyethylene frame for maximum rigidity and support
- Shoulder and mandible contours for increased stability
- Anterior and posterior flex tabs for individual comfort
- Latex free, non-toxic, open-cell liner helps prevent skin breakdown
- Five adult and five child sizes to ensure best patient fit
- Replacement pad set available upon request

Indications: Cervical sprains and strains; Cervical spinal stenosis; Cervical disc displacement; Cervical arthritis; Whiplash injury; Cervical spondylosis and spondylolisthesis; Degenerative disc disease; Cervicalgia

Size	Collar	Collar & Pad Set
Adult X-Short	11201000	11202000
Adult Short	11211000	11212000
Adult Regular	11221000	11222000
Adult Tall	11231000	11232000
Adult X-Tall	11241000	11242000
Infant Short	11251000	11252000
Infant Regular	11261000	11262000
Child Tot	11271000	11272000
Child Short	11281000	11282000
Child Regular	11291000	11292000

Capital™ is a trademark of DeRoyal Industries, Inc.



Warrior™ Spine Universal Back Brace powered by Boa® Technology

Head/Neck/Spine



Warrior™ Spine 650
HCPCS Code: L0650

Warrior™ is a trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Warrior™ Spine 648

HCPCS Code: L0648

- Universal size (24" to 60" waist)
- Only three (3) SKUs (159-01, 159-02, 159-03)
- Patent pending, single pull bi-lateral sizing adjustment
- Powered by Boa® Technology for a powerful and secure closure that is easy for the patient to tighten with the turn of a dial
- Enlarged Boa® dials allow micro-adjustability and on-the-fly compression and adjustment with fewer turns
- Integrated anterior panels and closure hand strap for easy donning and doffing
- Soft breathable material to provide a cool, comfortable fit
- Rounded anterior closure to accommodate varying patient anatomies (inferior and superior angulation)
- Hot/cold temperature therapy packs included with each device

Indications: Post-operative lumbar fusion; Post-operative lumbar immobilization; Laminectomy or discectomy; Lumbar compression fractures; Degenerative disc disease; Osteoporosis; Spondylolisthesis; Spondylolysis; Spinal stenosis; Single column instability immobilization; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; Spinal blocks

Size	Description	Belt Adjustment Length	Product #
Universal	Warrior™ Spine 642	24" to 60"	159-01
Universal	Warrior™ Spine 648	24" to 60"	159-02
Universal	Warrior™ Spine 650	24" to 60"	159-03



Warrior™ Spine 642

HCPCS Code: L0642



Head/Neck/Spine

Warrior™ is a trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.



HCPCS Code: L0642

Solace® Low Profile Back Brace

- Soft, breathable material provides cool, comfortable fit increasing patient compliance
- Adjustable pulley style closure system is easy to secure and provides maximum compression
- Anterior and posterior plastic inserts provide control and support
- Low profile, lightweight design fits easily under patient's clothing
- Plastic D-rings for ease of application and removal

Indications: Chronic low back pain; Lumbar muscle weakness; Lumbar sprains and strains; Mechanical or discogenic lumbar pain

Size	Waist Circumference	Product #
Small	26" - 34"	13400005
Medium	32" - 40"	13400006
Large	38" - 46"	13400007
X-Large	44" - 51"	13400008



HCPCS Code: L0627 or L0642

PROlign® Lumbar Orthosis powered by Boa® Technology

- Powered by The Boa® Closure System for a powerful and secure closure that is easy for the patient to tighten with the turn of a dial
- Dial allows micro-adjustability and on-the-fly adjustment
- Soft breathable material to provide a cool, comfortable fit
- Patented dual closure for adjustability
- Rigid anterior insert provides support to the intra-abdominal and viscera region helping to unload the spine under gravitational load
- Features posterior plastic inserts to help maintain neutral sagittal alignment while maximizing support of anatomical structures

Indications: Post-operative lumbar immobilization; Chronic low back pain; Lumbar muscle weakness; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; IDET procedure; Spinal blocks

Size	Waist Circumference	15°	25°
X-Small	25" - 34"	13502004	13503004
Small	30" - 34"	13502005	13503005
Medium	35" - 40"	13502006	13503006
Large	40" - 46"	13502007	13503007
X-Large	45" - 50"	13502008	13503008
XX-Large	50" - 55"	13502009	13503009

PROlign® & Solace® are registered trademarks of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroval.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





PROlign® EXT Lumbar Orthosis powered by Boa® Technology



HCPCS Code: L0631 or L0648

- Powered by The Boa® Closure System for a powerful and secure closure that is easy for the patient to tighten with the turn of a dial
- Dial allows micro-adjustability and on-the-fly adjustment
- Soft breathable material to provide a cool, comfortable fit
- Patented dual closure for adjustability
- Rigid anterior insert provides support to the intra-abdominal and viscera region helping to unload the spine under gravitational load
- Features posterior plastic inserts to help maintain neutral sagittal alignment while maximizing support of anatomical structures
- PROlign® EXT features extended posterior panel for coverage from sacrococcygeal junction to the T9 vertebrae

Indications: Post-operative lumbar immobilization; Chronic low back pain; Lumbar muscle weakness; Lumbar strains and sprains; Mechanical or discogenic disc disease; Postural support; Degenerative disc disease; IDET procedure; Spinal blocks

Size	Waist Circumference	15°	25°
X-Small	25" - 34"	158-01	157-01
Small	30" - 34"	158-02	157-02
Medium	35" - 40"	158-03	157-03
Large	40" - 46"	158-04	157-04
X-Large	45" - 50"	158-05	157-05
XX-Large	50" - 55"	158-06	157-06

PROlign® is a registered trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.



HCPCS Code:
L0637 or L0650



ULTRAlign® + Lumbar Orthosis powered by Boa® Technology

- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 0°, 15° or 25° lordosis
- Non-tapered (male) and tapered (female) versions available for best patient fit
- Mid and low profile version for short waisted patients
- LSO meets requirements for most T6 to L5 pathologies

Indications: Post-operative thoracic/lumbar fusion; Laminectomy or discectomy; Compression fractures; Degenerative disc disease; Osteoporosis; Single column instability immobilization; Facet syndrome

Size	15° Non-Tapered	15° Tapered	15° Low-Profile Tapered
Small	UPL3330-15	UPL4330-15	UPLP4330-15
Medium	UPL3430-15	UPL4430-15	UPLP4430-15
Large	UPL3530-15	UPL4530-15	UPLP4530-15
X-Large	UPL3630-15	UPL4630-15	UPLP4630-15
XX-Large	UPL3730-15	UPL4730-15	UPLP4730-15

Sizing Guide

Sizing determined by three circumferential measurements:

- 1" below tip of breast bone or xyphoid process
- Around navel or waist
- Around fullest area of trochanter or hips

Non-Tapered	Xyphoid	Waist	Trochanter
Small	31" - 35"	31" - 35"	31" - 35"
Medium	34" - 38"	32" - 38"	34" - 39"
Large	37" - 41"	35" - 41"	38" - 42"
X-Large	38" - 45"	38" - 45"	38" - 45"
XX-Large	45" - 53"	45" - 53"	45" - 53"

Tapered	Xyphoid	Waist	Trochanter
Small	25" - 31"	25" - 31"	27" - 33"
Medium	28" - 35"	29" - 36"	32" - 38"
Large	32" - 38"	32" - 38"	35" - 41"
X-Large	35" - 42"	36" - 42"	38" - 46"
XX-Large	40" - 46"	39" - 47"	43" - 49"

ULTRAlign® is a registered trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)



ULTRAlign® + Thoracic Lumbar Orthosis powered by Boa® Technology



HCPCS Code: L0454

- Constructed of hook receivable material with soft, breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closure system allows infinite adjustability and maximum compression
- Combines heat-moldable plastic inserts (anterior, posterior, lateral) for rigid support
- Available in 15° or 25° lordosis
- Non-tapered (male) and tapered (female) versions available for best patient fit
- DLE (Dorsal Lumbar Extension) or SPK (Sternal Pad Kit) available for maximum control
- Low profile version for short waisted patients
- Thoracic Lumbar Sacral Orthosis (TLSO) meets requirements for most T6 to L5 pathologies

Indications: Post-operative thoracic/lumbar fusion; Laminectomy or discectomy; Compression fractures; Degenerative disc disease; Osteoporosis; Single column instability immobilization; Facet syndrome

Size	15° Tapered	15° Non-Tapered
Small	UPT4333-15	UPT3333-15
Medium	UPT4433-15	UPT3433-15
Large	UPT4533-15	UPT3533-15
X-Large	UPT4633-15	UPT3633-15
XX-Large	UPT4733-15	UPT3733-15



HCPCS Code: L0456

OSO™ Osteoporosis Spinal Orthosis

- Patent pending Noble Hinge™ technology for patient specific fit, rigid support and progressive adjustment
- Constructed of a hook receivable material with soft breathable liner that provides a cooler, lighter fit
- Patented dual dynamic closures allow adjustment and maximum compression
- Additional flexible posterior panel extends from the sacrococcygeal junction to the scapular spine, providing additional support and comfort
- Padded shoulder straps can be easily adjusted by the patient for increased comfort and compliance in one of two configurations

Indications: Osteoporosis; Acute and chronic pain due to osteoporosis; Hyperkyphosis; Compression fractures; Pain relieved by thoracic extension; Spinal stenosis; Post-operative support; Vertebral collapse

Size	Waist	Short <5'7"	Regular 5'5" - 6'6"	Tall > 6'
Small	26" - 34"	OPB3333-S	OPB3333-R	OPB3333-T
Medium	34" - 41"	OPB3433-S	OPB3433-R	OPB3433-T
Large	41" - 48"	OPB3533-S	OPB3533-R	OPB3533-T

ULTRAlign® is a registered trademark of DeRoyal Industries, Inc.
Noble Hinge™ & OSO™ are trademarks of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.





Abduction/Immobilizers

Arm Slings

Elbow Braces

Sleeves

Tennis Supports

Wrist

Wrist/Thumb

Wrist/Forearm

Shoulder/Arm Elbow/Wrist



HCPCS Code: L3960

Warrior™ Wing Shoulder Abduction Device

- Universal Sling Design = one size fits most
- Aluminum waist frame is easily formed for accurate positioning & comfort
- Airplane hinge for optimal abduction and shoulder recovery
- Allows both abduction, rotation, and internal or external fixation
- Includes three triangle pillows for increased abduction and/or fixation (90°, 45°, 30°, 15°)
- Easy enough to be self-applied

Indications: For non-surgical or post-op shoulder conditions that require diagnosis-specific arm positioning; Rotator cuff repairs; Bankart repairs; SLAP lesions; Glenohumeral dislocations/subluxation; Posterior/anterior capsular reconstruction; Shoulder instabilities; Muscle and tendon repair; Joint reconstruction

Size	Description	Product #
Universal	Left	11631000
Universal	Right	11632000

Warrior™ is a trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Shoulder P.A.D. III

- Unique design allows sling length to be adjusted to fit most patients
- Constructed of hook receivable material for universal adjustability
- Soft, breathable liner wicks away moisture
- Foam block provides 35° of abduction
- Shoulder and waist straps can be trimmed to desired length
- Quick-release buckles on shoulder and waist straps for ease of application and removal
- Adjustable thumb loop to prevent wrist drop
- Arm can be positioned more anteriorly or posteriorly depending on post-operative requirements
- Exercise ball packaged with product
- Universal right/left and size

Indications: Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

HCPCS Code: L3670

Size	Product #
Universal	11620000



Shoulder P.A.D. II

- Constructed of hook receivable material for adjustability
- Soft, breathable liner wicks away moisture
- Foam block provides 35° of abduction
- Shoulder and waist straps can be trimmed to desired length
- Adjustable for open or closed pouch design
- Quick-release buckle on waist strap for ease of application and removal
- Adjustable thumb loop to prevent wrist drop
- Arm can be positioned more anteriorly or posteriorly depending on post-operative requirements
- Exercise ball packaged with product
- Universal right/left

HCPCS Code: L3670

Indications: Surgical rotator cuff repairs; Anterior glenohumeral capsular reconstruction; Dislocation/subluxation of G-H joint; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Pouch Size	Product #
Small	7" x 16"	11630005
Medium	7" x 17"	11630006
Large	7" x 18"	11630007
X-Large	7" x 19"	11630008



HCPCS Code: L3670

Shoulder Abduction Positioner

- Constructed of hook receivable foam material for adjustability
- Foam block provides 55° of abduction
- Adjustable wrist and humeral cuffs for proper positioning
- Foam shoulder pad for additional patient comfort
- Designed for abduction and immobilization of the shoulder
- Universal right/left

Indications: Acute shoulder instability; Surgical rotator cuff repair; Anterior glenohumeral capsular reconstruction; Acute rotator cuff tendinopathy; Dislocation/subluxation of G-H joint; Acromioclavicular joint dislocations; S/P glenohumeral acromioplasty; Bankart repairs; Soft tissue repairs; Shoulder sprains and strains

Size	Block Dimensions	Product #
Small	25¼" x 23¾"	11650005
Medium	29½" x 26¾"	11650006
Large	30½" x 36"	11650007



Premium Arm Sling

- Deep pocket pouch design constructed of black cotton material
- Infant and child sizes feature cartoon print
- Foam shoulder pad provides additional patient comfort
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with O-rings

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Infant	5" x 9"	1" x 21"	11690001
Child	6" x 11"	1" x 27"	11690002
Small	7" x 13"	2" x 34"	11690005
Medium	8" x 17"	2" x 37"	11690006
Large	9" x 20"	2" x 37"	11690007
X-Large	9" x 23"	2" x 42"	11690008



Specialty Arm Sling

- Durable navy canvas construction with deep pocket design
- Infant and child sizes feature cartoon print
- Choice of foam or vinyl shoulder pad for patient comfort
- Thumb loop helps prevent wrist drop
- Adjustable slide buckle on strap for best fit
- Hook and loop closure with O-rings

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #	Vinyl Pad
Infant	5" x 9"	1" x 21"	8004-01	8004-11
Child	6" x 11"	1" x 27"	8004-02	8004-12
Small	7" x 13"	2" x 34"	8004-03	8004-13
Medium	8" x 17"	2" x 37"	8004-04	8004-14
Large	9" x 20"	2" x 37"	8004-05	8004-15
X-Large	9" x 23"	2" x 42"	8004-06	8004-16



Universal Arm Sling

- Constructed of heavy-duty canvas material
- Features easy fold-back universal design
- Foam shoulder pad provides additional comfort
- Hook and loop closure with O-rings

Indications: Cast support; Injuries to the arm, wrist or hand; Acromioclavicular sprains; Anterior glenohumeral dislocations

Size	Pouch Size	Strap Length	Product #
Adult	8" x 19"	2" x 37"	8020-00
Child	5" x 9"	1" x 31"	8020-01



HCPCS Code: L3760

Warrior™ Elbow Post-Operative Brace

- Universal sizing (left & right specific)
- Single hinge lateral application for easy donning and doffing
- Hinge is easily adjusted in 10-degree increments (Flexion: 0° to 120°, Extension: 0° to 90°)
- Length can be adjusted from 16" to 24"
- Shoulder strap included for extra support and immobilization


Indications: Post-operative rehabilitation; Olecranon fractures; Radial head fractures; Hairline fracture of the radius, humerus and ulna; Ligamentous reconstruction procedures of the elbow

Size	Description	Product #
Universal	Left	EB6000-01
Universal	Right	EB6000-02

Warrior™ is a trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroval.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





HCPCS Code: L3710

Shoulder/Arm/Elbow/Wrist

Hypercontrol® Elbow Brace

- Constructed of breathable Tri-tex material with four way stretch
- Medial and lateral hinges control hyperextension of the elbow
- Criss-cross straps assist in controlling hyperextension
- Elbow pad protects bony prominence
- Proximal and distal straps provide secure fit and help prevent migration

Indications: Medial/lateral ligament instabilities; Hyperextension of the elbow; UCL sprains; Valgus Extension Overload (VEO); Joint capsular sprains

Size	Forearm Circumference	Product #
X-Small	8" - 9"	12850004
Small	9" - 10"	12850005
Medium	10" - 11"	12850006
Large	11" - 12"	12850007
X-Large	12" - 13"	12850008
XX-Large	13" - 14"	12850009
XXX-Large	14" - 15"	12850010
XXXX-Large	15" - 16"	12850011

*Measure largest circumference of forearm

Hypercontrol® is a registered trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





HCPCS Code: L3760

Transition Elbow Brace

- Breathable foam softgoods are trimmable for a customized fit
- Hinge can be easily adjusted in 10-degree increments without tools (Flexion: 120° to 0° and Extension: 90° to 0°)
- Loop lock strapping is trimmable to allow for superior fit
- Length can be adjusted from 19" to 16"
- Universal right/left

Indications: Post-operative rehabilitation

Size	Measurement	Product #
Small	6" - 8"	EB5000-01
Medium	8" - 10"	EB5000-02
Large	10" +	EB5000-03



Neoprene Elbow Sleeve

- 1/8" neoprene construction provides warmth, compression, and support
- Anatomically contours to the elbow
- Available with or without elbow strap

Indications: Tendonitis; Osteoarthritis; Arthritis; Olecranon bursitis; (Medial and lateral epicondylitis w/ strap)

Size	Forearm Circumference	w/o Strap	w/ Strap
Small	9" - 10"	NE7727-72	NE7728-72
Medium	10" - 11"	NE7727-73	NE7728-73
Large	11" - 12"	NE7727-74	NE7728-74
X-Large	12" - 14"	NE7727-75	NE7728-75
XX-Large	14" - 16"	NE7727-76	NE7728-76



Premium Tennis Elbow Support

- Made of lightweight, ventilated material
- Support is low profile increasing patient compliance
- Inner pad can be strategically placed for maximum pain relief
- D-ring closure for easy application
- Universal one size fits all

Indications: Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	2" x 16"	12870000



PneuGel® Tennis Elbow Support

- Constructed of soft foam material for patient comfort
- Concentration of pressure allows compression to be applied where needed while not constricting uninvolved muscles or blood flow
- Air/gel removable bladder can be placed in freezer to provide cold therapy
- Inflation bulb included to adjust bladder
- Universal one size fits all

Indications: Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	2" x 20"	12840000



Universal Neoprene Tennis Elbow Support

- 1/8" neoprene construction provides compression, warmth and support
- Adjustable felt buttress can be placed just below the epicondyle for pain relief
- Loop lock contact closure for adjustable compression and application
- Universal one size fits most patients

Indications: Medial and lateral epicondylitis

Size	Measurement	Product #
Universal	3" x 16"	NE7730-70



HCPCS Code: L3908

Premium Universal Wrist Splint

- Constructed of comfortable foam material that helps keep patient's skin cool and dry
- Dynamic closure allows easy adjustment for an intimate fit
- Removable palmar and dorsal stays make product universal and allow for proper support
- Universal right/left and size fits most patients

Indications: Mild sprains and strains of the wrist; Tendonitis, Tenosynovitis; Carpal tunnel syndrome; Overuse injuries

Size	Length	Product #
Universal	8"	5068-00



HCPCS Code: L3908

Premium Wrist Splint

- Breathable material with four way stretch provides patient comfort
- Removable/malleable palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
X-Small	8"	351XSR	351XSL
Small	8"	351SR	351SL
Medium	8"	351MR	351ML
Large	8"	351LR	351LL
X-Large	8"	351XLR	351XLL
XX-Large	8"	351XXLR	351XXLL



HCPCS Code: L3908

Universal Black Wrist Splint

- Constructed of comfortable black foam material
- Malleable palmar and dorsal stays help provide anatomical fit and rigid support
- Movable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

Indications: Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist; Carpal tunnel syndrome

Size	Length	Right	Left
Universal	7"	5067-01	5067-02



HCPCS Code: L3908

Black Wrist Splint

- Constructed of comfortable black foam material
- Malleable palmar stay provides anatomical fit and rigid support
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
Pediatric	4"	5073-11	5073-12
X-Small	8"	5073-01	5073-06
Small	8"	5073-02	5073-07
Medium	8"	5073-03	5073-08
Large	8"	5073-04	5073-09
X-Large	8"	5073-05	5073-10



HCPCS Code: L3908

Lace Up Suede Leatherette Wrist Splint

- Soft perforated suede leatherette material with nylon tricot lining
- Removable/malleable palmar stay
- Single pull lace closure for compression and easy adjustment
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
X-Small	8"	5023-01	5023-06
Small	8"	5023-02	5023-07
Medium	8"	5023-03	5023-08
Large	8"	5023-04	5023-09
X-Large	8"	5023-05	5023-10



HCPCS Code: L3908

Black Wrist/Forearm Splint

- Constructed of comfortable black foam material
- Malleable palmar stay provides anatomical fit and rigid support
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: Mild sprains and strains of the wrist; Tendonitis; Carpal tunnel syndrome; Overuse injuries; Tenosynovitis

Size	Length	Right	Left
X-Small	10"	5074-01	5074-06
Small	10"	5074-02	5074-07
Medium	10"	5074-03	5074-08
Large	10"	5074-04	5074-09
X-Large	10"	5074-05	5074-10



HCPCS Code: L3908

Suede Leatherette Wrist & Wrist/Forearm Splint

- Suede leatherette material reduces migration
- Unique wedge contour design ensures an intimate fit increasing patient comfort
- Extra-wide loop lock strap at wrist for compression and support
- Removable/malleable palmar and dorsal stays
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: Carpal tunnel syndrome; Sprains and strains of the wrist; Post-operative immobilization and rehabilitation; Stable fractures of the distal radius and carpal joints

Size	Wrist Length	Right	Left	Wrist/Forearm Length	Right	Left
Pediatric	—	—	—	6"	12991202	12991102
X-Small	—	—	—	9"	12992202	12992102
Small	8"	12982205	12982105	11½"	12993205	12993105
Medium	8"	12982206	12982106	11½"	12993206	12993106
Large	8"	12982207	12982107	11½"	12993207	12993107
X-Large	8"	12982208	12982108	11½"	12993208	12993108

HCPCS Code: L3908



Premium Wrist & Wrist/Forearm Splint

- Constructed of a unique odor resistant material with moisture transport properties to keep patient's skin cool and dry
- Malleable palmar and dorsal stays provide anatomical fit and proper support
- Movable dorsal stay makes product universal in size
- D-ring closures for easy application and removal

Indications: Sprains and strains of the wrist; Tendonitis; Tenosynovitis of the wrist and distal forearm

Size	Length	Right	Left
Universal Wrist	8"	5066-92	5066-93
Universal Wrist/Forearm	10"	5066-90	5066-91



HCPCS Code: L3809

Premium Universal Wrist & Thumb Splint

- Constructed of comfortable foam material that helps keep patient's skin cool and dry
- Removable palmer and dorsal stays make product universal and allows for proper support
- Dynamic closure allows easy adjustment for intimate fit
- Adjustable thumb stay ensures desired thumb positioning
- Universal right/left and size fits most patients

Indications: De Quervain's syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

Size	Length	Product #
Universal	8"	5069-00



HCPCS Code: L3809

Premium Wrist & Thumb Splint

- Breathable material with four way stretch provides patient comfort
- Removable/adjustable thumb and palmar stay
- Low profile and lightweight design
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: De Quervain's syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

Size	Length	Right	Left
X-Small	8"	350XSR	350XSL
Small	8"	350SR	350SL
Medium	8"	350MR	350ML
Large	8"	350LR	350LL
X-Large	8"	350XLR	350XLL
XX-Large	8"	350XXLR	350XXLL



HCPCS Code: L3809


Black Foam Wrist & Thumb Splint

- Constructed of comfortable black foam material
- Removable/adjustable thumb and palmer stay
- D-ring closures for easy application and removal
- Does not extend past distal palmar crease allowing full finger range of motion

Indications: De Quervain's syndrome; Scaphoid injuries; Carpometacarpal instabilities; Metacarpal fractures; Post cast healing; Gamekeeper's thumb; Cumulative trauma injuries

Size	Length	Right	Left
X-Small	8"	A125204	A125104
Small	8"	A125205	A125105
Medium	8"	A125206	A125106
Large	8"	A125207	A125107
X-Large	8"	A125208	A125108
XX-Large	8"	A125209	A125109





- OA Bracing
- ACL Bracing
- Post-Op
- Hinged Bracing
- Patellafemoral
- Supports & Sleeves

Knee/Leg



HCPCS Code:
L1843 or K0901



OA Single Upright

- Lightweight, low profile, and easy-to-use design
- Unique varus/valgus angle adjustment to provide compartmental off-loading and optimize pain relief
- Lateral application avoids tibial crest pressure and medial clearance issues
- Maximized 3 points of leverage for medial compartment
- 4 SKUs promote inventory reduction and minimal stocking levels
- Interchangeable condyle pads for increased comfort, stability and pain control
- Quick release buckles for easy donning and doffing
- Silicone gripped liners and strap pads to help prevent migration
- Optional suspension wrap (included) for optimal fixation
- Optional flexion/extension stops for range of motion control
- Measurement Guide: Standard - Under 6'0" and Tall - 6'0" and above

Indications: Mild to moderate unicompartmental osteoarthritis

Size	Left Valgus/Vargus	Right Valgus/Vargus
Standard	14514100	14514200
Tall	14514300	14514400

OA Double Upright

- Features a unique thigh cuff that is flexible yet durable; the three strapping system enables optimal suspension, protection and comfort by moving with the muscles throughout range of motion
- Thigh cuff accommodates changes due to muscle atrophy, hypertrophy, weight gain or loss, and can be adjusted to accommodate larger patients
- Strong and lightweight fiber-reinforced composite frame
- Interchangeable condyle pads for OA unloading, increased comfort, stability and pain control
- Numbered strapping system for intuitive donning and doffing
- Silicone gripped calf strap liner to help prevent migration
- Optional flexion/extension stops for range of motion control

Indications: Mild to moderate unicompartmental osteoarthritis

HCPCS Code:
L1845 or K0901



Size	6" Above Mid-Patella	Mid-Patella	6" Below Mid-Patella	Left Valgus/Vargus	Right Valgus/Vargus
Small	15½" - 18½"	11" - 13"	12¼" - 13¾"	14515105	14515205
Medium	18½" - 21"	13" - 15"	13¼" - 15"	14515106	14515206
Large	21" - 23½"	15" - 17"	14¼" - 15¾"	14515107	14515207
X-Large	23½" - 25"	17" - 19"	15" - 17"	14515108	14515208
XX-Large	25" - 28"	19" - 22"	17" - 19"	14515109	14515209



HCPCS Code:
L1845 or K0901

Low-Contact

Size	6" Above Mid-Patella	Mid-Patella	6" Below Mid-Patella	Left Valgus/Vargus	Right Valgus/Vargus
Small	14½" - 17"	11½" - 13"	11½" - 13"	14533105	14533205
Medium	17" - 19½"	13" - 15"	13" - 14½"	14533106	14533206
Large	19½" - 22½"	15" - 17"	14½" - 16½"	14533107	14533207
X-Large	22½" - 25½"	17" - 19"	16½" - 18½"	14533108	14533208
XX-Large	25½" - 30"	19" - 22"	18½" - 22"	14533109	14533209

Contact

Size	6" Above Mid-Patella	Mid-Patella	6" Below Mid-Patella	Left Valgus/Vargus	Right Valgus/Vargus
Small	14½" - 17"	11½" - 13"	11½" - 13"	14534105	14534205
Medium	17" - 19½"	13" - 15"	13" - 14½"	14534106	14534206
Large	19½" - 22½"	15" - 17"	14½" - 16½"	14534107	14534207
X-Large	22½" - 25½"	17" - 19"	16½" - 18½"	14534108	14534208
XX-Large	25½" - 30"	19" - 22"	18½" - 22"	14534109	14534209

HCPCS Code:
L1845 or K0901



Functional ACL

- Features a unique thigh cuff that is flexible yet durable and enables optimal suspension, protection and comfort by moving with the muscles throughout the range of motion (for LOW-CONTACT sports/use only)
- Aluminum frame is strong yet customizable
- Also comes in full-frame design (for CONTACT sports/use)
- Posterior calf cuff reduces tibial pressure
- Interchangeable condyle pads for increased comfort, stability, and protection
- Numbered strapping system for intuitive application and removal
- Silicone gripped calf strap liner to help prevent migration
- Optional flexion/extension stops for range of motion control

Indications: Mild to moderate ACL, PCL, MCL, and LCL sprains; Combined ligament instabilities



Warrior® Recovery Post-Operative Knee Brace

- Easy-to-use, lever-action, micro-adjustable telescoping struts for length adjustments from 18" to 29"
- Four convenient preset drop lock settings, 0°, 15°, 30°, 45°
- Enhanced range-of-motion (ROM) hinge interface for easy setting and patient compliance
- ROM is set in 10° increments
- Soft, circumferential and trimmable straps plus removable Y-tab feature optimizes individual patient fit and comfort
- Inner pad panel adjustment (2mm hex wrench, not included) to avoid incision sites and provide a more customized fit
- Customizable anti-migration strap pads for superior fixture and reduced slippage
- Integrated "bend zone" on lower hinge arms allows for customization and accommodation to varying patient anatomy
- Laminate reinforced padding to support prolonged use and maintain pad integrity
- Extended padding at malleoli to prevent irritation and improve comfort during use

Indications: Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

Size	Adjustment Length	Product #
Universal	18" - 29"	KB9000-01

HCPCS Code:
L1832 or L1833

Warrior® is a registered trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroval.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Transition Telescoping

- Telescoping struts allow for length adjustment
- Clip buckles for easy application and removal
- Soft circumferential straps for comfort and even compression
- Drop-lock design offers easy access for pin adjustment
- Range of motion can be set in 10° increments
- Dial can be locked to encourage patient compliance
- Softgoods available in full or cool foam options

Indications: Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Full Foam	27" - 18"	KB7000-00
Cool Foam	27" - 18"	KB7000-01

**Measure circumference 6" above mid-patella*

HCPCS Code:
L1832 or L1833



Transition


- Drop-lock design offers easy access for pin adjustment
- Range of motion can be set in 10° increments
- Superior and inferior score zones for length adjustment
- Dial can be locked to encourage patient compliance
- Circumferential loop lock straps allow for even compression
- Trimmable softgoods available in full or cool foam options

Indications: Post-operative stabilization and rehabilitation; Post total knee replacement; Post ORIF; Fractures of the femur, tibia, fibula or patella; Ligament reconstructive procedures; Meniscal repairs; Post-fracture management

Description	Length Adjustment	Product #
Full Foam	23" - 17"	KB5000-00
Cool Foam	23" - 17"	KB5000-01
Full Foam, X-Long	26" - 20"	KB5000-02

**Measure circumference 6" above mid-patella*

HCPCS Code:
L1832 or L1833



Warrior® II

- Classic pull-up design, wrap-around design, and new half-wrap design to optimize individual patient fit and comfort
- Integrated anti-migration strip (pull-up designs) for superior fixture and reduced slippage
- Flex-tab feature to prevent hinge digging or irritation
- Plaited stretch-nylon material (engineered for swimwear) helps wick away moisture and promote aeration behind the leg and popliteal region
- Waffled spacer fabric promotes breathability and a cooler environment
- Removable condyle pads for increased comfort and stability
- Easy-grip silicone pull tabs for easier donning and doffing
- Two detachable, medial/lateral half-moon shaped buttresses (pull-up style) and one detachable, inferior half-moon shaped buttress (half-wrap) for patellar control
- Optional flexion/extension stops (0° to 100° in 10° increments) for range-of-motion control
- Optional “hyper control” strapping to help prevent hyperextension (regular designs)

Indications: Mild to moderate ACL, PCL, MCL, LCL; Combined ligament instabilities; Patellar subluxation

HCPCS Code: L1832 or L1833

Size	Measurement	Pull-Up Long	Pull-Up Short	Wrap-Around Long	Wrap-Around Short	Half-Wrap Long	Half-Wrap Short
X-Small	13" - 15½"	14011004	14021004	14012004	14022004	14013004	14023004
Small	15½" - 18½"	14011005	14021005	14012005	14022005	14013005	14023005
Medium	18½" - 21"	14011006	14021006	14012006	14022006	14013006	14023006
Large	21" - 23½"	14011007	14021007	14012007	14022007	14013007	14023007
X-Large	23½" - 26½"	14011008	14021008	14012008	14022008	14013008	14023008
XX-Large	26½" - 29½"	14011009	14021009	14012009	14022009	14013009	14023009
XXX-Large	29½" - 32"	14011010	14021010	14012010	14022010	14013010	14023010

*Measure circumference 6" above mid-patella

Warrior® is a registered trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.derooyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Warrior®

- "Genucentric" hinge design tracks the knee's instantaneous center of motion throughout range of motion
- Breathable Tri-tex material helps keep skin cool by wicking away moisture
- Lined popliteal region helps prevent material from bunching
- Pull-up or wrap-around design available
- Two half-moon shaped buttresses can be aligned in desired position for patellar control (pull-up style only)
- Optional extension (0°, 5°, 10° and 15°) and flexion (0°, 25°, 45°, 75° and 100°) stops
- Extension tabs at thigh and calf region accommodate circumferential variances (pull-up style only)
- Features pull tabs to aid in donning

Indications: Mild to moderate ACL, PCL, MCL, LCL; Combined ligament instabilities; Patellar subluxation

HCPCS Code: L1832 or L1833

Size	Measurement	Pull-Up Regular	Wrap-Around Regular	Pull-Up Short	Wrap-Around Short
Small	15½" - 18"	14612005	14611005	14582005	14581005
Medium	18" - 20½"	14612006	14611006	14582006	14581006
Large	20½" - 23"	14612007	14611007	14582007	14581007
X-Large	23" - 25½"	14612008	14611008	14582008	14581008
XX-Large	25½" - 28"	14612009	14611009	14582009	14581009
XXX-Large	28" - 30½"	14612010	14611010	14582010	14581010

*Measure circumference 6" above mid-patella



Hypercontrol®

- Tri-tex construction serves as a breathable alternative to neoprene
- Anatomical design features unique strapping system at the back of the knee to assist in preventing hyperextension
- Medial and lateral hinges help prevent hyperextension of the knee
- Pull-up or wrap-around design available
- Dual adjustable buttresses for maximum patella stabilization (pull-up style only)
- Extension tabs at the thigh and calf accommodate circumferential variances (pull-up style only)
- Superior adjustability for patient comfort

Indications: Grade I or II collateral ligament sprain; Mild ACL/PCL sprain; Grade I varus and valgus instability

Size	Measurement	Pull-Up	Wrap-Around
X-Small	13½" - 15½"	14630004	—
Small	15½" - 18"	14630005	14910005
Medium	18" - 20½"	14630006	14910006
Large	20½" - 23"	14630007	14910007
X-Large	23" - 25½"	14630008	14910008
XX-Large	25½" - 28"	14630009	14910009
XXX-Large	28" - 30½"	14630010	14910010

HCPCS Code: L1812

*Measure circumference 6" above mid-patella

Warrior® & Hypercontrol® are registered trademarks of DeRoyal Industries, Inc.



HCPCS Code: L1820

Deluxe Hinged Knee Support

- ¼" neoprene construction provides warmth, compression and support
- Medial and lateral hinges with range of motion at 0°, 15°, 30° and 45° to limit extension
- Covered hinges allow brace to be worn during contact sports
- Elastic straps above and below knee help prevent migration
- Superior felt buttress helps control patella

Indications: Grade I collateral ligament sprain; Mild ACL/ PCL sprain; Grade I varus and valgus instability; Patellar subluxation; Osteoarthritis

Size	Measurement	Product #
Small	13" - 14"	NE7722-72
Medium	14" - 15"	NE7722-73
Large	15" - 17"	NE7722-74
X-Large	17" - 18"	NE7722-75
XX-Large	18" - 20"	NE7722-76
XXX-Large	20" - 22"	NE7722-77

**Measure circumference below patella*

Hinged Knee Support

- Neoprene construction provides warmth, compression and support
- Hinges provide medial and lateral support
- Pull-up or wrap-around design available
- Features open popliteal for breathability
- Condyle pads for additional patient comfort

Indications: Grade I or II collateral ligament sprain; Mild ACL/ PCL sprain; Grade I varus and valgus instability

Size	Measurement	Pull-Up	Wrap-Around
X-Small	13½" - 15½"	NE7711-71	NE7711-81
Small	15½" - 18"	NE7711-72	NE7711-82
Medium	18" - 20½"	NE7711-73	NE7711-83
Large	20½" - 23"	NE7711-74	NE7711-84
X-Large	23" - 25½"	NE7711-75	NE7711-85
XX-Large	25½" - 28"	NE7711-76	NE7711-86

**Measure circumference 6" above mid-patella*



HCPCS Code: L1820



Full Length Patella Stabilizer

- 1/8" neoprene construction provides warmth, compression and support
- Also available in breathable Tri-tex material
- Full-length design for maximum support and compression
- Unique strapping system helps pull patella into desired position
- Horseshoe-shaped pad can be placed on medial or lateral side of patella
- Medial and lateral hinges for additional ligament support (Tri-tex only)

Indications: Medial or lateral patellar subluxation; Chondromalacia; Patellofemoral tracking dysfunction

HCPCS Code: L1812

Size	Measurement	Neoprene Right	Neoprene Left	Tri-tex Right	Tri-tex Left
X-Small	13½" - 15½"	NEP779371R	NEP779371L	14840204	14840104
Small	15½" - 18"	NEP779372R	NEP779372L	14840205	14840105
Medium	18" - 20½"	NEP779373R	NEP779373L	14840206	14840106
Large	20½" - 23"	NEP779374R	NEP779374L	14840207	14840107
X-Large	23" - 25½"	NEP779375R	NEP779375L	14840208	14840108
XX-Large	25½" - 28"	NEP779376R	NEP779376L	14840209	14840109

*Measure circumference 6" above mid-patella

Concise Patella Stabilizer



- 1/8" neoprene construction provides warmth, compression and support
- Also available in breathable Tri-tex material
- Low-profile design allows for increased range of motion during activity
- Unique strapping system helps pull patella into desired position
- Universal right/left

Indications: Medial or lateral patellar subluxation; Chondromalacia; Patellofemoral tracking dysfunction

Size	Measurement	Neoprene	Tri-tex
X-Small	13½" - 15½"	NEP7790-71	14750004
Small	15½" - 18"	NEP7790-72	14750005
Medium	18" - 20½"	NEP7790-73	14750006
Large	20½" - 23"	NEP7790-74	14750007
X-Large	23" - 25½"	NEP7790-75	14750008
XX-Large	25½" - 28"	NEP7790-76	14750009

**Measure circumference 6" above mid-patella*

Patella Band



- Constructed of neoprene for durability
- Unique bi-lateral strapping allows for greater control of buttress tension over patellar tendon
- Pull-up design helps promote compliance and reduce migration
- Low-profile design provides patient comfort

Indications: Patella tendonitis; Osgood-Schlatter's disease; Jumper's knee

Size	Measurement	Product #
X-Small	11" - 12"	14640004
Small	12" - 13"	14640005
Medium	13" - 14"	14640006
Large	14" - 15"	14640007
X-Large	15" - 16"	14640008
XX-Large	16" - 17"	14640009

**Measure circumference 6" above mid-patella*



Deluxe Knee Support with Trimmable Buttress

- 1/8" neoprene construction provides warmth, compression and support
- Anterior strapping above and below knee helps prevent migration while medial/lateral stays help provide additional stability
- Enclosed tubular buttress can be trimmed for desired patella position
- Available in open or closed popliteal

Indications: Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Patellar dislocation; Quadricep tendonitis; Plica syndrome

Size	Measurement	Closed Popliteal	Open Politeal
X-Small	13½" - 15½"	NE7718-71	14882004
Small	15½" - 18"	NE7718-72	14882005
Medium	18" - 20½"	NE7718-73	14882006
Large	20½" - 23"	NE7718-74	14882007
X-Large	23" - 25½"	NE7718-75	14882008
XX-Large	25½" - 28"	NE7718-76	14882009
XXX-Large	28" - 30½"	NE7718-77	14882010

*Measure circumference 6" above mid-patella



Universal Knee Support

- 1/8" neoprene construction provides warmth, compression and support
- Anterior wrap around strapping above and below knee helps prevent migration
- Universal adult and pediatric sizes available

Indications: Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Quadricep tendonitis; Plica syndrome

Size	Measurement	Product #
Universal	up to 25"	NE7710-70
Pediatric	up to 12"	NE7710-71



Knee Support with Trimmable Buttness

- 1/8" neoprene construction provides warmth, compression and support
- Enclosed tubular buttness can be trimmed for desired patella position

Indications: Grade I collateral ligament sprain; Patellar subluxation; Chondromalacia; Patellar tendonitis; Patellofemoral tracking dysfunction; Quadriceps tendonitis; Plica syndrome

Size	Measurement	Product #
X-Small	13½" - 15½"	NE7717-71
Small	15½" - 18"	NE7717-72
Medium	18" - 20½"	NE7717-73
Large	20½" - 23"	NE7717-74
X-Large	23" - 25½"	NE7717-75
XX-Large	25½" - 28"	NE7717-76
XXX-Large	28" - 30½"	NE7717-77

**Measure circumference 6" above mid-patella*



Visco Elastic Knee Support

- Visco elastic construction with four way stretch for maximum compression
- Silicone buttness option features spiral stays
- Closed patella design

Indications: Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

Size	Measurement	w/o Buttness	w/ Buttness
X-Small	13½" - 15½"	14761004	14801004
Small	15½" - 18"	14761005	14801005
Medium	18" - 20½"	14761006	14801006
Large	20½" - 23"	14761007	14801007
X-Large	23" - 25½"	14761008	14801008
XX-Large	25½" - 28"	14761009	14801009
XXX-Large	28" - 30½"	14761010	14801010

**Measure circumference 6" above mid-patella*

Open Patella Knee Support



- 1/8" neoprene construction provides warmth, compression and support
- Available with oval pad for additional patella protection
- Also available in breathable Tri-tex construction (w/o pad)
- Open patella for patient comfort

Indications: Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

Size	Measurement	w/o Pad	w/ Pad	Tri-tex
X-Small	13½" - 15½"	NE7702-71	NE7702-01	14811004
Small	15½" - 18"	NE7702-72	NE7702-02	14811005
Medium	18" - 20½"	NE7702-73	NE7702-03	14811006
Large	20½" - 23"	NE7702-74	NE7702-04	14811007
X-Large	23" - 25½"	NE7702-75	NE7702-05	14811008
XX-Large	25½" - 28"	NE7702-76	NE7702-06	14811009
XXX-Large	28" - 30½"	NE7702-77	NE7702-07	14811010
XXXX-Large	30½" - 33"	--	NE7702-08	--

*Measure circumference 6" above mid-patella

Closed Patella Knee Support



- 1/8" neoprene construction provides warmth, compression and support
- Available with oval pad for additional patella protection

Indications: Grade I collateral ligament sprain; Knee sprains and strains; Osteoarthritis; Mild joint effusion; Baker's cyst; Bursitis; Tendonitis

Size	Measurement	w/o Pad	w/ Pad
X-Small	13½" - 15½"	NE7701-71	NE7701-01
Small	15½" - 18"	NE7701-72	NE7701-02
Medium	18" - 20½"	NE7701-73	NE7701-03
Large	20½" - 23"	NE7701-74	NE7701-04
X-Large	23" - 25½"	NE7701-75	NE7701-05
XX-Large	25½" - 28"	NE7701-76	NE7701-06
XXX-Large	28" - 30½"	NE7701-77	NE7701-07

*Measure circumference 6" above mid-patella

Thigh Sleeve



- 1/8" neoprene construction provides warmth, compression and support
- Anatomically contoured to provide compression and support during healing

Indications: Quadriceps, adductor and hamstring strains

Size	Mid-Thigh Circumference	Product #
X-Small	18" - 19"	NE7734-71
Small	19" - 20"	NE7734-72
Medium	20" - 22"	NE7734-73
Large	22" - 25"	NE7734-74
X-Large	25" - 28"	NE7734-75
XX-Large	28" - 30"	NE7734-76

Calf Sleeve




- 1/8" neoprene construction provides warmth, compression and support
- Support adequately covers calf area for all patients
- Provides gentle compressive support to lower leg injuries

Indications: Calf (Gastrocnemius and Soleus) strains, Medial tibial stress syndrome (Shin splints), Peroneal compartment pain

Size	Calf Circumference	Product #
X-Small	12" - 13"	NE7744-71
Small	13" - 14"	NE7744-72
Medium	14" - 15"	NE7744-73
Large	15" - 17"	NE7744-74
X-Large	17" - 18"	NE7744-75
XX-Large	18" - 20"	NE7744-76
XXX-Large	20" - 22"	NE7744-77

**Measure circumference 6" above mid-patella*



Knee / Leg





Ankle Bracing

Ankle Stirrups

AFO

Walkers

Ankle/Foot



HCPCS Code: L1902



DeRoyal® Sports Orthosis Powered by The Boa® Closure System

- Powered by The Boa® Closure System for a secure fit
- Dial allows for micro-adjustability and on-the-fly adjustment
- Low friction lacing allows for even closure and a glove like fit
- Removable medial and lateral plastic stays provide rigid support
- Constructed of durable nylon for long-lasting wear
- Figure 8 strapping mimics athletic taping
- Low-profile design allows use in athletic shoe
- No seam under foot
- Universal right/left

Indications: Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Size	Ankle Circumference	w/o Stays	w/ Stays
X-Small	10" - 11"	AB3000-10	AB3100-10
Small	11" - 12"	AB3000-12	AB3100-12
Medium	12" - 13"	AB3000-14	AB3100-14
Large	13" - 14"	AB3000-16	AB3100-16
X-Large	14" - 15"	AB3000-18	AB3100-18
XX-Large	15" - 16"	AB3000-20	AB3100-20

*Measure circumference from base of heel around the ankle

DeRoyal® Sports Orthosis

HCPCS Code: L1902

- Constructed of durable nylon for long-lasting wear
- Figure 8 strapping mimics athletic taping
- Lace-up design for a secure fit
- Elastic cuff closure for added support and compression
- Low profile design allows use in athletic shoe
- Universal right/left

Indications: Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Description	Ankle Circumference	Product #
X-Small	10" - 11"	AB2900-10
Small	11" - 12"	AB2900-12
Medium	12" - 13"	AB2900-14
Large	13" - 14"	AB2900-16
X-Large	14" - 15"	AB2900-18
XX-Large	15" - 16"	AB2900-20

*Measure circumference from base of heel around the ankle



DeRoyal® is a registered trademark of DeRoyal Industries, Inc.
Boa® is a registered trademark of Boa Technology, Inc.

DeRoyal® Sports Ankle Brace 2



HCPCS Code: L1902

- Constructed of durable nylon for long-lasting wear
- Universal stabilization strap forces ankle into inversion or eversion
- Cuff closure for added support and compression
- Lace-up design for secure fit and primary stabilization
- Low profile design allows use in athletic shoe
- Universal right/left

Indications: Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Description	Ankle Circumference	Product #
X-Small	10" - 11"	15600004
Small	11" - 12"	15600005
Medium	12" - 13"	15600006
Large	13" - 14"	15600007
X-Large	14" - 15"	15600008
XX-Large	15" - 16"	15600009

*Measure circumference from base of heel around the ankle

DeRoyal® Sports Ankle Brace



HCPCS Code: L1902

- Constructed of durable nylon for long-lasting wear
- Figure 8 strapping configuration, in conjunction with medial to lateral stirrup straps, places foot in slightly everted position for maximum ankle stability
- Removable medial and lateral plastic stays with malleoli cutout provides comfort and support
- Lace-up design for a secure fit
- Low profile design allows use in athletic shoe

Indications: Chronic ankle instability; Prophylactic talocrural support; Prevention or treatment of acute ankle sprains; Subtalar joint sprain; Syndesmosis sprain

Size	Male Shoe Size	Female Shoe Size	Right	Left
X-Small	4 - 6	5 - 7	AB2801-20	AB2801-10
Small	6½ - 8	7½ - 9	AB2801-22	AB2801-12
Medium	8½ - 10	9½ - 11	AB2801-24	AB2801-14
Large	10½ - 12	11½ - 13	AB2801-26	AB2801-16
X-Large	12½ +	13½ +	AB2801-28	AB2801-18

DeRoyal® is a registered trademark of DeRoyal Industries, Inc.



HCPCS Code: L1902

Size	Male Shoe Size	Female Shoe Size	Right	Left
Small	< 8	< 9½	AB2150-10	AB2150-16
Medium	8 - 12	9½ - 13½	AB2150-12	AB2150-18
Large	12 +	13½ +	AB2150-14	AB2150-20



Element[®] Sport

Powered by The Boa[®] Closure System

- Patented heel control strapping system holds the calcaneus under the talus controlling subtalar inversion and talar rotation
- Powered by The Boa[®] Closure System for a secure glove-like fit that is easy for the patient to apply
- Dial allows on-the-fly adjustment
- Rigid medial and lateral uprights help prevent abnormal eversion and inversion
- Functional hinge design allows for utilization of the brace prophylactically
- Low profile design for a comfortable fit inside of shoes promoting patient compliance

Indications: Acute and chronic ankle sprains; Chronic ankle instability; Posterior tibial tendonitis; Peroneal tendonitis; Osteoarthritis



Element[®] is a registered trademark of DeRoyal Industries, Inc.
Boa[®] is a registered trademark of Boa Technology, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





HCPCS Code: L1902

Element[®] Sport Ankle Brace

- Patented heel control strapping system holds the calcaneus under the talus controlling subtalar inversion and talar rotation
- Rigid medial and lateral uprights help prevent abnormal eversion and inversion
- Functional hinge design allows for utilization of the brace prophylactically
- User friendly dynamic closure strap provides circumferential compression for an intimate fit
- Low profile design helps promote patient compliance
- Low profile foot plate for a comfortable fit inside of shoes

Indications: Acute and chronic ankle sprains; Chronic ankle instability; Posterior tibial tendonitis; Peroneal tendonitis; Osteoarthritis

Size	Male Shoe Size	Female Shoe Size	Right	Left
Small	< 8	< 9½	AB2140-10	AB2140-16
Medium	8 - 12	9½ - 13½	AB2140-12	AB2140-18
Large	12 +	13½ +	AB2140-14	AB2140-20



HCPCS Code: L1906

Element[®] Ankle Brace

- Patented heel control strapping system holds the calcaneus under the talus controlling subtalar inversion and talar rotation
- Rigid medial and lateral uprights prevent abnormal eversion and inversion and provide protection during activity
- Functional hinge design allows for utilization of the brace prophylactically, after acute injuries and during rehabilitation
- Rear foot stabilization is enhanced with heat moldable footplate to control pronation during loading and mid-stance phases of gait
- Pulley system strapping design provides optimal compression enhancing proprioception
- Wide base of support at impact reduces ankle inversion
- Unique design provides syndesmosis stabilization
- Concise version is 1½ shorter for lower profile fit
- Available in standard and wide calf cuff

Indications: Acute and chronic ankle sprains; Chronic ankle instability; Syndesmosis sprain; Posterior tibial tendonitis; Peroneal tendonitis; Osteoarthritis

Standard	Male Shoe Size	Female Shoe Size	Right	Left
Small	< 8	< 9½	AB2110-10	AB2110-16
Medium	8 - 12	9½ - 13½	AB2110-12	AB2110-18
Large	12 +	13½ +	AB2110-14	AB2110-20

Concise	Male Shoe Size	Female Shoe Size	Right	Left
Small	< 8	< 9½	AB2130-10	AB2130-16
Medium	8 - 12	9½ - 13½	AB2130-12	AB2130-18
Large	12 +	13½ +	AB2130-14	AB2130-20

Wide	Male Shoe Size	Female Shoe Size	Right	Left
Small	< 8	< 9½	AB2120-10	AB2120-16
Medium	8 - 12	9½ - 13½	AB2120-12	AB2120-18
Large	12 +	13½ +	AB2120-14	AB2120-20

*Measure calf circumference 8" from ground; Standard <10", Wide > 10"



HCPCS Code: L1906

DeRoyal® Functional Ankle Brace

- Anatomically designed uprights for a secure fit
- Hinged stirrup helps maintain support throughout range of motion
- Padded foam inlay for high energy absorption during activity
- Polyurethane foam padding for patient comfort and shock absorption
- Stretch loop straps for maximum support and edema control
- Universal right/left

Indications: Acute and chronic ankle sprains; Chronic instability; Prophylactic talocrural support

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	--	5 - 7	AB2000-10
Small	5 - 8	7½ - 10	AB2000-12
Medium	8½ - 12	10½ - 13	AB2000-14
Large	12½ - 16	13 +	AB2000-16



Neoprene Ankle Support

- ¼" neoprene construction provides warmth, compression and support
- Pull-on design
- Universal right/left

Indications: Mild ankle sprains and strains; Chronic ankle instability; Ankle edema; Osteoarthritis; Joint effusion

Size	Male Shoe Size	Female Shoe Size	Product #
Small	4½ - 6	5½ - 7	NE7732-72
Medium	6½ - 8	7½ - 9	NE7732-73
Large	8½ - 10	9½ - 11	NE7732-74
X-Large	10½ - 12	11½ - 13	NE7732-75
XX-Large	12½ +	13½ +	NE7732-76

DeRoyal® is a registered trademark of DeRoyal Industries, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Confor[®] Ankle Support

- Confor[®] open cell polyurethane foam provides a superior fit
- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Stretch loop straps provide compression to help reduce edema
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits in shoe
- Universal right/left
- Available in each or 6/case quantity

Indications: Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Each	Case
Standard	10"	AB2340-00	AB2340-06
Small	8½"	AB2341-00	AB2341-06

HCPCS Code: L4350

Ankle/Foot

Confor[®] is a registered trademark of Cabot Safety Intermediate, Inc.

Air/Gel Ankle Stirrup



HCPCS Code: L4350

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air/gel bladder can be placed in freezer for cold therapy
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits in shoe
- Universal right/left

Indications: Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2350-00
Small	8½"	AB2351-00

Air Ankle Stirrup



HCPCS Code: L4350

- Rigid medial and lateral shells minimize ankle rotation and inversion/eversion
- Air bladder allows for adjustment of air with inflation bulb
- Swivel straps allow even compression when applied in opposite directions
- Adjustable outside heel strap for maximum fit and comfort
- Low profile design fits in shoe
- Universal right/left

Indications: Acute and chronic ankle sprains; Osteoarthritis; Syndesmosis sprain

Size	Height	Product #
Standard	10"	AB2360-00
Small	8½"	AB2361-00
Pediatric	6"	AB2362-00



Deluxe Ankle Foot Orthosis

- Closed heel support provides additional security
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits in shoe
- Available in right/left design for best fit

Indications: Patients suffering from neurological or orthopedic conditions; Lower leg instabilities

Size	Male Shoe Size	Female Shoe Size	Right	Left
Small	6 – 7½	7½ - 9	AF1001-01	AF1001-04
Medium	8 – 9 ½	9½ - 11	AF1001-02	AF1001-05
Large	10 – 11½	11½ - 13	AF1001-03	AF1001-06
X-Large	12 +	13½ +	AF1001-07	AF1001-08



Ankle Foot Orthosis

- Open back and heel design
- Can be modified with heat gun and trimmed for customized fit
- Hook and loop closure for easy application
- Low profile design fits in shoe
- Universal in size

Indications: Patients suffering from neurological or orthopedic conditions; Lower leg instabilities

Size	Male	Female
Right	AF1000-01	AF1000-03
Left	AF1000-02	AF1000-04



Full Shell Walker 2

- Built in medial and lateral air pumps allow easy operation and instant air bladder control
- Air bladders accommodate swelling and bulky dressings
- Wide foot bed accommodates varying patient foot anatomies
- Textured tread, non-skid sole provides maximum shock absorption and grip
- Rocker bottom mimics natural gait
- Unique adjustable toe guard for increased protection
- Additional pads included for extra comfort where needed

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

Size	Male Shoe Size	Female Shoe Size	Standard	Low Top
X-Small	up to 4	up to 5	15500041	15500046
Small	4 - 7	5 - 8	15500042	15500047
Medium	7 - 10	8 - 11	15500043	15500048
Large	10 - 13	11 - 15	15500044	15500049
X-Large	13 +	15 +	15500045	15500050

HCPCS Code: L4361

Inline Pneumatic Walker 2



HCPCS Code: L4361

- High quality rigid outer skeleton for superior protection
- Dense foam liner to provide comfort and support
- Shock absorbing embedded sole tread improves patient comfort and helps prevent slips
- Rocker bottom mimics natural gait
- Plastic uprights help reduce ankle rotation and range of motion
- Large air bladder creates a larger surface for contact pressure to help control swelling and ensure a snug, comfortable fit

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

Size	Male Shoe Size	Female Shoe Size	Standard	Short
X-Small	up to 4½	up to 5½	15500031	15500036
Small	4½ - 7½	5½ - 8	15500032	15500037
Medium	7½ - 10½	8 - 11½	15500033	15500038
Large	10½ - 12½	11½ - 13½	15500034	15500039
X-Large	12½ +	13½ +	15500035	15500040

Pacesetter™ Air Walker

- Inflatable pump design for compression and stabilization
- Unique non-welded air bladder reduces chance of over-inflation to the point of bursting or tearing
- Wide foot plate for patient comfort and stability
- Rocker bottom mimics natural gait
- Offered in both standard and low top versions
- Heel strap helps prevent foot migration

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

Size	Male Shoe Size	Female Shoe Size	Standard Open Heel	Low Top Open Heel
Small	3 - 6½	4 - 7½	15780005	15800005
Medium	7 - 8½	8 - 10	15780006	15800006
Large	9 - 10½	10½ - 12½	15780067	15800067
X-Large	11 - 16	13 +	15780008	15800008
Youth	--	--	15780004	--



HCPCS Code: L4361

Ankle/Foot

Pacesetter™ is a trademark of DeRoyal Industries, Inc.



Pacesetter™ II Walker

- Rocker sole helps promote rolling gait pattern
- Double laminated foam wrap is breathable to prevent skin breakdown
- Wide foot plate for patient comfort and stability
- Available in open or closed heel design
- Offered in standard and low top versions

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

HCPCS Code: L4387

Size	Male Shoe Size	Female Shoe Size	Standard Open Heel	Standard Closed Heel	Low Top Open Heel
Small	3 - 6½	4 - 7½	CSG53510-A	CSG53550-A	CSG53520-A
Medium	7 - 8½	8 - 10	CSG53510-B	CSG53550-B	CSG53520-B
Large	9 - 10½	10½ - 12½	CSG53510-D	CSG53550-D	CSG53520-D
X-Large	11 - 16	13 +	CSG53510-C	CSG53550-C	CSG53520-C



Tracker EX Walker

- Lightweight design for patient comfort
- Low profile rocker bottom sole helps promote natural gait
- Wide foot plate for patient comfort and stability
- Closed heel design helps prevent heel migration
- Extra long wrap provides full coverage
- Available in standard and low top versions

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

HCPCS Code: L4387

Size	Male Shoe Size	Female Shoe Size	Standard	Low Top
X-Small	up to 5½	up to 4½	15590002	15580002
Small	6 - 8½	5 - 8½	15590005	15580005
Medium	9 - 10½	9 - 11½	15590006	15580006
Large	11 - 13	12 - 13½	15590067	15580067
X-Large	13½ +	14 +	15590008	15580008

Pacesetter™ is a trademark of DeRoyal Industries, Inc.



HCPCS Code: L4387

Grid Iron XR® Walker

- Features rigid metal uprights for maximum stability
- Low profile rocker bottom sole helps promote natural gait
- Wide foot plate for patient comfort and stability
- Closed heel design helps prevent heel migration
- Available in standard and low top versions

Indications: Moderate to severe ankle sprains and strains; Stable foot and/or ankle fractures; Syndesmosis sprain; Post ORIF; Achilles tendonitis

Size	Male Shoe Size	Female Shoe Size	Standard	Low Top
X-Small	up to 3½	up to 4½	15570004	15560004
Small	4 - 7	5 - 8½	15570005	15560005
Medium	7½ - 10	9 - 11	15570006	15560006
Large	10½ - 12	11½ - 13	15570007	15560007
X-Large	12½ +	13½ +	15570008	—



HCPCS Code: L4387

Grid Iron XR® R.O.M. Walker

- Range of motion is adjustable from fixed to 45°
- Rigid metal uprights for maximum stability
- Low profile rocker bottom sole helps promote natural gait
- Wide foot plate for patient comfort and stability
- Closed heel design helps prevent heel migration

Indications: Achilles tendon repair; Moderate to severe ankle sprains and strains; Stable foot and ankle fractures; Post ORIF; Achilles tendonitis

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	up to 3½	up to 4½	15510004
Small	4 - 7	5 - 8½	15510005
Medium	7½ - 10	9 - 11	15510006
Large	10½ - 12	11½ - 13	15510007
X-Large	12½ +	13½ +	15510008

Grid Iron XR® is a registered trademark of DeRoyal Industries, Inc.



JETSTREAM

QUICK START GUIDE

1. Remove the shipping container from the shipping box.
2. Remove the shipping container from the shipping box.
3. Remove the shipping container from the shipping box.
4. Attach the handle and lock.
5. Attach the handle and lock.
6. Attach the handle and lock.
7. Attach the handle and lock.
8. Attach the handle and lock.
9. Attach the handle and lock.
10. Attach the handle and lock.



DeRoyal
Improving Care.
Improving Business.

DeRoyal
Improving Care.
Improving Business.



Temperature Therapy Units
Therapy Blankets/Wraps
Thermal Supports
Accessories

Temperature Therapy



JetStream® Hot/Cold Temperature Therapy

- Patented ATC™ measures and regulates temperature to provide therapeutic hot or cold therapy without manual adjustment
- Intuitive LED light lets user know when in therapeutic zone and when to add ice or hot water
- Portable system suitable for any medical setting or treatment at home
- Lightweight, low-profile design and Quick Start Guide are patient friendly
- Both Hot and Cold Therapy in one unit
- Therapy Blankets optimize temperature distribution, insulation and comfort
- Blanket Wraps increase insulation and ensure precise placement of blankets

T700 Combo Units	Non-Sterile	Sterile
Chest/Lumbar	T725NS	T725ST
Knee/Shoulder	T735NS	T735ST
Knee/Shoulder (w/attached straps)	T736NS	T736ST
Shoulder/Knee (w/attached straps)	T737NS	T737ST

ATC™ is a trademark of DeRoyal Industries, Inc.
 JetStream® is a registered trademark of DeRoyal Industries, Inc.



CoolJet™ Cold Temperature Therapy

- Patented ATC™ measures and regulates temperature to provide therapeutic cold therapy without manual adjustment
- Intuitive LED light lets user know when in therapeutic zone and when to add ice
- Portable system suitable for any medical setting or treatment at home
- Therapy Blankets optimize temperature distribution, insulation and comfort
- Blanket Wraps increase insulation and ensure precise placement of blankets
- Intuitive screw top opening

T800 Combo Units	Non-Sterile	Sterile
Chest/Lumbar	T825NS	T825ST
Knee/Shoulder	T835NS	T835ST
Knee/Shoulder (w/attached straps)	T836NS	T836ST
Shoulder/Knee (w/attached straps)	T837NS	T837ST


ATC™ & JetStream™ are trademarks of DeRoyal Industries, Inc.



Blankets & Wraps

Blankets	Non-Sterile	Sterile
Arm/Leg	T5015NS	T5015ST
Hip/Lumbar	T5025NS	T5025ST
Knee/Shoulder	T5035NS T5035NS-10 (10 per/case)	T5035ST T5035ST-10 (10 per/case)
Knee/Shoulder-Large	T5035NS-L T5035NS-L-10 (10 per/case)	T5035ST-L T5035ST-L-10 (10 per/case)
Knee/Shoulder (w/ attached straps)	T5036NS T5036NS-10 (10 per/case)	T5036ST T5036ST-10 (10 per/case)
Knee/Shoulder-Large (w/attached straps)	T5036NS-L T5036NS-L-10 (10 per/case)	T5036ST-L T5036ST-L-10 (10 per/case)
Shoulder/Knee (w/ attached straps)	T5037NS T5037NS-10 (10 per/case)	T5037ST T5037ST-10 (10 per/case)
Shoulder/Knee - Large (w/attached straps)	T5037NS-L T5037NS-L-10 (10 per/case)	T5037ST-L T5037ST-L-10 (10 per/case)
Hand/Wrist	T5055NS	T5055ST
Cervical	T5075NS	T5075ST
Foot/Ankle	T5085NS	T5085ST

Foam Wraps	Product #
Chest/Lumbar	T6100
Knee	T6110
Shoulder	T6120
Hips	T6160



Neoprene Hot/Cold Therapy Wraps

- Constructed of comfortable neoprene material for an ergonomic fit
- Packaged with 2 gel packs to provide hot or cold therapy
- Packs can be uniquely positioned within wraps to provide more targeted treatment and improved compression
- Can be worn without gel packs for compression, support and warmth
- Shoulder and knee wraps available in each or 10/case quantity

Indications: Pain and soft tissue swelling as a result of surgery or trauma



Shoulder	Chest Circumference	Case	Product #
Universal	< 75"	9361-10	9361-00
Knee	Thigh Circumference	Case	Product #
Universal	< 25"	9363-10	9363-00
Back	Waist Circumference	Case	Product #
Universal	< 58"	—	9364-00
Ankle	Ankle Circumference	Case	Product #
Universal	< 20"	—	9365-00
Wrist		Case	Product #
Universal	-----	—	9366-00

Temperature Therapy

ActiveWrap® Thermal Supports



- Constructed of comfortable Active Medical Plush™ material for best ergonomic fit
- Gel packs provide hot or cold therapy
- Packs can be removed and uniquely positioned within wraps to provide more targeted treatment and improved compression
- Can be used without thermal packs for compression, support and warmth
- Post-op knee/leg and replacement packs are designed to stay cold for 2 plus hours (for cold therapy only)

Indications: Pain and soft tissue swelling as a result of surgery or trauma

Shoulder	Chest Circumference	Packs	Product #	
Small/Medium	< 38"	2	9371-00	
Large/X-Large	38" +	2	9371-01	
Back	Waist Circumference	Packs	Product #	
Small/Medium	< 35"	1	9374-00	
Large/X-Large	35" - 54"	1	9374-01	
Wrist/Hand		Packs	Product #	
Universal	-----	2	9376-00	
Knee/Leg	Thigh Circumference	Packs	Product #	
Small/Medium	< 16"	2	9373-00	
Large/X-Large	16" +	2	9373-01	
Post-Op Knee/Leg	Thigh Circumference	Packs	Product #	
Small/Medium	< 16"	2	9373-02	
Large/X-Large	16" +	2	9373-03	
Small/Medium	< 16"	4	9373-04	
Large/X-Large	16" +	4	9373-05	
X-Wide	16" +	3	9373-06	
Post-Op Replacement Pack		Qty	Product #	
Small/Medium	-----	48/Case	9378-01	
Large/X-Large	-----	24/Case	9378-02	
Large/X-Large	-----	12/Case	9378-03	
Foot/Ankle	Male Shoe Size	Female Shoe Size	Packs	Product #
Small/Medium	4 - 10	4 - 11	2	9375-00
Large/X-Large	10 +	11 +	2	9375-01
Hip		Packs	Product #	
Universal	-----	2	9372-00	

ActiveWrap® is a registered trademark of ActiveWrap, Inc.
Active Medical Plush™ is a trademark of ActiveWrap, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroval.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





Elastic Bandages

- Durable elastic construction
- Available in standard length (5 yds) and double length (10 yds)
- Packaged individually
- Non-sterile
- Available in hook and loop or clip closure
- Premium option features hook closure on both ends of roll

Size	Quantity	Premium Clip Closure	Premium Hook & Loop Closure	Econo Clip Closure	Econo Hook & Loop Closure
2" X 5 yds	10 rolls/case	EB2	EBV2	DEB2	DEBV2
3" X 5 yds	10 rolls/case	EB3	EBV3	DEB3	DEBV3
4" X 5 yds	10 rolls/case	EB4	EBV4	DEB4	DEBV4
6" X 5 yds	10 rolls/case	EB6	EBV6	DEB6	DEBV6
4" X 10 yds	6 rolls/case	HEB4	—	—	—
6" X 10 yds	6 rolls/case	HEB6	—	—	—



Chill-n-Wrap Ice Bag

- Constructed of soft outer cover with moisture absorbing middle layer
- Leak-proof plastic bladder prevents patient exposure
- Plastic closure clip is completely enclosed
- Conformable opening holds shape while ice is being added
- Packaged 10/case

Indications: Pain and soft tissue swelling as a result of surgery or trauma

Size	Pouch Dimensions	Strap Length	Product #
Large	6" x 16"	16½"	8726-50



Hot/Cold Pack

- Reusable for hot or cold therapy
- Packaged each or 20/case

Indications: Pain and soft tissue swelling as a result of surgery or trauma

Measurement	Each	Case
4½" x 10"	9380-00	9380-20

Coldstar™ is a trademark of Coldstar International.





Fracture Bracing

Post-Op/Cast Shoes

Miscellaneous Ambulatory Aids

Ambulatory Aids



Fracture Bracing

- Orthopedic grade thermoplastic construction
- Prefabricated and anatomically designed for best fit
- Durable, lightweight, adjustable and easy to remove for proper hygiene
- Allows for movement of muscles surrounding fracture site
- Hook and loop closures help maintain soft tissue compression

Indications: Treatment and management of diaphyseal fractures; Post surgical immobilization



Ulnar	Distal Circumference	Proximal Circumference	Product #	
Small	6½" - 7½"	7½" - 9"	FB4000-10	
Medium	7½" - 8½"	9" - 11"	FB4000-12	
Large	8½" - 9½"	11" - 13"	FB4000-14	
Hand/Thumb	Palm Width	Left	Right	
Small	2½" - 3"	FB5000-10	FB5000-16	
Medium	3" - 3½"	FB5000-12	FB5000-18	
Large	3½" - 4"	FB5000-14	FB5000-20	
Wrist	Palm Width	Left	Right	
Small	2½" - 3"	FB5000-30	FB5000-36	
Medium	3" - 3½"	FB5000-32	FB5000-38	
Large	3½" - 4"	FB5000-34	FB5000-40	
Humeral	Distal Circ.	Proximal Circ.	Regular	Extended
Small	7" - 9"	9" - 11"	FB1000-10	FB2000-10
Medium	9" - 11"	11" - 13"	FB1000-12	FB2000-12
Large	11" - 13"	13" - 15"	FB1000-14	FB2000-14
X-Large	13" - 15"	15" - 17"	FB1000-16	FB2000-16
Tibial	Calf Circ.	Ankle Circ.	Left	Right
Small	12½"	9"	FB3000-10	FB3000-20
Medium	14"	10"	FB3000-12	FB3000-22
Medium/Long	15"	10"	FB3000-14	FB3000-24
Large	16"	11"	FB3000-16	FB3000-26



Original Darby Shoe



- Semi-rigid metatarsal shank provides controlled forefoot motion
- Reinforced padded heel counter to help reduce slippage
- Padded upper for patient comfort
- Loop-lock closure for easy application
- Universal right/left

Indications: Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	DS1014-10	4 - 6	DS1013-10
Medium	8½ - 10	DS1014-12	6½ - 8	DS1013-12
Large	10½ - 12	DS1014-14	8½ - 10	DS1013-14
X-Large	12½ - 14	DS1014-16	—	—
Pediatric	12 - 1 (Child)	DS1012-00	—	—

Darco® MedSurg™ Shoe

- MetaShank™ protection provides rigid control under metatarsal heads
- Strapless forefoot closure expands to accommodate bulky dressings
- Square toe design serves as a bumper for additional protection
- Rocker sole reduces plantar pressure on forefoot and heel
- Ankle strap helps eliminate heel migration
- Universal right/left



Indications: Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	DS1016-10	4 - 6	DS1015-10
Medium	8½ - 10	DS1016-12	6½ - 8	DS1015-12
Large	10½ - 12	DS1016-14	8½ - 10	DS1015-14
X-Large	12½ - 14	DS1016-16	—	—

*Darco® is a registered trademark of Darco International, Inc.
MedSurg™ & MetaShank™ are trademarks of Darco International, Inc.*



Darco® OrthoWedge™ Shoe

- Reduces weight bearing pressure from the forefoot with 15° (dorsi-angle)
- Removable 10 mm insole allows customization for off-loading specific areas
- Square toe design serves as a bumper for additional protection
- Hook and loop closure for easy application
- Universal right/left

Indications: Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	4 - 7	HS1005-08
Small	6 - 8	7½ - 10	HS1005-10
Medium	8½ - 10	10½ - 13	HS1005-12
Large	10½ - 12	13½ +	HS1005-14
X-Large	12½ - 14	—	HS1005-16



Darco® SlimLine Cast Shoe

- Square toe design serves as a bumper for protection and provides room for larger casts
- Stretchable upper material helps boot conform to cast allowing precise fit
- Rocker sole promotes rolling gait pattern
- Padded insole to cushion patient's foot
- Strapless forefoot closure adjusts easily to cast size

Indications: Lower extremity cast protection

Size	Male Shoe Size	Female Shoe Size	Product #
X-Small	—	3 - 5	2039-01
Small	—	5½ - 8	2039-02
Medium	6 - 8	8½ -10	2039-03
Large	8½ - 11	10½ -13	2039-04
X-Large	11½ - 14	13+	2039-05

Size	Child Shoe Size	Product #
Child/Small	8 - 10½	2039-98
Child/Medium	11 - 1½	2039-99
Child/Large	2 - 4	2039-00



Post-Op Shoe

- Semi-rigid rubber sole for support
- Black nylon upper for comfort
- Loop-lock closure for easy application
- Universal right/left

Indications: Post-operative foot protection and support; Minor foot fractures

Size	Male Shoe Size	Product #	Female Shoe Size	Product #
Small	6 - 8	A152305	4 - 6	A152405
Medium	8½ - 10	A152306	6½ - 8	A152406
Large	10½ - 12	A152307	8½ - 10	A152407
X-Large	12½ - 14	A152308	—	—

Darco® is a registered trademark of Darco International, Inc.
Orthowedge™ is a trademark of Darco International, Inc.

DeRoyal Industries, Inc. | 200 DeBusk Lane | Powell, TN 37849 | www.deroyal.com | 1.800.DeRoyal (U.S.) | 865.938.7828 (Outside U.S.)





DeRoyal® Night Splint

- Dorsiflexion design applies stretch for effective pain relief
- Lightweight and low profile design for patient comfort
- Adjustable straps to vary amount of stretch
- Removable foam liner for laundering
- Optional foam wedges provided for additional stretch
- Kit features: Night splint, ice pack, exercise band, heel cups, foot control strap and self-management guide

Indications: Plantar fasciitis; Achilles tendonitis; Tendonosis

Size	Male Shoe Size	Female Shoe Size	Splint	Kit
Small	4 - 6	5 - 7	NS2000-10	SC2000-10
Medium	6½ - 9	7½ - 10	NS2000-12	SC2000-12
Large	9½ +	10½ +	NS2000-14	SC2000-14




HCPCS Code: L4397

Dorsal Night Splint

- Thermoplastic design holds foot in 90° neutral position
- Neoprene straps with hook and loop closure
- Plastic free bottom and non-skid transfer surface
- Low profile design allows for a cooler, more comfortable fit

Indications: Plantar fasciitis; Achilles tendonitis; Tendonosis

Size	Male Shoe Size	Female Shoe Size	Product #
Small/Medium	4 - 9	5 - 10	NS2000-20
Large/X-Large	9½ +	10½ +	NS2000-21



Silicone Heel Cups

- Constructed of medical grade silicone for long-lasting support and comfort
- Soft blue spot helps relieve heel strike fatigue and discomfort through shock absorption
- Packaged 1/Pair

Indications: Heel spurs; Bruises; Callosities; Plantar fasciitis and apophyseal injuries; Retrocalcaneal bursitis; Fat pad atrophy

Size	Product #
Small	HC2000-10
Medium	HC2000-12
Large	HC2000-14

DeRoyal® is a registered trademark of DeRoyal Industries, Inc.

Crutches



- Lightweight, sturdy aluminum construction
- Push-button height adjustment
- Rubber tips for maximum surface traction
- Padded underarm rest and hand grips for patient comfort
- Packaged 10 Pair/case

Size	Height	Product #
Short	4'8" - 5'2"	CR4000-LS
Regular	5'2" - 5'10"	CR2000-LS
Tall	5'10" - 6'6"	CR3000-LS

Folding Walkers



- Lightweight, sturdy aluminum construction
- Push-button height adjustment
- Push-button folding for easy storage and transport
- Rubber tips for maximum surface traction
- Padded hand grips for patient comfort
- Available with or without 3" wheels
- Packaged 2/case

Size	w/o Wheels	w/ Wheels
Universal	WLK02	WLK02-W

Walker Canes



- Lightweight and sturdy aluminum construction
- Push-button height adjustable from 30" - 40"
- Rubber tip for maximum surface traction
- Padded hand grip for patient comfort
- Packaged 20/case

Size	Product #
Universal	CA298

Patents

Element® - #7,828,758, #8,100,845, #8,317, 735

Ultralign® – #5,853,378; #5,967,998 and #6,478,759

Prolign® – #6,666,838

Information related to L-Codes or ICD-9 Codes for use in determining the proper Medicare or Medicaid reimbursement status for any DeRoyal product listed in this catalog, represents the best efforts of DeRoyal Industries, Inc. in determining relevant L-Codes and corresponding diagnoses for certain of DeRoyal's products. This information is provided as a no-cost service to customers of DeRoyal. The information is compiled from the descriptors and functions associated with the L-Codes and a review of ICD-9 Codes from publicly available information and anecdotal information obtained from sources familiar with the billing process as indicated above. However, only the appropriate Medicare authority can ultimately determine the appropriateness of reimbursement of a particular product, for a particular patient, given the appropriate diagnosis. DeRoyal strongly cautions the users of this information to verify the applicability of a given L-Code and/or ICD-9 Code in a particular situation prior to seeking reimbursement for a DeRoyal product.

Product Index

11201000.....	5	12850010.....	19	14022010.....	32	14640004.....	36
11202000.....	5	12850011.....	19	14023004.....	32	14640005.....	36
11211000.....	5	12870000.....	21	14023005.....	32	14640006.....	36
11212000.....	5	12982105.....	24	14023006.....	32	14640007.....	36
11221000.....	5	12982106.....	24	14023007.....	32	14640008.....	36
11222000.....	5	12982107.....	24	14023008.....	32	14640009.....	36
11231000.....	5	12982108.....	24	14023009.....	32	14750004.....	36
11232000.....	5	12982205.....	24	14023010.....	32	14750005.....	36
11241000.....	5	12982206.....	24	14514100.....	28	14750006.....	36
11242000.....	5	12982207.....	24	14514200.....	28	14750007.....	36
11251000.....	5	12982208.....	24	14514300.....	28	14750008.....	36
11252000.....	5	12991102.....	24	14514400.....	28	14750009.....	36
11261000.....	5	12991202.....	24	14515105.....	28	14761004.....	38
11262000.....	5	12992102.....	24	14515106.....	28	14761005.....	38
11271000.....	5	12992202.....	24	14515107.....	28	14761006.....	38
11272000.....	5	12993105.....	24	14515108.....	28	14761007.....	38
11281000.....	5	12993106.....	24	14515109.....	28	14761008.....	38
11282000.....	5	12993107.....	24	14515205.....	28	14761009.....	38
11291000.....	5	12993108.....	24	14515206.....	28	14761010.....	38
11292000.....	5	12993205.....	24	14515207.....	28	14801004.....	38
11311000.....	4	12993206.....	24	14515208.....	28	14801005.....	38
11312000.....	4	12993207.....	24	14515209.....	28	14801006.....	38
11313000.....	4	12993208.....	24	14533105.....	29	14801007.....	38
11421000.....	5	13400005.....	8	14533106.....	29	14801008.....	38
11422000.....	5	13400006.....	8	14533107.....	29	14801009.....	38
11423000.....	5	13400007.....	8	14533108.....	29	14801010.....	38
11431000.....	5	13400008.....	8	14533109.....	29	14811004.....	39
11432000.....	5	13502004.....	8	14533205.....	29	14811005.....	39
11433000.....	5	13502005.....	8	14533206.....	29	14811006.....	39
11441000.....	5	13502006.....	8	14533207.....	29	14811007.....	39
11442000.....	5	13502007.....	8	14533208.....	29	14811008.....	39
11443000.....	5	13502008.....	8	14533209.....	29	14811009.....	39
11451000.....	5	13502009.....	8	14534105.....	29	14811010.....	39
11452000.....	5	13503004.....	8	14534106.....	29	14840104.....	35
11461000.....	5	13503005.....	8	14534107.....	29	14840105.....	35
11462000.....	5	13503006.....	8	14534108.....	29	14840106.....	35
11471000.....	5	13503007.....	8	14534109.....	29	14840107.....	35
11472000.....	5	13503008.....	8	14534205.....	29	14840108.....	35
11481000.....	5	13503009.....	8	14534206.....	29	14840109.....	35
11481001.....	5	14011004.....	32	14534207.....	29	14840204.....	35
11481002.....	5	14011005.....	32	14534208.....	29	14840205.....	35
11481003.....	5	14011006.....	32	14534209.....	29	14840206.....	35
11482000.....	5	14011007.....	32	14581005.....	33	14840207.....	35
11482001.....	5	14011008.....	32	14581006.....	33	14840208.....	35
11482002.....	5	14011009.....	32	14581007.....	33	14840209.....	35
11482003.....	5	14011010.....	32	14581008.....	33	14882004.....	37
11483000.....	5	14012004.....	32	14581009.....	33	14882005.....	37
11483001.....	5	14012005.....	32	14581010.....	33	14882006.....	37
11483002.....	5	14012006.....	32	14582005.....	33	14882007.....	37
11483003.....	5	14012007.....	32	14582006.....	33	14882008.....	37
11620000.....	15	14012008.....	32	14582007.....	33	14882009.....	37
11630005.....	15	14012009.....	32	14582008.....	33	14882010.....	37
11630006.....	15	14012010.....	32	14582009.....	33	14910005.....	33
11630007.....	15	14013004.....	32	14582010.....	33	14910006.....	33
11630008.....	15	14013005.....	32	14611005.....	33	14910007.....	33
11631000.....	14	14013006.....	32	14611006.....	33	14910008.....	33
11632000.....	14	14013007.....	32	14611007.....	33	14910009.....	33
11650005.....	16	14013008.....	32	14611008.....	33	14910010.....	33
11650006.....	16	14013009.....	32	14611009.....	33	15500031.....	53
11650007.....	16	14013010.....	32	14611010.....	33	15500032.....	53
11690001.....	16	14021004.....	32	14612005.....	33	15500033.....	53
11690002.....	16	14021005.....	32	14612006.....	33	15500034.....	53
11690005.....	16	14021006.....	32	14612007.....	33	15500035.....	53
11690006.....	16	14021007.....	32	14612008.....	33	15500036.....	53
11690007.....	16	14021008.....	32	14612009.....	33	15500037.....	53
11690008.....	16	14021009.....	32	14612010.....	33	15500038.....	53
12840000.....	21	14021010.....	32	14630004.....	33	15500039.....	53
12850004.....	19	14022004.....	32	14630005.....	33	15500040.....	53
12850005.....	19	14022005.....	32	14630006.....	33	15500041.....	52
12850006.....	19	14022006.....	32	14630007.....	33	15500042.....	52
12850007.....	19	14022007.....	32	14630008.....	33	15500043.....	52
12850008.....	19	14022008.....	32	14630009.....	33	15500044.....	52
12850009.....	19	14022009.....	32	14630010.....	33	15500045.....	52

15500046.....	52	350XLL.....	25	9361-00.....	61	AB2140-20.....	47
15500047.....	52	350XLR.....	25	9361-10.....	61	AB2150-10.....	46
15500048.....	52	350XSL.....	25	9363-00.....	61	AB2150-12.....	46
15500049.....	52	350XSR.....	25	9363-10.....	61	AB2150-14.....	46
15500050.....	52	350XXLL.....	25	9364-00.....	61	AB2150-16.....	46
15510004.....	55	350XXLR.....	25	9365-00.....	61	AB2150-18.....	46
15510005.....	55	351LL.....	22	9366-00.....	61	AB2150-20.....	46
15510006.....	55	351LR.....	22	9371-00.....	62	AB2340-00.....	49
15510007.....	55	351ML.....	22	9371-01.....	62	AB2340-06.....	49
15510008.....	55	351MR.....	22	9372-00.....	62	AB2341-00.....	49
15560004.....	55	351SL.....	22	9373-00.....	62	AB2341-06.....	49
15560005.....	55	351SR.....	22	9373-01.....	62	AB2350-00.....	50
15560006.....	55	351XLL.....	22	9373-02.....	62	AB2351-00.....	50
15560007.....	55	351XLR.....	22	9373-03.....	62	AB2360-00.....	50
15570004.....	55	351XSL.....	22	9373-04.....	62	AB2361-00.....	50
15570005.....	55	351XSR.....	22	9373-05.....	62	AB2362-00.....	50
15570006.....	55	351XXLL.....	22	9373-06.....	62	AB2801-10.....	45
15570007.....	55	351XXLR.....	22	9374-00.....	62	AB2801-12.....	45
15570008.....	55	5023-01.....	23	9374-01.....	62	AB2801-14.....	45
15580002.....	54	5023-02.....	23	9375-00.....	62	AB2801-16.....	45
15580005.....	54	5023-03.....	23	9375-01.....	62	AB2801-18.....	45
15580006.....	54	5023-04.....	23	9376-00.....	62	AB2801-20.....	45
15580008.....	54	5023-05.....	23	9378-01.....	62	AB2801-22.....	45
15580067.....	54	5023-06.....	23	9378-02.....	62	AB2801-24.....	45
15590002.....	54	5023-07.....	23	9378-03.....	62	AB2801-26.....	45
15590005.....	54	5023-08.....	23	9380-00.....	63	AB2801-28.....	45
15590006.....	54	5023-09.....	23	9380-20.....	63	AB2900-10.....	44
15590008.....	54	5023-10.....	23	A125104.....	25	AB2900-12.....	44
15590067.....	54	5066-90.....	24	A125105.....	25	AB2900-14.....	44
15600004.....	45	5066-91.....	24	A125106.....	25	AB2900-16.....	44
15600005.....	45	5066-92.....	24	A125107.....	25	AB2900-18.....	44
15600006.....	45	5066-93.....	24	A125108.....	25	AB2900-20.....	44
15600007.....	45	5067-01.....	23	A125109.....	25	AB3000-10.....	44
15600008.....	45	5067-02.....	23	A125204.....	25	AB3000-12.....	44
15600009.....	45	5068-00.....	22	A125205.....	25	AB3000-14.....	44
157-01.....	9	5069-00.....	25	A125206.....	25	AB3000-16.....	44
157-02.....	9	5073-01.....	23	A125207.....	25	AB3000-18.....	44
157-03.....	9	5073-02.....	23	A125208.....	25	AB3000-20.....	44
157-04.....	9	5073-03.....	23	A125209.....	25	AB3100-10.....	44
157-05.....	9	5073-04.....	23	A152305.....	68	AB3100-12.....	44
157-06.....	9	5073-05.....	23	A152306.....	68	AB3100-14.....	44
15780004.....	53	5073-06.....	23	A152307.....	68	AB3100-16.....	44
15780005.....	53	5073-07.....	23	A152308.....	68	AB3100-18.....	44
15780006.....	53	5073-08.....	23	A152405.....	68	AB3100-20.....	44
15780008.....	53	5073-09.....	23	A152406.....	68	AF1000-01.....	51
15780067.....	53	5073-10.....	23	A152407.....	68	AF1000-02.....	51
15800005.....	53	5073-11.....	23	AB2000-10.....	48	AF1000-03.....	51
15800006.....	53	5073-12.....	23	AB2000-12.....	48	AF1000-04.....	51
15800008.....	53	5074-01.....	24	AB2000-14.....	48	AF1001-01.....	51
15800067.....	53	5074-02.....	24	AB2000-16.....	48	AF1001-02.....	51
158-01.....	9	5074-03.....	24	AB2110-10.....	47	AF1001-03.....	51
158-02.....	9	5074-04.....	24	AB2110-12.....	47	AF1001-04.....	51
158-03.....	9	5074-05.....	24	AB2110-14.....	47	AF1001-05.....	51
158-04.....	9	5074-06.....	24	AB2110-16.....	47	AF1001-06.....	51
158-05.....	9	5074-07.....	24	AB2110-18.....	47	AF1001-07.....	51
158-06.....	9	5074-08.....	24	AB2110-20.....	47	AF1001-08.....	51
159-01.....	7	5074-09.....	24	AB2120-10.....	47	CA298.....	70
159-02.....	7	5074-10.....	24	AB2120-12.....	47	CR2000-LS.....	70
159-03.....	7	8004-01.....	17	AB2120-14.....	47	CR3000-LS.....	70
2039-00.....	68	8004-02.....	17	AB2120-16.....	47	CR4000-LS.....	70
2039-01.....	68	8004-03.....	17	AB2120-18.....	47	CSG53510-A.....	54
2039-02.....	68	8004-04.....	17	AB2120-20.....	47	CSG53510-B.....	54
2039-03.....	68	8004-05.....	17	AB2130-10.....	47	CSG53510-C.....	54
2039-04.....	68	8004-06.....	17	AB2130-12.....	47	CSG53510-D.....	54
2039-05.....	68	8004-11.....	17	AB2130-14.....	47	CSG53520-A.....	54
2039-98.....	68	8004-12.....	17	AB2130-16.....	47	CSG53520-B.....	54
2039-99.....	68	8004-13.....	17	AB2130-18.....	47	CSG53520-C.....	54
350LL.....	25	8004-14.....	17	AB2130-20.....	47	CSG53520-D.....	54
350LR.....	25	8004-15.....	17	AB2140-10.....	47	CSG53550-A.....	54
350ML.....	25	8004-16.....	17	AB2140-12.....	47	CSG53550-B.....	54
350MR.....	25	8020-00.....	17	AB2140-14.....	47	CSG53550-C.....	54
350SL.....	25	8020-01.....	17	AB2140-16.....	47	CSG53550-D.....	54
350SR.....	25	8726-50.....	63	AB2140-18.....	47	DEB2.....	63


Product Index

DEB3.....	63	HS1005-12	68	NE7727-73.....	20	T5035ST-L	60
DEB4.....	63	HS1005-14	68	NE7727-74.....	20	T5035ST-L-10.....	60
DEB6.....	63	HS1005-16	68	NE7727-75.....	20	T5036NS.....	60
DEBV2	63	KB5000-00	31	NE7727-76.....	20	T5036NS-10	60
DEBV3	63	KB5000-01	31	NE7728-72.....	20	T5036NS-L	60
DEBV4	63	KB5000-02	31	NE7728-73.....	20	T5036NS-L-10.....	60
DEBV6	63	KB7000-00	31	NE7728-74.....	20	T5036ST	60
DS1012-00	67	KB7000-01	31	NE7728-75.....	20	T5036ST-10	60
DS1013-10	67	KB9000-01	30	NE7728-76.....	20	T5036ST-L	60
DS1013-12	67	NE7701-01.....	39	NE7730-70.....	21	T5036ST-L-10.....	60
DS1013-14	67	NE7701-02.....	39	NE7732-72.....	48	T5037NS.....	60
DS1014-10	67	NE7701-03.....	39	NE7732-73.....	48	T5037NS-10	60
DS1014-12	67	NE7701-04.....	39	NE7732-74.....	48	T5037NS-L	60
DS1014-14	67	NE7701-05.....	39	NE7732-75.....	48	T5037NS-L-10	60
DS1014-16	67	NE7701-06.....	39	NE7732-76.....	48	T5037ST	60
DS1015-10	67	NE7701-07.....	39	NE7734-71.....	40	T5037ST-10	60
DS1015-12	67	NE7701-71.....	39	NE7734-72.....	40	T5037ST-L	60
DS1015-14	67	NE7701-72.....	39	NE7734-73.....	40	T5037ST-L-10.....	60
DS1016-10	67	NE7701-73.....	39	NE7734-74.....	40	T5055NS.....	60
DS1016-12	67	NE7701-74.....	39	NE7734-75.....	40	T5055ST	60
DS1016-14	67	NE7701-75.....	39	NE7734-76.....	40	T5075NS.....	60
DS1016-16	67	NE7701-76.....	39	NE7744-71.....	40	T5075ST	60
EB2.....	63	NE7701-77.....	39	NE7744-72.....	40	T5085NS.....	60
EB3.....	63	NE7702-01.....	39	NE7744-73.....	40	T5085ST	60
EB4.....	63	NE7702-02.....	39	NE7744-74.....	40	T6100.....	60
EB5000-01.....	20	NE7702-03.....	39	NE7744-75.....	40	T6110.....	60
EB5000-02.....	20	NE7702-04.....	39	NE7744-76.....	40	T6120.....	60
EB5000-03.....	20	NE7702-05.....	39	NE7744-77.....	40	T6160.....	60
EB6.....	63	NE7702-06.....	39	NEP7790-71	36	T700.....	58
EB6000-01.....	18	NE7702-07.....	39	NEP7790-72	36	T725NS.....	58
EB6000-02.....	18	NE7702-08.....	39	NEP7790-73	36	T725ST	58
EBV2	63	NE7702-71.....	39	NEP7790-74	36	T735NS.....	58
EBV3	63	NE7702-72.....	39	NEP7790-75	36	T735ST	58
EBV4	63	NE7702-73.....	39	NEP7790-76	36	T736NS.....	58
EBV6	63	NE7702-74.....	39	NEP779371L.....	35	T736ST	58
FB1000-10.....	66	NE7702-75.....	39	NEP779371R.....	35	T737NS.....	58
FB1000-12.....	66	NE7702-76.....	39	NEP779372L.....	35	T737ST	58
FB1000-14.....	66	NE7702-77.....	39	NEP779372R.....	35	T800.....	59
FB1000-16.....	66	NE7710-70.....	37	NEP779373L.....	35	T825NS.....	59
FB2000-10.....	66	NE7710-71.....	37	NEP779373R.....	35	T825ST	59
FB2000-12.....	66	NE7711-71.....	34	NEP779374L.....	35	T835NS.....	59
FB2000-14.....	66	NE7711-72.....	34	NEP779374R.....	35	T835ST	59
FB2000-16.....	66	NE7711-73.....	34	NEP779375L.....	35	T836NS.....	59
FB3000-10.....	66	NE7711-74.....	34	NEP779375R.....	35	T836ST	59
FB3000-12.....	66	NE7711-75.....	34	NEP779376L.....	35	T837NS.....	59
FB3000-14.....	66	NE7711-76.....	34	NEP779376R.....	35	T837ST	59
FB3000-16.....	66	NE7711-81.....	34	NS2000-10	69	UPL3330-15.....	10
FB3000-20.....	66	NE7711-82.....	34	NS2000-12	69	UPL3430-15.....	10
FB3000-22.....	66	NE7711-83.....	34	NS2000-14	69	UPL3530-15.....	10
FB3000-24.....	66	NE7711-84.....	34	NS2000-20	69	UPL3630-15.....	10
FB3000-26.....	66	NE7711-85.....	34	NS2000-21	69	UPL3730-15.....	10
FB4000-10.....	66	NE7711-86.....	34	OPB3333-R	11	UPL4330-15.....	10
FB4000-12.....	66	NE7717-71.....	38	OPB3333-S.....	11	UPL4430-15.....	10
FB4000-14.....	66	NE7717-72.....	38	OPB3333-T.....	11	UPL4530-15.....	10
FB5000-10.....	66	NE7717-73.....	38	OPB3433-R	11	UPL4630-15.....	10
FB5000-12.....	66	NE7717-74.....	38	OPB3433-S.....	11	UPL4730-15.....	10
FB5000-14.....	66	NE7717-75.....	38	OPB3433-T.....	11	UPLP4330-15	10
FB5000-16.....	66	NE7717-76.....	38	OPB3533-R	11	UPLP4430-15	10
FB5000-18.....	66	NE7717-77.....	38	OPB3533-S.....	11	UPLP4530-15	10
FB5000-20.....	66	NE7718-71.....	37	OPB3533-T.....	11	UPLP4630-15	10
FB5000-30.....	66	NE7718-72.....	37	SC2000-10	69	UPLP4730-15	10
FB5000-32.....	66	NE7718-73.....	37	SC2000-12	69	UPT3333-15	11
FB5000-34.....	66	NE7718-74.....	37	SC2000-14	69	UPT3433-15	11
FB5000-36.....	66	NE7718-75.....	37	T5015NS.....	60	UPT3533-15	11
FB5000-38.....	66	NE7718-76.....	37	T5015ST	60	UPT3633-15	11
FB5000-40.....	66	NE7718-77.....	37	T5025NS.....	60	UPT3733-15	11
HC2000-10	69	NE7722-72.....	34	T5025ST	60	UPT4333-15	11
HC2000-12	69	NE7722-73.....	34	T5035NS.....	60	UPT4433-15	11
HC2000-14	69	NE7722-74.....	34	T5035NS-10	60	UPT4533-15	11
HEB4	63	NE7722-75.....	34	T5035NS-L	60	UPT4633-15	11
HEB6	63	NE7722-76.....	34	T5035NS-L-10	60	UPT4733-15	11
HS1005-08	68	NE7722-77.....	34	T5035ST	60	WLK02	70
HS1005-10	68	NE7727-72.....	20	T5035ST-10	60	WLK02-W	70






Head/Neck/Spine



Shoulder/Arm/Elbow/Wrist



Knee/Leg



Ankle/Foot



Temperature Therapy



Ambulatory Aids